

BD Product Name: **Care Coordination Engine**

Date of Critical or Security Patches: Jan 2022

Abstract: Critical or Security Patches – Jan 2022

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for Jan 2022.

The patches include updates that may adversely affect BD CCE. These patches should not be installed CCE Windows Server 2008 standard.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

Patch Name	Description	Patch ID	Notes
Security Only Quality Update for Windows Server 2008 for x64-based Systems.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5009601	1
Security Monthly Quality Rollup for Windows Server 2008 for x64-based Systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5009627	1

Security and Quality Rollup for .NET Framework 2.0, 3.0 for Windows Server 2008 SP2 for x64	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5008866	1
Security Only Update for .NET Framework 2.0, 3.0 for Windows Server 2008 SP2 for x64	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5008889	1
Security Only Update for .NET Framework 4.5.2 for Windows Server 2008 R2 for x64	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5008887	1
Security and Quality Rollup for .NET Framework 3.5, 4.5.2, 4.6, 4.6.1, 4.6.2, 4.7, 4.7.1, 4.7.2, 4.8	This security update addresses an issue where an unauthenticated attacker could cause a denial of service on an affected system.	KB5008868	1

- Note 1: This available update is not able to install due to an expired license issue on DoD 1.4.1, CCER 1.3.1, CCE 1.4.1 Windows 2008 Server 64-bit operating system.