

BD Product Name: **Care Coordination Engine**

Date of Critical or Security Patches: Oct 2021

Abstract: Critical or Security Patches – Oct 2021

### **Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for Oct 2021.

The patches include updates that may adversely affect BD CCE. These patches should not be installed CCE Windows Server 2008 standard.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

<b>Patch Name</b>	<b>Description</b>	<b>Patch ID</b>	<b>Notes</b>
Servicing Stack Update for Windows Server 2008 for x64-based Systems.	Install this update to resolve issues in Windows. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article for more information. After you install this item, you may have to restart your computer.	<a href="#">KB5006750</a>	1
Cumulative Security Update for Internet Explorer 11 for Windows Server 2008 R2 for x64-based systems.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	<a href="#">KB5006671</a>	1
Security Only Quality Update for Windows Server 2008 for x64-based Systems.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a	<a href="#">KB5006715</a>	1



	complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.		
Security Monthly Quality Rollup for Windows Server 2008 for x64-based Systems.	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	<a href="#">KB5006736</a>	1

- Note 1: This available update is not able to install due to an expired license issue on DoD 1.4.1, CCER 1.3.1, CCE 1.4.1 Windows 2008 Server 64-bit operating system.