

BD EpiCenter™ on Windows 10

August 2020

Critical and/or Security Patches January 2020 - April 2020.

**Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related and include updates through April 2020.

The patches include updates that may adversely affect BD EpiCenter. These patches will not be and should not be installed on BD EpiCenter at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

Patch Name	Patch Description	Patch ID	Month
------------	-------------------	----------	-------

**Notes**



BD EpiCenter™ on Windows 10

August 2020

Critical and/or Security Patches for January 2020-April 2020

### Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related and include updates through April 2020. These patches were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Month
Service Stack Update for Windows 10 IoT for x64-based Systems	This update makes quality improvements to the servicing stack, which is the component that installs Windows updates.	<a href="#">4550994</a>	April 2020
Cumulative Security Update for Windows 10 IoT for x64-based Systems	This security update includes quality improvements. Key changes include: <ul style="list-style-type: none"><li>Addresses an issue that causes <b>winspool.driv</b> to stop working and generates the error 0xc0000005 when users attempt to print from certain client applications.</li><li>Security updates to the Microsoft Scripting Engine, Windows Kernel, Windows App Platform and Frameworks, Microsoft Graphics Component, Windows Media, Windows Cloud Infrastructure, Windows Fundamentals, Windows Core Networking, Windows Update Stack, and the Microsoft JET Database Engine.</li></ul>	<a href="#">4550929</a>	April 2020
Security Update for Adobe Flash Player	This security update resolves vulnerabilities in Adobe Flash Player that is installed on any of the Windows operating systems	<a href="#">4537759</a>	February 2020

### Notes

1. Service Stack update 4550994 must be installed prior to other updates listed.



Advancing the world of health