

BD EpiCenter™ on Windows 7

June 2020

Critical and/or Security Patches January 2020 – May 2020.

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related and include updates through May 2020.

The patches include updates that may adversely affect BD EpiCenter. These patches will not be and should not be installed on BD EpiCenter at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

| Patch Name | Patch Description | Patch ID | Month |
|------------|-------------------|----------|-------|
|------------|-------------------|----------|-------|

Notes



BD EpiCenter™ on Windows 7

June 2020

Critical and/or Security Patches for January 2020 – May 2020

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related and include updates through May 2020. These patches were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

| Patch Name | Patch Description | Patch ID | Month |
|--|--|-------------------------|----------------|
| Security update for Windows 7 for x64-based Systems | This update corrects an issue in the Windows Update Agent that has the ZH-HK language pack installed. | 4474419 | September 2019 |
| Service Stack Update for Windows 7 for x64-based Systems | This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. | 4490628 | March 2019 |
| Security Update for Windows 7 for x64-based Systems | This update provides an additional set of licensing changes to enable installation of the ESU add-on key. This is one of the steps to prepare for installation of Extended Security Updates | 4538483 | May 2020 |
| Service Stack Update for Windows 7 for x64-based Systems | This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. | 4550735 | March 2020 |
| Service Stack Update for Windows 7 for x64-based Systems | This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. | 4550738 | April 2020 |
| Security monthly rollup for Windows 7. | This security update includes improvements and fixes that were a part of update KB4540688 (released March 10, 2020) and addresses the following issues: <ul style="list-style-type: none">• Security updates to Internet Explorer, the Microsoft Scripting Engine, Windows Kernel, Windows App Platform and Frameworks, Windows Graphics, Windows Media, Windows Cloud Infrastructure, Windows Fundamentals, Windows Core Networking, and the Microsoft JET Database Engine. | 4550964 | April 2020 |



Advancing the world of health

| | | | |
|----------------------|---|-------------------------|--------------|
| Security only update | <p>This security update includes quality improvements. Key changes include:</p> <ul style="list-style-type: none"> • Security updates to Windows Kernel, Windows App Platform and Frameworks, Windows Graphics, Windows Media, Windows Cloud Infrastructure, Windows Fundamentals, Windows Core Networking, and the Microsoft JET Database Engine. | 4550965 | January 2020 |
|----------------------|---|-------------------------|--------------|

Notes

1. Service Stack update 4490628, 4550735 and 4550738 must be installed prior to other updates listed.