

Security Patches: BD EpiCenter™

May 2025

BD has identified patches from Microsoft® that have been identified as critical or security related through January 2025 for BD EpiCenter™ on Windows LTSC 1809. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Microsoft® patches

Patch name	Description	Patch ID	Notes
Cumulative security update for Windows 10 1809	This update addresses security issues for the Windows 10 1809 operating system	KB5050008	
Servicing stack update for Windows 10 1809	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	KB5005112	This update must be installed prior to all other updates
Cumulative security update for .NET framework 3.5 and 4.7.2 for Win 10 1809.	This update addresses security issues for the .NET framework in Windows.	KB5049608	
SQL Server 2019 Security Update.	This update addresses security issues for the SQL Server database.	KB5046859	
SQL Server 2019 Security Update.	This update addresses security issues for the SQL Server database.	KB5046860	
SQL Server 2019 Security Update.	This driver update addresses security issues for the SQL Server database.	KB5037570	
SQL Server 2019 Security Update.	This driver update addresses security issues for the SQL Server database.	KB5040711	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

