

# Security Patches: BD FACSCanto™ 10-Color System

May 2023

## BD Biosciences - BD FACSCanto™ 10-Color Cell Analyzers (IVD and RUO)

### Microsoft™ and Third-Party Patches: March 2023 – May 2023

BD has identified patches from Microsoft™ that have been identified as critical or security related for March 2023 – May 2023. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

#### Microsoft™ patches: Windows™ 10 Enterprise LTSC

| Patch name                            | Description   | Patch ID                  | Notes |
|---------------------------------------|---|---------------------------|-------|
| March 2023 Monthly Cumulative Updates | Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v1809 | <a href="#">KB5023702</a> | 1, 2  |
| April 2023 Monthly Cumulative Updates | Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v1809 | <a href="#">KB5025229</a> | 1, 3  |
| May 2023 Monthly Cumulative Updates   | Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v1809 | <a href="#">KB5026362</a> | 1, 4  |



## Microsoft™ patches: Windows™ 10 Professional

| Patch name                            | Description  | Patch ID                  | Notes |
|---------------------------------------|--|---------------------------|-------|
| March 2023 Monthly Cumulative Updates | Security and quality updates for Microsoft™ Windows™ 10 Professional v22H2 | <a href="#">KB5023696</a> | 1, 2  |
| April 2023 Monthly Cumulative Updates | Security and quality updates for Microsoft™ Windows™ 10 Professional v22H2 | <a href="#">KB5025221</a> | 1, 3  |
| May 2023 Monthly Cumulative Updates   | Security and quality updates for Microsoft™ Windows™ 10 Professional v22H2 | <a href="#">KB5026361</a> | 1, 4  |

### Notes

1. Patches may require the system to reboot.
2. For further information from Microsoft regarding CVEs covered by the March Security updates, please refer to the following link – <https://msrc.microsoft.com/update-guide/releaseNote/2023-Mar>
3. For further information from Microsoft regarding CVEs covered by the April Security updates, please refer to the following link – <https://msrc.microsoft.com/update-guide/releaseNote/2023-Apr>
4. For further information from Microsoft regarding CVEs covered by the May Security updates, please refer to the following link – <https://msrc.microsoft.com/update-guide/releaseNote/2023-May>

BD FACSCanto™ Flow Cytometers are Class 1 Laser Products.

The BD FACSCanto™ Flow Cytometer is for In Vitro Diagnostic Use for up to six colors. Seven to ten colors are for Research Use Only.

BD Biosciences, San Jose, CA, 95131, U.S.  
877.232.8995

**bd.com**

BD, the BD Logo and FACSCanto are trademarks of Becton, Dickinson and Company or its affiliates.  
All other trademarks are the property of their respective owners. © 2023 BD. All rights reserved. BD-92316 (v1.0) 0623

