

BD FACSCanto™ with BD FACSLink™

Date of Critical or Security Patches: March 2019

### Microsoft® & Third-Party Patches

BD has reviewed patches from Microsoft® that have been identified as critical or security related for March 2019.

The March patches were confirmed to install without adversely affecting BD FACSCanto™ and BD FACSLink™ products. BD recommends installing the March patches to resolve the SMBv1 vulnerability.

Customers that maintain patches independent of BD delivery should ensure these actions are performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Notes
Monthly patches for Windows 7 SP1 and Windows Server 2008 R2 SP1	This security update resolves various vulnerabilities in Windows 7 SP1 and Windows Server 2008 R2 SP1. Note: Installing update 4489885 installs only security fixes that are listed here.	<a href="#">KB4489878</a>	1

### Notes

1. These patches contain an updated SMBv1 patch. January patches were confirmed to disable SMBv1 file shares created on Windows 7 SP1 x86 (32-bit) operating systems. Such file shares may be used for transferring data between FACSCanto instruments and FACSLink middleware

