

BD Product Security - Security Patch Recommendations

BD Biosciences Product Name: BD FACSLytic™ Cell Analyzers (IVD and RUO)

October / November / December 2020

Critical and/or Security Patches: The October 2020 security rollup/cumulative update was released on October 13, 2020. The November 2020 security rollup/cumulative update was released on November 10, 2020. The December 2020 security rollup/cumulative update was released on December 8, 2020.

Microsoft™ and Third-Party Patches

BD has identified patches from Microsoft that have been indicated as security related for the months of October, November, and December. The BD products listed above have been tested with security patches from the October, November, and December rollup/cumulative releases installed. The patches in these updates were not found to adversely affect operation of the listed BD products.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s). For products that are supported remotely by BD Service, patches will be applied according to customers' service agreements.

Windows™ 10 Enterprise:

Patch Name	Patch Description	Patch ID	Notes
October 2020 Monthly Cumulative Updates	Security and quality updates for Microsoft Windows 10 Enterprise (LTSC) v1809	KB4577668	1,2
November 2020 Monthly Cumulative Updates	Security and quality updates for Microsoft Windows 10 Enterprise (LTSC) v1809	KB4586793	1,3
December 2020 Monthly Cumulative Updates	Security and quality updates for Microsoft Windows 10 Enterprise (LTSC) v1809	KB4592440	1,4



Windows™ 10 Professional:

Patch Name	Patch Description	Patch ID	Notes
October 2020 Monthly Cumulative Updates	Security and quality updates for Microsoft Windows 10 Professional (SAC) v2004*	KB4579311	1,2
November 2020 Monthly Cumulative Updates	Security and quality updates for Microsoft Windows 10 Professional (SAC) v2004*	KB4586781	1,3
December 2020 Monthly Cumulative Updates	Security and quality updates for Microsoft Windows 10 Professional (SAC) v20H2**	KB4592438	1,4

*Windows 10 version 2004 released on May 27, 2020

**Windows 10 version 20H2 released on October 20, 2020

Notes

1. Patches may require the system to reboot.
2. For further information from Microsoft regarding CVEs covered by the October Security updates, please refer to the following link –
<https://msrc.microsoft.com/update-guide/releaseNote/2020-Oct>
3. For further information from Microsoft regarding CVEs covered by the November Security updates, please refer to the following link –
<https://msrc.microsoft.com/update-guide/releaseNote/2020-Nov>
4. For further information from Microsoft regarding CVEs covered by the December Security updates, please refer to the following link –
<https://msrc.microsoft.com/update-guide/releaseNote/2020-Dec>

BD FACSLyric™ Flow Cytometers are Class 1 Laser Products.

For BD FACSLyric CE-IVD, the BD FACSLyric™ Flow Cytometer with the BD FACSuite™ Clinical and BD FACSuite™ Applications are CE marked in compliance with the European In Vitro Diagnostic Medical Device Directive 98/79/EC.

For BD FACSLyric US-IVD, the BD FACSLyric™ Flow Cytometer is for In Vitro Diagnostic Use with BD FACSuite™ Clinical Application for up to six colors. The BD FACSLyric™ Flow Cytometer is for Research Use Only with BD FACSuite™ Application for up to 12 colors. Not for use in diagnostic or therapeutic procedures.

BD, the BD Logo, FACSLyric and FACSuite are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2021 BD. All rights reserved. BD-31978 (v1.0) 0721

