

# Security Patches: BD FACSLyric™ System

October 2025

## BD Biosciences - BD FACSLyric™ Cell Analyzers (IVD and RUO)

### Microsoft™ and Third-Party Patches: August 2025 – October 2025

BD has identified patches from Microsoft™ that have been identified as critical or security related for August 2025 – October 2025. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

**Microsoft™ patches: Windows™ 10 Professional GAC (Note: Windows™ 10 GAC reached Microsoft™ end of support on Oct 14, 2025)**

Patch name	Description	Patch ID	Notes
August 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Professional (GAC) v22H2	<a href="#">KB5063709</a>	1, 2
September 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Professional (GAC) v22H2	<a href="#">KB5065429</a>	1, 3
October 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Professional (GAC) v22H2	<a href="#">KB5066791</a>	1, 4

BD Biosciences, Milpitas, CA, 95035, U.S.  
877.232.8995

**bd.com**

BD, the BD Logo, BD FACSLyric and BD FACSuite are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2026 BD. All rights reserved. BD-168971 (v1.0) 0226



## Microsoft™ patches: Windows™ 10 Enterprise LTSC

Patch name	Description	Patch ID	Notes
August 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v1809	<a href="#">KB5063877</a>	1, 2
August 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v21H2	<a href="#">KB5063709</a>	1, 2
September 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v1809	<a href="#">KB5065428</a>	1, 3
September 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v21H2	<a href="#">KB5065429</a>	1, 3
October 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v1809	<a href="#">KB5066586</a>	1, 4
October 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v21H2	<a href="#">KB5066791</a>	1, 4

## Third-Party patches: Adobe Acrobat® Reader® and Google Chrome™ Browser

Third-Party Application	Version Description	URL	Notes
Adobe Acrobat® Reader®	Version 2025.001.20813 64-bit	<a href="#">Release Notes</a>	N/A
Google Chrome™ Browser	Version 141.0.7390.123 (Official Build) (64-bit)	<a href="#">Release Notes</a>	N/A

BD Biosciences, Milpitas, CA, 95035, U.S.  
877.232.8995

**bd.com**

BD, the BD Logo, BD FACSLyric and BD FACSuite are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2026 BD. All rights reserved. BD-168971 (v1.0) 0226



## Notes

1. Patches may require the system to reboot.
2. For further information from Microsoft regarding CVEs covered by the August Security updates, please refer to the following link – <https://msrc.microsoft.com/update-guide/releaseNote/2025-Aug>
3. For further information from Microsoft regarding CVEs covered by the September Security updates, please refer to the following link – <https://msrc.microsoft.com/update-guide/releaseNote/2025-Sep>
4. For further information from Microsoft regarding CVEs covered by the October Security updates, please refer to the following link – <https://msrc.microsoft.com/update-guide/releaseNote/2025-Oct>

The BD FACSLyric™ Flow Cytometer is a Class 1 Laser Product.

In the U.S., the BD FACSLyric™ Flow Cytometer is for In Vitro Diagnostic Use with BD FACSuite™ Clinical Application for up to six colors.

In the U.S., the BD FACSLyric™ Flow Cytometer is for Research Use Only with BD FACSuite™ Application for up to 12 colors. Not for use in diagnostic or therapeutic procedures.

In the E.U., the BD FACSLyric™ Flow Cytometer with the BD FACSuite™ Clinical and BD FACSuite™ Applications is an in vitro diagnostic medical device bearing a CE mark.

BD Biosciences, Milpitas, CA, 95035, U.S.  
877.232.8995

**bd.com**

BD, the BD Logo, BD FACSLyric and BD FACSuite are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2026 BD. All rights reserved. BD-168971 (v1.0) 0226

