

Security Patches: BD FACSymphony™ A3/A5 Systems

June 2025

BD Biosciences - BD FACSymphony™ A3/A5 Cell Analyzers

Microsoft™ and Third-Party Patches: April 2025 – June 2025

BD has identified patches from Microsoft™ that have been identified as critical or security related for April 2025 – June 2025. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Microsoft™ patches: Windows™ 10 Professional GAC

Patch name	Description	Patch ID	Notes
April 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Professional (GAC) v22H2	KB5055518	1, 2
May 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Professional (GAC) v22H2	KB5058379	1, 3
June 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Professional (GAC) v22H2	KB5060533	1, 4

BD Biosciences, Milpitas, CA, 95035, U.S.
877.232.8995

bd.com

BD, the BD Logo and BD FACSymphony are trademarks of Becton, Dickinson and Company or its affiliates.
All other trademarks are the property of their respective owners. © 2025 BD. All rights reserved. BD-155601 (v1.0) 0725



Microsoft™ patches: Windows™ 10 Enterprise LTSC

Patch name	Description	Patch ID	Notes
April 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v1809	KB5055519	1, 2
April 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v21H2	KB5055518	1, 2
May 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v1809	KB5058392	1, 3
May 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v21H2	KB5058379	1, 3
June 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v1809	KB5060531	1, 4
June 2025 Monthly Cumulative Updates	Security and quality updates for Microsoft™ Windows™ 10 Enterprise (LTSC) v21H2	KB5060533	1, 4

BD Biosciences, Milpitas, CA, 95035, U.S.
877.232.8995

bd.com

BD, the BD Logo and BD FACSymphony are trademarks of Becton, Dickinson and Company or its affiliates.
All other trademarks are the property of their respective owners. © 2025 BD. All rights reserved. BD-155601 (v1.0) 0725



Notes

1. Patches may require the system to reboot.
2. For further information from Microsoft regarding CVEs covered by the April Security updates, please refer to the following link – <https://msrc.microsoft.com/update-guide/releaseNote/2025-Apr>
3. For further information from Microsoft regarding CVEs covered by the May Security updates, please refer to the following link – <https://msrc.microsoft.com/update-guide/releaseNote/2025-May>
4. For further information from Microsoft regarding CVEs covered by the June Security updates, please refer to the following link – <https://msrc.microsoft.com/update-guide/releaseNote/2025-Jun>

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Class 1 Laser Product.

BD Biosciences, Milpitas, CA, 95035, U.S.
877.232.8995

bd.com

BD, the BD Logo and BD FACSymphony are trademarks of Becton, Dickinson and Company or its affiliates.
All other trademarks are the property of their respective owners. © 2025 BD. All rights reserved. BD-155601 (v1.0) 0725

