

Security Patches: BD FocalPoint™ Large Lab Server

December 2023

BD has identified patches from Microsoft® that have been identified as critical or security related up to November 2023. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

Microsoft® patches

Patch name	Description	Patch ID	Notes
Servicing stack update for Windows Server 2012	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	KB5032308	This update must be installed prior to all other updates
October 2023 rollup for Windows Server 2012 R2	Fixes various security issues for Windows Server 2012 R2.	KB5031407	
November 2023 rollup for Windows Server 2012 R2	Fixes various security issues for Windows Server 2012 R2.	KB5032249	
Update for extended security updates licensing.	This update provides an additional set of licensing changes to enable installation of the ESU add on key.	KB5017220	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.





BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.