

BD Kiestra™ InoqUA

Critical and/or Security Patches: October – December 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for October – December 2021. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Month
Security Update for Microsoft Windows 10 (1809) for X64-based Systems	Security improvements and fixes published by Microsoft for Windows 10 (1809).	KB5006672	OCT-2021
		KB5007206	NOV-2021
		KB5008218	DEC-2021
Service Stack Update (SSU) for Microsoft Windows 10 (1809) for X64-based Systems	Servicing stack updates fixes and improves the reliability of the update process to mitigate potential issues while installing.	KB5005112	OCT-2021
		KB5005112	NOV-2021
		KB5005112	DEC-2021

Notes

1. Customers using the System Control Unit (SCU) server system:
Your regional Customer Support Center will contact you to make an appointment to install the patches.

