

# Security Patches: BD Kiestra™ Inoqua +

March 2024

## BD Kiestra™ Inoqua +

Critical and/or Security Patches: January – March 2024.

BD has identified patches from Microsoft® that have been identified as critical or security related for January – March 2024. These patches were not found to adversely affect BD products and will be applied according to customers’ service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

### Microsoft® patches

Patch Name	Patch Description	Patch ID	Month
Security Update for Microsoft Windows 10 (1809) for X64-based Systems	Security improvements and fixes published by Microsoft for Windows 10 (1809).	<a href="#">KB5034127</a>	JAN-2024
		<a href="#">KB5034768</a>	FEB-2024
		<a href="#">KB5035849</a>	MAR-2024
Service Stack Update (SSU) for Microsoft Windows 10 (1809) for X64-based Systems	Servicing stack updates fixes and improves the reliability of the update process to mitigate potential issues while installing.	<a href="#">KB5005112</a>	JAN-2024
		<a href="#">KB5005112</a>	FEB-2024
		<a href="#">KB5005112</a>	MAR-2024

### Notes

1. Customers using the System Control Unit (SCU) server system:  
Your regional Customer Support Center will contact you to make an appointment to install the patches.
2. In order to improve the reliability of the update process, the Service Stack Update (SSU) should be installed first.
3. Security Updates are cumulative, therefore all previous security updates are included.

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

**bd.com**

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

