

BD Kiestra™ Inoqua+

January 2020

Critical and/or Security Patches for April – June 2019

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for July 2019. These patches were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Month
Security Update for Microsoft Windows Windows 7 SP1 for X86-based Systems	Security improvements and fixes published by Microsoft between April and June 2019 for Windows 7.	4503292	June
		4503269	
		4499164	May
		4499175	
		4493472	April
		4493448	
Security Update for Microsoft Windows Windows Server 2008 R2 SP1 for x64-based Systems	Security improvements and fixes published by Microsoft between April and June 2019 for Windows Server 2008R2.	4503292	June
		4503269	
		4499164	May
		4499175	
		4493472	April
		4493448	
Security Update for Microsoft Windows Windows Server 2012 R2 for x64-based Systems	Security improvements and fixes published by Microsoft between April and June 2019 for Windows Server 2012R2.	4503276	June
		4503290	
		4499151	May
		4499165	
		4493446	April
		4493467	

Notes

1. Customers using the System Control Unit (SCU) server system: Your regional Customer Support Center will contact you to make an appointment to install the patches



Advancing the world of health