

BD Kiestra™ Inoqua+

January 2020

Critical and/or Security Patches for July – September 2019

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for October 2019. These patches were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Month
Security Update for Microsoft Windows 7 SP1 for X86-based Systems	Security improvements and fixes published by Microsoft between July and September 2019 for Windows 7.	4516065	September
		4516033	
		4512506	August
		4512486	
		4507449	July
		4507456	
Security Update for Microsoft Windows Server 2008 R2 SP1 for x64-based Systems	Security improvements and fixes published by Microsoft between July and September 2019 for Windows Server 2008R2.	4516065	September
		4516033	
		4512506	August
		4512486	
		4507449	July
		4507456	
Security Update for Microsoft Windows Server 2012 R2 for x64-based Systems	Security improvements and fixes published by Microsoft between July and September 2019 for Windows Server 2012R2.	4516067	September
		4516064	
		4512488	August
		4512489	
		4507448	July
		4507457	

Notes

1. Customers using the System Control Unit (SCU) server system: Your regional Customer Support Center will contact you to make an appointment to install the patches

