

Security Patches: BD Kiestra™ ReadA

September 2022

BD Kiestra™ ReadA

Critical and/or Security Patches: July – September 2022.

BD has identified patches from Microsoft® that have been identified as critical or security related for July – September 2022. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Microsoft® patches

Patch Name	Patch Description	Patch ID	Month
Security Update for Microsoft Windows 10 (1809) for X64-based Systems	Security improvements and fixes published by Microsoft for Windows 10 (1809).	KB5015811	JULY-2022
		KB5016623	AUG-2022
		KB5017315	SEPT-2022
Service Stack Update (SSU) for Microsoft Windows 10 (1809) for X64-based Systems	Servicing stack updates fixes and improves the reliability of the update process to mitigate potential issues while installing.	KB5005112	JULY-2022
		KB5005112	AUG-2022
		KB5005112	SEPT-2022
Security Update for Microsoft Windows Server 2016 (1607) for x64-based Systems	Security improvements and fixes published by Microsoft for Windows Server 2016.	KB5015808	JULY-2022
		KB5016622	AUG-2022
		KB5017305	SEPT-2022
Service Stack Update (SSU) for Windows Server 2016 (1607) for x64-based Systems	Servicing stack updates fixes and improves the reliability of the update process to mitigate potential issues while installing.	KB5016058	JULY-2022
		KB5017095	AUG-2022
		KB5017396	SEPT-2022

Notes

1. Customers using the System Control Unit (SCU) server system:
Your regional Customer Support Center will contact you to make an appointment to install the patches.

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

