

## BD Kiestra™ TLA and BD Kiestra™ WCA

Critical and/or Security Patches: October – December 2021

### Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for October – December 2021. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Month
Security Update for Microsoft Windows 10 (1809) for X64-based Systems	Security improvements and fixes published by Microsoft for Windows 10 (1809).	<a href="#">KB5006672</a>	OCT-2021
		<a href="#">KB5007206</a>	NOV-2021
		<a href="#">KB5008218</a>	DEC-2021
Security Update for Microsoft Windows Server 2012 R2 for x64-based Systems	Security improvements and fixes published by Microsoft for Windows Server 2012 R2.	<a href="#">KB5006729</a>	OCT-2021
		<a href="#">KB5007255</a>	NOV-2021
		<a href="#">KB5008285</a>	DEC-2021
Monthly Rollup for Microsoft Windows Server 2012 R2 for x64-based Systems	This update package includes not only a given month's security updates, but also all past updates, going back to at least October 2016.	<a href="#">KB5006714</a>	OCT-2021
		<a href="#">KB5007247</a>	NOV-2021
		<a href="#">KB5008263</a>	DEC-2021
Security Update for Microsoft Windows Server 2016 (1607) for x64-based Systems	Security improvements and fixes published by Microsoft for Windows Server 2016.	<a href="#">KB5006669</a>	OCT-2021
		<a href="#">KB5007192</a>	NOV-2021
		<a href="#">KB5008207</a>	DEC-2021
Service Stack Update (SSU) for Microsoft Windows 10 (1809) for X64-based Systems	Servicing stack updates fixes and improves the reliability of the update process to mitigate potential issues while installing.	<a href="#">KB5005112</a>	OCT-2021
		<a href="#">KB5005112</a>	NOV-2021
		<a href="#">KB5005112</a>	DEC-2021
Service Stack Update (SSU) for Windows Server 2012 R2 for x64-based Systems	Servicing stack updates fixes and improves the reliability of the update process to mitigate potential issues while installing.	<a href="#">KB5005698</a>	OCT-2021
		<a href="#">KB5005698</a>	NOV-2021
		<a href="#">KB5005698</a>	DEC-2021
Service Stack Update (SSU) for Windows Server 2016 (1607) for x64-based Systems	Servicing stack updates fixes and improves the reliability of the update process to mitigate potential issues while installing.	<a href="#">KB5001403</a>	OCT-2021
		<a href="#">KB5001403</a>	NOV-2021
		<a href="#">KB5001403</a>	DEC-2021

Additionally, in this quarter a workaround has been released by VMware for the Apache Log4j vulnerability in VMware Vcenter software, which also is used by this product. The workaround will be applied by your regional Service organization during deployment of this patch cycle. For more information, please refer to the Apache Log4J Product Security bulletin on the [BD Trust Center](#) website.



## Notes

1. Customers using the System Control Unit (SCU) server system:  
Your regional Customer Support Center will contact you to make an appointment to install the patches.

