

BD MAX™

Date of Critical or Security Patches: April 2021

Critical Security Patches: October 2020-March 2021

Microsoft® & Third-Party Patches for Windows 10

BD has identified patches from Microsoft® that have been identified as critical or important security related updates thru March 2021. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Notes
Security Update for Secure Boot DBX	This security update makes improvements to Secure Boot DBX	KB4535680	
Adobe Flash Player	This security update resolves vulnerabilities in Adobe Flash Player that is installed	KB4580325	
Cumulative Security Update	<ul style="list-style-type: none">▪ Updates security for the Windows user interface.▪ Updates to improve security when Windows performs basic operations.▪ Updates to improve security when using Microsoft Office products.	KB5000803	
Servicing Stack Update	This update makes quality improvements to the servicing stack, which is the component that installs Windows updates.	KB5001078	1

Notes:

Servicing Stack update KB5001078 must be installed prior to other updates.

