

BD Multiprocessor Remote Station™

January 2019

### Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for October 2018.

- These patches were not found to adversely affect BD products. Customers that maintain patches independent of BD automated delivery should ensure these actions are performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest version of the BD offering

Patch Name	Patch Description	Patch ID	Notes
10-2018 Security Monthly Quality Rollup for <b>Windows 7</b> for <b>x64</b> -based Systems	Security updates to Windows Media Player, Windows Graphics, Microsoft Graphics Component, Windows Storage and Filesystems, Windows Kernel, and the Microsoft JET Database Engine.	<a href="#">KB4462923</a>	1

### Notes

1. Tested against BD Totalys Multiprocessor Workstation configured with Windows 7 Professional 64-bit 2009/6.1 (Build 7601: Service Pack 1).



Advancing the world of health