

# Security Patches: BD Phoenix™ M50

March 2026

BD has identified patches from Microsoft® that have been identified as critical or security related through December 2025. These patches were not found to adversely affect the BD Phoenix M50 and can be applied according to customers' service agreements. These patches will impact on the onboard operator's manual of the M50, making it illegible. The operator's manual is available online through the BD public facing website.

## Microsoft® patches

Patch name	Description	Patch ID	Notes
Security Update for SQL Server 2019.	This security update contains fixes and resolves vulnerabilities for the windows operating system	<a href="#">kb5068404</a>	
Servicing Stack Update	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	<a href="#">kb5031539</a>	Must be installed before other updates.
Cumulative Update for .NET Framework 3.5 and 4.8 for Windows 10 version 21h2	The primary purpose of this KB5066135 update is to address a security vulnerability in the .NET Framework	<a href="#">kb5066135</a>	

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

**bd.com**

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

