

BD Product Name: BD Pyxis™ Pyxis Anesthesia Station

Date of Critical or Security Patches: September 2019

Abstract: Security Update for SQL Server 2012 Service Pack 3 GDR (KB4057115) issue on BD Pyxis Pyxis Anesthesia Station

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for August 2019.

The patches include updates that may adversely affect BD Pyxis Anesthesia Station ES, These patches will not be and should not be installed on BD Pyxis Pyxis Anesthesia Station ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

Patch Name	Patch Description	Patch ID	Notes
4057115	Security Update for SQL Server 2012 Service Pack 3 GDR (KB4057115)	KB 4057115	1

Notes

1. KB 4057115 is declined for installation, findings appeared on ES Pyxis Anesthesia Station 1.4.4.1.

