

# Security Patches: BD Pyxis™ MedStation™ 3500

April 2023

BD has identified patches from Microsoft® that have been identified as critical or security related for April 2023 will not install on BD Pyxis™ MedStation™ 3500.

The patches include updates that may adversely affect BD Pyxis™ MedStation™ 3500. These patches will not be and should not be installed on BD Pyxis™ MedStation™ 3500 at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

## Microsoft® patches

Patch name	Description	Patch ID	Notes
2023-04 Security Only Quality Update for Windows Embedded Standard 7 for x86-based Systems	This update makes improvements for the following issues: <ul style="list-style-type: none"><li>By the March 1, 2023, order of the Arab Republic of Egypt, daylight savings time (DST) will resume on April 28, 2023, and end on October 27, 2023. This update applies to Egypt Standard Time - (UTC+02:00) Cairo.</li><li>Known compatibility issues exist with certain printer models which feature GDI printer drivers that do not completely adhere to GDI specifications.</li></ul>	<a href="#">KB5025277</a>	1
2023-04 Security Monthly Quality Rollup for Windows Embedded Standard 7 for x86-based Systems	This cumulative security update contains improvements that are part of update KB5023769 (released March 14, 2023). This update also makes improvements for the following issues: <ul style="list-style-type: none"><li>By the March 1, 2023, order of the Arab Republic of Egypt, daylight savings time (DST) will resume on April 28, 2023, and end on October 27, 2023. This update applies to Egypt Standard Time - (UTC+02:00) Cairo.</li><li>Known compatibility issues exist with certain printer models which feature GDI printer drivers that do not completely adhere to GDI specifications.</li></ul>	<a href="#">KB5025279</a>	1

Note:-

1. Due to environment issue patches were unable to be validated this month.

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

[bd.com](#)

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.



# Security Patches: BD Pyxis™ MedStation™ 4000

April 2023

BD has identified patches from Microsoft® that have been identified as critical or security related for April 2023 will not install on BD Pyxis™ MedStation™ 4000.

The patches include updates that may adversely affect BD Pyxis™ MedStation™ 4000. These patches will not be and should not be installed on BD Pyxis™ MedStation™ 4000 at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

## Microsoft® patches

Patch name	Description	Patch ID	Notes
Windows Malicious Software Removal Tool - v5.112	The Windows Malicious Software Removal Tool (MSRT) helps remove malicious software from computers.	<a href="#">KB890830</a>	1
2023-04 Security Only Quality Update for Windows Server 2008 for x86-based Systems	This update makes improvements for the following issues: <ul style="list-style-type: none"><li>By the March 1, 2023, order of the Arab Republic of Egypt, daylight savings time (DST) will resume on April 28, 2023, and end on October 27, 2023. This update applies to Egypt Standard Time - (UTC+02:00) Cairo.</li></ul>	<a href="#">KB5025273</a>	1
2023-04 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems	This cumulative security update contains improvements that are part of update KB5023755 (released March 14, 2023). This update also makes improvements for the following issues: <ul style="list-style-type: none"><li>By the March 1, 2023, order of the Arab Republic of Egypt, daylight savings time (DST) will resume on April 28, 2023, and end on October 27, 2023. This update applies to Egypt Standard Time - (UTC+02:00) Cairo.</li></ul>	<a href="#">KB5025271</a>	1

Note:-

1. Due to environment issue patches were unable to be validated this month.

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

[bd.com](#)

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

