

# Security Patches: BD Pyxis™ Medstation ES

July 2022

BD has identified patches from Microsoft® that have been identified as critical or security related for July 2022 will not install on BD Pyxis Medstation ES. The patch includes updates that may adversely affect BD Pyxis Medstation ES. This patch will not be and should not be installed on BD Pyxis™ Medstation ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: This patch only apply to the latest supported version(s) of the BD offering.

## Microsoft® patches

Patch name	Description	Patch ID	Notes
2022-07 Security Only Quality Update for Windows Embedded Standard 7 for x86-based Systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	<a href="#">KB5015862</a>	1

Note:

1. Application does not launch on Windows 7 after installing this patch.

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

**bd.com**

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.



# Security Patches: BD Pyxis™ Anesthesia ES

July 2022

BD has identified patches from Microsoft® that have been identified as critical or security related for July 2022 will not install on BD Pyxis Anesthesia ES. The patch includes updates that may adversely affect BD Pyxis Anesthesia ES. This patch will not be and should not be installed on BD Pyxis™ Anesthesia ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: This patch only apply to the latest supported version(s) of the BD offering.

## Microsoft® patches

Patch name	Description	Patch ID	Notes
2022-07 Security Only Quality Update for Windows Embedded Standard 7 for x86-based Systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	<a href="#">KB5015862</a>	1

Note:

1. Application does not launch on Windows 7 after installing this patch.

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

**bd.com**

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

