

BD Product Name: **BD Pyxis™ CII Safe ES**

Date of Critical or Security Patches: June 2021

Abstract: Critical or Security Patches – June 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for June 2021. The patches include updates that may adversely affect BD Pyxis™ CII Safe ES. These patches will not be and should not be installed on BD Pyxis™ CII Safe ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

Patch Name	Description	Patch ID	Notes
Security Update for SQL Server 2016 Service Pack 2 CU	This security update fixes the following issue: KB4583468 - Microsoft SQL Server elevation of privilege vulnerability	KB4583461	1
SQL Server 2016 Service Pack 2 Cumulative Update (CU) 16	CU16 for SQL Server 2016 Service Pack 2 upgraded all SQL Server 2016 Service Pack 2 instances and components installed through the SQL Server setup. CU16 can upgrade all editions and servicing levels of SQL Server 2016 Service Pack 2 to the CU16 level.	KB5000645	1
SQL Server 2016 Service Pack 2 Cumulative Update (CU) 15	CU15 for SQL Server 2016 Service Pack 2 upgraded all SQL Server 2016 Service Pack 2 instances and components installed through the SQL Server setup. CU15 can upgrade all editions and servicing levels of SQL Server 2016 Service Pack 2 to the CU15 level. For customers in need of additional installation	KB4577775	1



	options, please visit the Microsoft Download Center to download the latest Cumulative Update (https://support.microsoft.com/en-us/kb/957826)options, please visit the Microsoft Download Center to download the latest Cumulative Update (https://support.microsoft.com/en-us/kb/957826)		
--	--	--	--

Note:

1. After installation of this KB. DB is going into recovery mode. Bug raised for same. Bug Id: 1171472

BD Product Name: **BD Pyxis™ Supply**

Date of Critical or Security Patches: June 2021

Abstract: Critical or Security Patches – June 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for June 2021 will not install on BD Pyxis™ Supply.

The patch includes updates that may adversely affect Pyxis™ Supply. This patch will not be and should not be installed on BD Pyxis™ Supply at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: This patch only apply to the latest supported version(s) of the BD offering.

Patch Name	Description	Patch ID	Notes
2021-06 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems	This security update includes improvements and fixes that were a part of update KB5003210 (released May 11, 2021) and addresses the following issues: Security updates to Windows App Platform and Frameworks, Windows Cloud Infrastructure, Windows Authentication, Windows Fundamentals, Windows Storage and Filesystems, Windows HTML Platform, and Microsoft Scripting Engine.	KB5003661	1
2021-06 Cumulative Security Update for Internet Explorer 9 for Windows Server 2008 for x86-based systems	This security update resolves vulnerabilities in Internet Explorer.	KB5003636	1
2021-06 Security Only Quality Update for Windows Server 2008 for x86-based Systems	This security update includes quality improvements. Key changes include: Security updates to Windows App Platform and Frameworks, Windows Cloud Infrastructure, Windows Authentication, Windows Fundamentals, and Windows Storage and Filesystems.	KB5003695	1



2021-06 Security and Quality Rollup for .NET Framework 2.0, 3.0, 4.5.2, 4.6 and 4.6.2 for Windows Server 2008 SP2	This update is included in the Security and Quality Rollup that's dated June 8, 2021. Parts of this update were previously released in the rollup that's dated October 13, 2020.	KB5003782	1
---	--	---------------------------	---

Note:

1. Due to environment issue patch was unable to be validated this month.

BD Product Name: **BD Pyxis™ CIISafe™**

Date of Critical or Security Patches: June 2021

Abstract: Critical or Security Patches – June 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for June 2021. The patches include updates that may adversely affect BD Pyxis™ CII Safe ES. These patches will not be and should not be installed on BD Pyxis™ CII Safe ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

Patch Name	Description	Patch ID	Notes
Windows Malicious Software Removal Tool - v5.88 (KB890830)	After the download, this tool runs one time to check your computer for infection by specific, prevalent malicious software (including Blaster, Sasser, and Mydoom) and helps remove any infection that is found. If an infection is found, the tool will display a status report the next time that you start your computer. A new version of the tool will be offered every month. If you want to manually run the tool on your computer, you can download a copy from the Microsoft Download Center, or you can run an online version from microsoft.com. This tool is not a replacement for an antivirus product. To help protect your computer, you should use an antivirus product.	KB890830	1
2021-06 Cumulative Update for .NET Framework 4.8 for Windows 10 Version	Install this update to resolve issues in Windows. For a complete listing of the issues that are included in this	KB5003542	1



1607 for x64 (KB5003542)	update, see the associated Microsoft Knowledge Base article for more information. After you install this item, you may have to restart your computer.		
-----------------------------	---	--	--

Note:

1. System not installing the June cycle patches on the CIISafe Win 10 device and they are getting failed. (KB890830 and KB5003542). Bug raised for same.
Bug Id: 1221374

BD Product Name: **BD Pyxis™ Anesthesia™ ES**
 Date of Critical or Security Patches: June 2021
 Abstract: Critical or Security Patches – June 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for June 2021. The patches include updates that may adversely affect BD Pyxis™ MedStation™ ES. These patches will not be and should not be installed on BD Pyxis™ MedStation™ ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

Patch Name	Description	Patch ID	Notes
Cumulative update for windows10 version 1809for x64 based systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5003646	1
Cumulative update for .net framework 3.5 and 4.7.2 and 4.8 for windows 10 version 1809 for x64	Install this update to resolve issues in Windows. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article for more information. After you install this item, you may have to restart your computer.	KB5003778	1

Note:

1. System not installing the June cycle patches on the PAS-CHAOS-LTSC Win 10 device. Bug raised for same.
Bug Id: 1210252



BD Product Name: **BD Pyxis™ Medstation™ ES**
Date of Critical or Security Patches: June 2021
Abstract: Critical or Security Patches – June 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for June 2021. The patches include updates that may adversely affect BD Pyxis™ MedStation™ ES. These patches will not be and should not be installed on BD Pyxis™ MedStation™ ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

Patch Name	Description	Patch ID	Notes
Cumulative update for windows10 version 1809for x64 based systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5003646	1

Note:

1. System not installing the June cycle patches on the MedSTN-CHAOS-LTSC Win 10 device. Bug raised for same.
Bug Id: 1210252

