

BD Product Name: **BD Pyxis™ Supply**

Date of Critical or Security Patches: December 2020

Abstract: Critical or Security Patches – December 2020

### **Microsoft® & Third-Party Patches**

BD has identified patches from Microsoft® that have been identified as critical or security related for December 2020 will not install on BD Pyxis™ Supply.

The patch includes updates that may adversely affect Pyxis™ Supply. This patch will not be and should not be installed on BD Pyxis™ Supply at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: This patch only apply to the latest supported version(s) of the BD offering.

<b>Patch Name</b>	<b>Description</b>	<b>Patch ID</b>	<b>Notes</b>
2020-12 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems	This security update includes improvements and fixes that were a part of update KB4586807 (released November 10, 2020) and addresses the following issues: Addresses an issue in which PDF24 Creator version 9.1.1 cannot open .txt files. Security updates to Windows Graphics and Windows Peripherals.	<a href="#">KB4592498</a>	1
2020-12 Security Only Quality Update for Windows Server 2008 for x86-based Systems	This security update includes Security updates to Windows Graphics and Windows Peripherals.	<a href="#">KB4592504</a>	1

Note:

1. Due to environment issue patch was unable to be validated this month.

