

BD Product Name: BD Pyxis Medstation ES

Date of Critical or Security Patches: August 2019

Abstract: Compatablitiy issue with KB4512517

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for August 2019.

The patches include updates that may adversely affect BD Pyxis MedStation ES. These patches will not be and should not be installed on BD Pyxis MedStation ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

Patch Name	Patch Description	Patch ID	Notes
4512517	OS build 14393.3144	KB 4512517	1

Notes

1. Update will not install on ES 1.6.0 Chaos images. Patch is not approved for installation.

