

Security Patches: BD Synapsys™

August 2023

BD has identified patches from Microsoft® that have been identified as critical or security related through June 2023. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

The following patches have been validated for Windows Server 2016.

Patch name	Description	Patch ID	Notes
Servicing stack update for Windows Server 2016	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	KB5023788	This update must be installed prior to all other updates
Windows Server 2016 Cumulative security update.	Addresses several security issues for the Windows operating system.	KB5027219	
Windows Server 2016 .NET Core 3.1 security update.	Addresses several security issues for the .NET Core package	KB5021953	
Windows Server 2016 Security update for .NET framework 4.8.	This security update includes cumulative reliability improvements for .NET framework 4.8.	KB5027123	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.





The following patches have been validated for Windows Server 2019

Patch name	Description	Patch ID	Notes
Servicing stack update for Windows Server 2019	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	KB5005112	This update must be installed prior to all other updates
Windows Server 2019 Cumulative security update.	Addresses several security issues for the Windows operating system.	KB5027222	
Windows Server 2016 .NET Core 3.1 security update.	Addresses several security issues for the .NET Core package	KB5021953	
Windows Server 2019 Security update for .NET framework 4.8	This security update includes cumulative reliability improvements for .NET framework 4.8.	KB5027124	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.