

Security Patches: BD Synapsys™

December 2025

BD has identified patches from Microsoft® that have been identified as critical or security related through October 2025. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

The following patches have been validated for Windows Server 2016.

Patch name	Description	Patch ID	Notes
Servicing stack update for Windows Server 2016	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	KB5066584	This update must be installed prior to all other updates
Cumulative Update for .NET Framework 3.5 & 4.8 for Windows Server 2016	Addresses several security issues for the Windows operating system.	KB5066136	
Security update for SQL Server	Addresses several fixes and resolves vulnerabilities for the windows operating system	KB5065227	
Windows Server 2016 Cumulative Security update.	Addresses several security issues for the Windows operating system.	KB5070882	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.





The following patches have been validated for Windows Server 2019

Patch name	Description	Patch ID	Notes
Servicing stack update for Windows Server 2019	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	KB5005112	This update must be installed prior to all other updates
Windows Server 2019 Cumulative security update.	Addresses several security issues for the Windows operating system.	KB5070883	
Security update for SQL Server	Addresses several fixes and resolves vulnerabilities for the windows operating system	KB5065222	
Cumulative Update for .NET Framework 3.5 & 4.8 for Windows Server 2019	Addresses several security issues for the Windows operating system.	KB5066137	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.