

BD Product Name: **Identity provider management (IDM)**

Date of Critical or Security Patches: December 2020

Abstract: Critical or Security Patches – December 2020

### Microsoft® & Third-Party Patches

BD has identified patches from Microsoft that have been identified as critical or security related for December 2020. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Description	Patch ID	Notes
2020-12 Security Only Quality Update for Windows Server 2012 R2 for x64-based Systems (KB4592495)	This security update includes quality improvements. Key changes include: <ul style="list-style-type: none"><li>Security updates to Windows Graphics, Windows Peripherals, and Windows Core Networking.</li></ul>	<a href="#">KB4592495</a>	None
2020-12 Security Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems (KB4592484)	This security update includes improvements and fixes that were a part of update KB4586845 (released November 10, 2020) and addresses the following issues: <ul style="list-style-type: none"><li>Addresses an issue in which PDF24 Creator version 9.1.1 cannot open .txt files.</li><li>Security updates to Windows Graphics, Windows Peripherals, and Windows Core Networking.</li></ul>	<a href="#">KB4592484</a>	None

<p>2020-12 Cumulative Update for Windows Server 2016 for x64-based Systems (KB4593226)</p>	<p>This security update includes quality improvements. Key changes include:</p> <ul style="list-style-type: none"><li>• Addresses an issue that prevents the PDF24 app, version 9.1.1, from opening .txt files.</li><li>• Security updates to the Microsoft Graphics Component, Windows Media, Windows Fundamentals, and Windows Virtualization.</li></ul>	<p><a href="#">KB4593226</a></p>	<p>None</p>
--	--	----------------------------------	-------------