

Security Patches: Identity Management

March 2024

BD has identified patches from Microsoft® that have been identified as critical or security related for March 2024. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity to maintain the correct security posture of the system(s).

Microsoft® patches

Patch name	Description	Patch ID	Notes
2024-03 Cumulative Update for Windows Server 2016 for x64-based Systems (KB5035855)	<p>This security update includes quality improvements. When you install this KB:</p> <p>This update makes miscellaneous security improvements to internal OS functionality. No additional issues were documented for this release.</p>	KB5035855	N/A
2024-03 Servicing Stack Update for Windows Server 2016 for x64-based Systems (KB5035962)	<p>This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. Additionally, reliability improvements to the secure boot servicing stack were made to help with safe roll out of upcoming security updates for Secure Boot.</p>	KB5035962	N/A
Remove specific prevalent malware with Windows Malicious Software Removal Tool (KB890830)	<p>This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. Servicing stack updates (SSU) makes sure that you have a robust and reliable servicing stack so that your devices can receive and install Microsoft updates.</p>	KB890830	N/A

2024-03 Servicing Stack Update for Windows Server 2012 R2 for x64-based Systems (KB5035968)

This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. Additionally, reliability improvements to the secure boot servicing stack were made to help with safe roll out of upcoming security updates for Secure Boot.

[KB5035968](#)

N/A

2024-03 Security Monthly Quality Rollup for Windows Server 2012 R2 for x64-based Systems (KB5035885)

This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. Servicing stack updates (SSU) makes sure that you have a robust and reliable servicing stack so that your devices can receive and install Microsoft updates

[KB5035885](#)

N/A

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

