

Security Patches: Identity Management

November 2023

BD has identified patches from Microsoft® that have been identified as critical or security related for November 2023. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity to maintain the correct security posture of the system(s).

Microsoft® patches

Patch name	Description	Patch ID	Notes
Remove specific prevalent malware with Windows Malicious Software Removal Tool (KB890830)	Microsoft releases the MSRT on a monthly cadence as part of Windows Update or as a standalone tool. Use this tool to find and remove specific prevalent threats and reverse the changes they have made (see covered malware families). For comprehensive malware detection and removal, consider using Windows Defender Offline or Microsoft Safety Scanner.	KB890830	N/A
Servicing stack update for Windows Server 2016: November 14, 2023(KB5032391)	This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. Servicing stack updates (SSU) makes sure that you have a robust and reliable servicing stack so that your devices can receive and install Microsoft updates.	KB5032391	N/A
OS Build 14393.6452(KB5032197)	This update addresses security issues for Windows operating system. Microsoft is not currently aware of any issues with this update.	KB5032197	N/A

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

