

## Products in Scope

BD has provided the list below in order to better help our customers identify any BD products with a computer processing unit (CPU) that has the potential to be vulnerable to the Meltdown and Spectre hard-ware based vulnerabilities. The list below of BD products is currently dynamic and will be updated as we complete analysis of products in scope\*.

BD uses the [Common Vulnerability Scoring System v3.0](#) (CVSS) to rate the characteristics of a vulnerability and produce a base numerical score reflecting the severity.

## Meltdown

The following BD products were determined to have a [CVSS rating of 2.5](#) during our internal evaluation.

- Accuri C6 Plus
- FACS Aria Fusion
- FACS Aria I/II/III
- FACSCalibur
- FACSCanto 10-color
- FACSCanto 10-color clinical
- FACSCanto II
- FACSCanto II clinical
- FACSCelesta
- FACSCount
- FACSDuet Sample Prep (ASaP)
- FACSJazz
- FACSLink Interface
- FACS Lyric
- FACSMelody
- FACSSample Prep Assistant III
- FACSVerse
- FACS Via
- GenCell CliC
- Influx
- LSR II
- LSRFortessa
- LSRFortessa X-20
- Lyse Wash Assistant
- Panel Designer
- BD Assurity Linc™
- BD BACTEC™ 9120/9240
- BD BACTEC™ FX
- BD BACTEC™ FX40
- BD MAX
- BD Data Innovations™
- BD EpiCenter™
- BD Kiestra™ TLA/WCA
- BD Phoenix™ 100
- BD Phoenix™ M50
- BD Totalys™ SlidePrep
- BD Viper LT™
- BD Viper XTR™
- Pyxis® Anesthesia ES
- Pyxis® Anesthesia System 3500
- Pyxis® CathRack
- Pyxis® CIISafe - Workstation
- Pyxis® CUBIE Replenishment Station
- Pyxis® DuoStation
- Pyxis® EcoStation System
- Pyxis® MedStation 3500
- Pyxis® MedStation 4000
- Pyxis® MedStation ES
- Pyxis® ParAssist System
- Pyxis® Parx
- Pyxis® Parx handheld
- Pyxis® ProcedureStation
- Pyxis® ScrubStation System
- Pyxis® StockStation System
- Pyxis® Supply Roller
- Pyxis® SupplyStation
- Pyxis® Infant Care Verification
- Pyxis® Medication Administration
- Pyxis® Nursing Data Collection
- Pyxis® Specimen Collection Verification
- Pyxis® Transfusion Verification
- Rowa Dose
- Rowa Smart
- Rowa Vmax System

## Spectre

The following product was determined to have a [CVSS Rating of 1.9](#) during our internal evaluation.

- BD Phoenix™ AP

The following products were determined to have a [CVSS rating of 2.5](#) during our internal evaluation.

- BD Phoenix™ 100
- BD Phoenix™ M50

The following products were determined to have a [CVSS rating of 3.0](#) during our internal evaluation.

- Accuri C6 Plus
- FACS Aria Fusion
- FACS Aria I/II/III
- FACSCalibur
- FACSCanto 10-color
- FACSCanto 10-color clinical
- FACSCanto II
- FACSCanto II clinical
- FACSCelesta
- FACSCount
- FACSDuet Sample Prep (ASaP)
- FACSJazz
- FACSMelody
- FACSPresto
- FACSSample Prep Assistant (SPA) III
- FACSVerse
- FACSVia
- GenCell CliC
- Influx
- LSR II
- LSRFortessa
- LSRFortessa X-20
- Lyse Wash Assistant
- BD Assurity Linc™
- BD BACTEC™ 9050
- BD BACTEC™ 9120/9240
- BD BACTEC™ FX
- BD BACTEC™ FX40
- BD Data Innovations™
- Focal Point - Linux
- Focal Point - Solaris
- Innova
- Kiestra InoquLA Standalone
- BD Kiestra™ TLA/WCA
- BD PrepStain
- Totalys Multiprocessor
- BD Totalys™ SlidePrep
- BD Veritor™ Plus System
- BD Viper LT™
- BD Viper XTR™

The following products were determined to have a [CVSS rating of 3.9](#) during our internal evaluation.

- Alaris Systems Manager
- FACSLink Interface
- FACSLyric
- Panel Designer
- Pyxis Anesthesia ES
- Pyxis Anesthesia System 3500
- Pyxis® CathRack v8
- Pyxis® CIISafe - Workstation
- Pyxis® CUBIE Replenishment Station
- Pyxis® DuoStation
- Pyxis® EcoStation System
- Pyxis® MedStation 3500
- Pyxis® MedStation 4000 Console
- Pyxis® MedStation ES
- Pyxis® ParAssist System
- Pyxis® Parx
- Pyxis® Parx handheld
- Pyxis® ProcedureStation
- Pyxis® ScrubStation System
- Pyxis® StockStation System
- Pyxis® Supply Roller
- Pyxis® SupplyStation
- Pyxis® Infant Care Verification



- Pyxis® Medication Administration
- Pyxis® Nursing Data Collection
- Pyxis® Specimen Collection Verification
- Pyxis® Transfusion Verification
- Rowa Dose
- Rowa Smart
- Rowa Vmax System

Note: This list provided does not indicate the patch or device status. The intended use of these products does not include email and/or internet browsing. BD is determining compensating controls and prioritizing patch validation and impact for all hosted solutions, including MedMinded and Knowledge Portal.

\*Some products experience processor specification changes during the manufacturing lifecycle. While specific versions of a product may not be affected BD recommends applying BD mitigations or compensating controls across all affected product lines.