

Catheters to aid in clot management

Fresh, soft emboli/thrombi

Fogarty™ Arterial Embolectomy Catheter

Removal of fresh, soft emboli and thrombi from vessels in the peripheral arterial system.

- Hand-tied using a recessed winding technique to secure balloon to catheter
- Rounded tip promotes easy insertion for reduced trauma

Fogarty™ Fortis Arterial Embolectomy Catheter

Removal of fresh, soft emboli and thrombi from vessels of the peripheral arterial system.

- Can withstand greater pull-force than the Fogarty™ Arterial Embolectomy Catheter
- The addition of barium to the catheter tubing provides radiopacity

Fresh, soft emboli/thrombi, temporary occlusion, infusion

Fogarty™ Thru-Lumen Embolectomy Catheter

Multi-use guidewire-compatible design enables removal of fresh, soft emboli and thrombi from vessels in the arterial system as well as temporary occlusion of blood vessels, infusion of fluids, and blood sampling

- Fogarty™ Balloon conforms to the vessel to maximize clot removal while minimizing trauma
- Provides the solution for distal thrombus removal in the arterial system
- The addition of barium to the catheter shaft material and stainless steel markers proximal and distal to the balloon provide radiopacity*

* 3F Fogarty™ Thru-Lumen Embolectomy Catheter body does not include barium

Emboli/thrombi

Fogarty Corkscrew™ Catheter

Removal of emboli and thrombi from either native vessels or synthetic grafts in the arterial system




- Corkscrew design increases the surface area for entrapping fibrous material
- Radiopaque and flexible

Organized thrombus/ adherent non-thrombotic material

Fogarty™ Graft Thrombectomy Catheter

Designed for removing organized thrombus and adherent non-thrombotic material from synthetic bypass grafts, specifically, AV fistulas for dialysis access and synthetic grafts for iliofemoral bypass

- Provides a larger diameter range than the Fogarty Corkscrew™ Catheter
- The pulling force can be varied by adjusting the diameter



To learn more, visit bd.com/fogarty
or call 1-800-FOGARTY.

CAUTION: Federal (United States) law restricts these devices to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

BD, Irvine, CA 92618 USA

bd.com

BD, the BD Logo, Fogarty, Fogarty Corkscrew, Fogarty Fortis and 1-800-FOGARTY are trademarks of Becton, Dickinson and company or its affiliates. All other trademarks are the property of their respective owners. © 2026 BD. All rights reserved.
BD-168754 v1.0 US

