# BD Eswab® liquid based collection and transport system for microbiology swab samples

BD Eswab® system elutes over 90% of patient specimen into the Liquid Amies transport medium. Multiple investigations can be performed from the same sample:



Multiple culture plates



Rapid antigen tests\*



Automation



Molecular testing\*





Gram stains



Flocked Swab

# Why Flocked?







1 mL Liquid Amies

# BD Eswab® helps busy laboratories do more with less



### Better diagnostics

Flocked swabs paired with a liquid based system give way to quantitative, measurable and consistent sample transfer. Evidence-based research shows that samples collected using Flocked swabs improve test sensitivity, eluting more than 90% of the specimen.




#### Reduce costs

A broad range of testing applications eliminates costs associated with stocking and management of numerous collection devices.



#### Patient comfort

Clinicians report better patient comfort due to ergonomic, anatomic swab design and softer texture. Additionally, because one sample can be used for multiple tests, fewer samples are collected from the patient.



# Automation ready

Liquid Based Microbiology systems are easily processed on automated specimen processors and liquid handling pipetting systems, minimizing manual handling.



Eswab® combines a flocked swab with 1 mL of Liquid Amies in a plastic self standing tube. The system elutes over 90% of patient specimen into the liquid medium.



BD Life Sciences, 7 Loveton Circle, Sparks, MD 21152-0999, U.S. 800.638.8663

#### bd.com



## Change Management



Ready to make the switch to better microbiology, improved patient care and laboratory cost savings? BD provides hands-on expertise to facilitate new product implementation. We are available to guide you with with training, verification guidance, and more!

