



BD Infusion Set Catalog

2024

Master Catalog (U.S.)

Extension, Gravity, Secondary and Pump Infusion Sets



BD is a global leader in infusion therapy across the healthcare continuum, developing, manufacturing and delivering innovative and clinically relevant products and services to the healthcare industry.

BD offers a broad portfolio of infusion sets, extension sets and accessories to meet your specific clinical needs.

The configuration and components used in an infusion set will depend on the specific needs of the patient and the type of fluid or medication being administered. Healthcare providers should carefully consider the options available and choose the most appropriate configuration and components for each patient and situation to help ensure safe and effective administration of infusion therapy.

Visit [bd.com](https://www.bd.com) for complete catalog

Sales and Customer Service

 **1.844.8BD LIFE (1.844.823.5433)**

Clinical Support

 **1.800.555.7422 Ext 5**

 **medical.information@bd.com**



Contents

Click on the titles below to view respective pages

Infusion Practices	2
Extension Sets	8
BD MaxZero™ Extension Sets	9
BD MaxZero™ Extension Sets	10
BD MaxZero™ Extension Sets Pressure Rated	10
BD MaxZero™ T-connector Extension Sets	12
BD MaxZero™ T-connector Extension Sets Pressure Rated	13
BD MaxZero™ Y-site Extension Sets Pressure Rated	14
BD MaxZero™ Bifuse Extension Sets	15
BD MaxZero™ Bifuse Extension Sets Pressure Rated	15
BD MaxZero™ Bifuse T-connector Extension Sets	16
BD MaxZero™ Bifuse Extension Sets with Back Check Valves	17
BD MaxZero™ Trifuse Extension Sets	17
BD MaxZero™ Trifuse T-connector Extension Sets	18
BD MaxZero™ Trifuse Extension Sets with Back Check Valves	18
BD MaxZero™ Quadfuse Extension Sets	19
BD MaxZero™ Quadfuse Extension Sets with Back Check Valves	19
BD MaxZero™ In-Line Filter Extension Sets	20
BD MaxZero™ Bifuse Filter Extension Sets	20
BD MaxZero™ Bifuse Filter Extension Sets with Back Check Valves	21
BD MaxZero™ Trifuse Filter Extension Sets	21
BD MaxZero™ Quadfuse Filter Extension Sets	23
BD MaxPlus™ Extension Sets	24
BD MaxPlus™ Extension Sets	25
BD MaxPlus™ Pressure Rated Extension Sets	25
BD MaxPlus™ T-connector Extension Sets	27
BD MaxPlus™ Pressure Rated T-connector Extension Sets	27
BD MaxPlus™ Pressure Rated Y-site Extension Sets	28
BD MaxPlus™ Y-site Extension Sets with BCVs	29
BD MaxPlus™ Bifuse Extension Sets	29

Contents (continued)

Click on the titles below to view respective pages

BD MaxPlus™ Bifuse Extension Sets Pressure Rated	29
BD MaxPlus™ Trifuse Extension Sets	30
BD Q-Syte™ Extension Sets	31
BD Q-Syte™ Extension Sets	32
BD Q-Syte™ Y-site Extension Sets	33
BD Q-Syte™ Bifuse Extension Sets	33
BD Q-Syte™ Trifuse Extension Sets	33
BD SmartSite™ Extension Sets	35
BD SmartSite™ Extension Sets	36
BD SmartSite™ Pressure Rated Extension Sets	37
BD SmartSite™ T- connector Extension Sets	38
BD SmartSite™ Extension Sets with Texium™ Male Luer	38
BD SmartSite™ Bifuse Extension Sets	38
BD SmartSite™ Trifuse Extension Sets	39
BD SmartSite™ Trifuse Extension Sets with Back Check Valves	40
BD SmartSite™ Multiport Extension Sets	40
BD SmartSite™ Filter Extension Sets	41
Extension Sets	42
Extension Sets	43
Pressure Rated Extension Sets	46
T-connector Extension Sets	47
T-connector Extension Sets Pressure Rated	47
J-loop Extension Sets	48
Extension Sets with Anti-Siphon Valve	49
Extension Sets with Anti-Siphon Valves	50
Extension Sets with Y-site Access Ports	51
Extension Sets with Y-sites	52
Extension Sets with Y-sites Pressure Rated	54
Extension Sets with Y-sites and Back Check Valves	54
Extension Sets with Sideport	55

Contents (continued)

Click on the titles below to view respective pages

Bifuse Extension Sets	55
Trifuse Extension Sets	56
Quadfuse Extension Sets	56
Specialty Tubing Extension Sets	57
PE Lined PVC Extension Sets	58
Filter Extension Sets	59
0.2 Micron Filter Extension Sets	60
0.2 Micron Filter PE Lined PVC Extension Sets	61
0.2 Micron Filter PE Lined PVC Extension Sets with Texium™ Male Luer	61
1.2 Micron Filter Extension Sets	61
1.2 Micron Filter PE Lined PVC Extension Sets	62
Stopcock Extension Sets	63
Stopcock Extension Sets	64
Stopcock Extension Sets with Y-sites	64
In-line Flow Controller Extension Sets	67
Flow Controller Single Scale with Y-sites	68
Gravity Sets	69
Primary Gravity Sets	71
BD Primary Gravity Sets	71
BD Primary Gravity Set Y Access Basic	71
BD Primary Gravity Sets with In-Line Filter	75
BD Primary Gravity Flow Controller Sets	75
BD Primary Gravity Burette (<i>Add On</i>) Sets	76
BD Primary Gravity Stopcock Sets	77
BD Primary Gravity Manifold Sets	80
Secondary Administration Sets	83
BD Secondary Gravity Sets	84
BD Secondary Texium™ Sets	85
BD Secondary Transfer Set	86
BD Secondary Blood Infusion Set	86

Contents (continued)

Click on the titles below to view respective pages

BD Alaris™ Infusion System - Infusion Sets	87
BD Alaris™ Pump Infusion Sets	88
BD Alaris™ Pump Infusion Sets	89
BD Alaris™ Pump 0.2µm Filter Infusion Sets	92
BD Alaris™ Pump PE Lined PVC 0.2µm Filter Infusion Sets	92
BD Alaris™ Pump 1.2µm Filter Infusion Sets	93
BD Alaris™ Pump PE Lined PVC 1.2µm Filter Infusion Sets	93
BD Alaris™ Pump Burette Infusion Sets	94
BD Alaris™ Pump PE Lined PVC Burette 0.2µm Filter Infusion Sets	94
BD Alaris™ Pump Stopcock/Manifold Infusion Sets	94
BD Alaris™ Pump Blood Infusion Sets	96
BD Alaris™ Pump Blood Dual Spike Infusion Sets	96
BD Alaris™ Pump Infusion Sets with SmartSite™ Bag Access Port and Texium™ Closed Male Luer	97
BD Alaris™ Pump Infusion Sets PE Lined PVC with SmartSite™ Bag Access Port and Texium™ Closed Male Luer	97
BD Alaris™ Pump Infusion Sets with Texium™ Closed Male Luer	98
BD Alaris™ Pump Infusion Sets PE Lined PVC with Texium™ Closed Male Luer	98
BD Alaris™ Pump Infusion Sets with SmartSite™ Bag Access Port	98
BD Alaris™ Pump Infusion Sets PE Lined PVC with SmartSite™ Bag Access Port	99
BD Alaris™ Pump Infusion Sets PE Lined PVC Filter Infusion Sets with SmartSite™ Bag Access Port and Texium™ Closed Male Luer	99
BD Alaris™ Pump Infusion Sets PE Lined PVC Filter Infusion Sets with SmartSite™ Bag Access	

Contents

Click on the titles below to view respective pages

Port	100
BD Alaris™ Pump Syringe Adapter Infusion Sets	100
BD Alaris™ Pump Infusion Half Sets	100
BD Alaris™ Syringe Infusion Sets	101
BD Alaris™ Syringe Infusion Sets with Pressure Sensing Disc	102
BD Alaris™ Syringe PE Lined PVC Infusion Sets with Pressure Sensing Disc	102
BD Alaris™ Syringe Filter Infusion Sets with Pressure Sensing Disc	102
BD Alaris™ Syringe Infusion Sets	103
BD Alaris™ Syringe Infusion Set	103
BD Alaris™ PCA Infusion Sets	104
Alaris™ PCA Infusion Set Kit	105
Alaris™ PCA Infusion Set with Open Y-site	105
Glossary	106
Symbols	109
Components details	111

Extension Sets


Sales and Customer Service

 **1-844-8BD LIFE (1-844-823-5433)**

Clinical Support

 **1-800-555-7422 Ext 5**


 **medical.information@bd.com**



BD MaxZero™ Extension Sets




BD MaxZero™ Extension Sets

Microbore Extension Set				
Part no. MZ9277	IVA			 0.5mm
Length 15 cm (6 in) Priming vol. 0.2 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxZero™ Extension Set Microbore 1 MaxZero™ Bonded 1 SC RLL	US	CE	STERILE
Part no. MZ9288	IVA			 0.5mm
Length 15 cm (6 in) Priming vol. 0.4 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxZero™ Extension Set Microbore 1 MaxZero™ Removable 1 SC RLL	US	CE	STERILE
Part no. MZ9267	IVA			 0.5mm
Length 152 cm (60 in) Priming vol. 0.6 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxZero™ Extension Set Microbore 1 MaxZero™ Bonded 1 SC RLL	US	CE	STERILE

BD MaxZero™ Extension Sets Pressure Rated

Smallbore/Minibore Pressure Rated (325 psi) Extension Set				
Part no. MZ5306	IVA			 1mm
Length 18 cm (7 in) Priming vol. 0.4 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxZero™ Extension Set Minibore 1 MaxZero™ Bonded Pressure Rated 325 PSI RLL	US	CE	STERILE



Part no. MZ5301		IVA			 1mm
Length	18 cm (7 in)				
Priming vol.	0.4 mL				
Ctn. qty.	50	Product description BD MaxZero™ Extension Set Minibore 1 MaxZero™ Bonded Pressure Rated 325 PSI 1 SC RLL		US CE	STERILE
Material	PVC*				

Part no. MZ5302		IVA			 1mm
Length	18 cm (7 in)				
Priming vol.	0.4 mL				
Ctn. qty.	50	Product description BD MaxZero™ Extension Set Minibore 1 MaxZero™ Removable Pressure Rated 325 PSI 1 SC RLL		US CE	STERILE
Material	PVC*				


Standardbore Pressure Rated (325 psi) Extension Set

Part no. MZ5305		IVA			 1.5mm
Length	18 cm (7 in)				
Priming vol.	0.5 mL				
Ctn. qty.	50	Product description BD MaxZero™ Extension Set Standardbore 1 MaxZero™ Bonded Pressure Rated 325 PSI RLL		US CE	STERILE
Material	PVC*				

Part no. MZ5309		IVA			 1.5mm
Length	18 cm (7 in)				
Priming vol.	0.5 mL				
Ctn. qty.	50	Product description BD MaxZero™ Extension Set Standardbore 1 MaxZero™ Bonded Pressure Rated 325 PSI 1 SC RLL		US CE	STERILE
Material	PVC*				


Part no. MZ5303		IVA			 1.5mm
Length	18 cm (7 in)				
Priming vol.	0.5 mL				
Ctn. qty.	50	Product description BD MaxZero™ Extension Set Standardbore 1 MaxZero™ Bonded Pressure Rated 325 PSI 1 PC RLL		US CE	STERILE
Material	PVC*				




Part no. MZ5310	IVA	 1.5mm		
Length 18 cm (7 in) Priming vol. 0.5 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxZero™ Extension Set Standardbore 1 MaxZero™ Removable Pressure Rated 325 PSI 1 SC RLL		US	CE	STERILE



Part no. MZ5304	IVA	 1.5mm		
Length 18 cm (7 in) Priming vol. 0.5 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxZero™ Extension Set Standardbore 1 MaxZero™ Removable Pressure Rated 325 PSI 1 PC RLL		US	CE	STERILE

BD MaxZero™ T-connector Extension Sets



Microbore T-connector Extension Set				
Part no. MZT1004	IVA	 0.5mm		
Length 15 cm (6 in) Priming vol. 0.3 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxZero™ Extension Set Microbore 1 MaxZero™ Bonded 1 SC INJ/RLL T-connector		US	CE	STERILE

Part no. MZXT9004	IVA	 0.5mm		
Length 15 cm (6 in) Priming vol. 0.4 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxZero™ Extension Set Microbore 1 MaxZero™ Bonded 1 SC Nac/LS T-connector		US	CE	STERILE




Part no. MZXT9001		IVA	 1mm
Length 15 cm (6 in)			
Priming vol. 0.2 mL			
Ctn. qty. 50			
Material PVC*			
Product description BD MaxZero™ Extension Set Microbore 1 MaxZero™ Bonded w FL ADPT 1 SC Nac/RLL T-connector		US CE	STERILE
Part no. MZT1003		IVA	 0.5mm
Length 15 cm (6 in)			
Priming vol. 0.4 mL			
Ctn. qty. 50			
Material PVC*			
Product description BD MaxZero™ Extension Set Microbore 1 MaxZero™ Removable 1 SC INJ/RLL T-connector		US CE	STERILE




BD MaxZero™ T-connector Extension Sets Pressure Rated

Smallbore/Minibore Pressure Rated (325 psi) T-connector Extension Set			
Part no. MZXT5307		IVA	 1mm
Length 15 cm (6 in)			
Priming vol. 0.6 mL			
Ctn. qty. 50			
Material PVC*			
Product description BD MaxZero™ Extension Set Minibore 1 MaxZero™ Removable Pressure Rated 325 PSI 1 SC Nac/RLL T-connector		US CE	STERILE
Standardbore Pressure Rated (325 psi) T-connector Extension Set			
Part no. MZXT5304		IVA	 1.5mm
Length 15 cm (6 in)			
Priming vol. 0.7 mL			
Ctn. qty. 50			
Material PVC*			
Product description BD MaxZero™ Extension Set Standardbore 1 MaxZero™ Bonded Pressure Rated 325 PSI 1 SC Nac/RLL T-connector		US CE	STERILE



Part no. MZXT5306	IVA	 1.5mm
Length 16 cm (6 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*		
Product description BD MaxZero™ Extension Set Standardbore 1 MaxZero™ Removable Pressure Rated 325 PSI 1 SC Nac/RLL T-connector		US CE STERILE

BD MaxZero™ Y-site Extension Sets Pressure Rated

Smallbore/Minibore Pressure Rated (325 psi) Y-site Extension Set		
Part no. MZX5301	IVA	 1mm
Length 20 cm (8 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*		
Product description BD MaxZero™ Extension Set Minibore 1 MaxZero™ Bonded Pressure Rated 325 PSI 1 SC 1 NAC-Y RLL		US CE STERILE
Part no. MZX5302	IVA	 1mm
Length 20 cm (8 in) Priming vol. 1 mL Ctn. qty. 50 Material PVC*		
Product description BD MaxZero™ Extension Set Minibore 1 MaxZero™ Removable Pressure Rated 325 PSI 1 SC 1 NAC-Y RLL		US CE STERILE
Standardbore Pressure Rated (325 psi) Y-site Extension Set		
Part no. MZX5305	IVA	 1.5mm
Length 20 cm (8 in) Priming vol. 1 mL Ctn. qty. 50 Material PVC*		
Product description BD MaxZero™ Extension Set Standardbore 1 MaxZero™ Bonded Pressure Rated 325 PSI 1 SC 1 NAC-Y RLL		US CE STERILE



Part no. MZ5306	IVA	
Length 20 cm (8 in) Priming vol. 1 mL Ctn. qty. 50 Material PVC*		
Product description BD MaxZero™ Extension Set Standardbore 1 MaxZero™ Removable Pressure Rated 325 PSI 1 SC 1 NAC-Y RLL		US CE STERILE

BD MaxZero™ Bifuse Extension Sets

Microbore Bifuse Extension Set		
Part no. MZ9265	IVA	
Length 15 cm (6 in) Priming vol. 0.4 mL Ctn. qty. 50 Material PVC*		
Product description BD MaxZero™ Bifuse Extension Set Microbore 2 MaxZero™ Bonded 2 SC RLL		US CE STERILE
Part no. MZ9290	IVA	
Length 15 cm (6 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*		
Product description BD MaxZero™ Bifuse Extension Set Microbore 2 MaxZero™ Removable 2 SC RLL		US CE STERILE

BD MaxZero™ Bifuse Extension Sets Pressure Rated

Smallbore/Minibore Bifuse Pressure Rated (325 psi) Extension Set		
Part no. MZ5307	IVA	
Length 18 cm (7 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*		
Product description BD MaxZero™ Bifuse Extension Set Minibore 2 MaxZero™ Removable Pressure Rated 325 PSI 2 SC RLL		US CE STERILE






Standardbore Bifuse Pressure Rated (325 psi) Extension Set				
Part no. MZ5312	IVA			 1.5mm
Length 18 cm (7 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxZero™ Bifuse Extension Set Standardbore 2 MaxZero™ Bonded Pressure Rated 325 PSI RLL			US CE	STERILE
Part no. MZ5316	IVA			 1.5mm
Length 18 cm (7 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxZero™ Bifuse Extension Set Standardbore 2 MaxZero™ Bonded Pressure Rated 325 PSI 2 SC RLL			US CE	STERILE
Part no. MZ5313	IVA			 1.5mm
Length 18 cm (8 in) Priming vol. 1 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxZero™ Bifuse Extension Set Standardbore 2 MaxZero™ Removable Pressure Rated 325 PSI 2 SC RLL			US CE	STERILE

BD MaxZero™ Bifuse T-connector Extension Sets


Smallbore/Minibore Bifuse T-connector Pressure Rated (325 psi) Extension Set				
Part no. MZXT5303	IVA			 1mm
Length 15 cm (6 in) Priming vol. 1 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxZero™ Bifuse Extension Set PVC Smallbore 2 MaxZero™ Bonded Pressure Rated 325 PSI 3 SC Nac/RLL T-connector			US CE	STERILE



BD MaxZero™ Bifuse Extension Sets with Back Check Valves

Microbore Bifuse Extension Set with BCV			
Part no. MZ9295	IVA		 0.5mm
Length 18 cm (7 in) Priming vol. 0.5 mL Ctn. qty. 50 Material PVC*			
	Product description BD MaxZero™ Bifuse Extension Set Microbore 2 MaxZero™ Bonded 1 BCV 2 SC RLL	US CE	STERILE
Microbore Bifuse Extension Set with BCV			
Part no. MZ9292	IVA		 0.5mm
Length 18 cm (7 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*			
	Product description BD MaxZero™ Bifuse Extension Set Microbore 2 MaxZero™ Removable 1 BCV 2 SC RLL	US CE	STERILE
Smallbore/Minibore Bifuse Extension Set			
Part no. MZ9279	IVA		 1mm
Length 18 cm (7 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*			
	Product description BD MaxZero™ Bifuse Extension Set Smallbore 2 MaxZero™ Bonded 2 BCV 2 SC RLL	US CE	STERILE

BD MaxZero™ Trifuse Extension Sets

Microbore Trifuse Extension Set			
Part no. MZ9266	IVA		 0.5mm
Length 10 cm (4 in) Priming vol. 0.6 mL Ctn. qty. 50 Material PVC*			
	Product description BD MaxZero™ Trifuse Extension Set Microbore 3 MaxZero™ Bonded 3 SC RLL	US CE	STERILE



Microbore Trifuse Extension Set			
Part no. MZ9327	IVA		 0.5mm
Length 18 cm (7 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Trifuse Extension Set Microbore 2 Open Female Luer 1 MaxZero™ Bonded 3 SC RLL		US CE	STERILE

Smallbore/Minibore Trifuse Extension Set			
Part no. MZ9326	IVA		 1mm
Length 18 cm (7 in) Priming vol. 0.6 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Trifuse Extension Set Minibore 2 Open Female Luer 1 MaxZero™ Bonded 3 SC RLL		US CE	STERILE

BD MaxZero™ Trifuse T-connector Extension Sets

Microbore Trifuse T-connector Extension Set			
Part no. MZXT9003	IVA		 0.5mm
Length 20 cm (8 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Trifuse Extension Set Microbore 3 MaxZero™ Bonded 4 SC Nac/RLL T-connector		US CE	STERILE

BD MaxZero™ Trifuse Extension Sets with Back Check Valves

Microbore Trifuse Extension Set with BCV			
Part no. MZ9278	IVA		 0.5mm
Length 15 cm (6 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Trifuse Extension Set Microbore 3 MaxZero™ Bonded 3 BCV 3 SC RLL		US CE	STERILE



Part no. MZ9313	IVA		0.5mm
Length 20 cm (8 in) Priming vol. 0.9 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Trifuse Extension Set Microbore 3 MaxZero™ Removable 3 SC RLL		US CE	STERILE

Smallbore/Minibore Trifuse Extension Set with 3 bonded MaxZero™ Connectors			
Part no. MZ9281	IVA		1mm
Length 15 cm (5 in) Priming vol. 1 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Trifuse Extension Set Smallbore 3 MaxZero™ Bonded 3 BCV 3 SC RLL		US CE	STERILE

BD MaxZero™ Quadfuse Extension Sets

Microbore Quadfuse Extension Set			
Part no. MZ9275	IVA		0.5mm
Length 15 cm (6 in) Priming vol. 0.9 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Quadfuse Extension Set Microbore 4 MaxZero™ Bonded 4 SC RLL		US CE	STERILE

BD MaxZero™ Quadfuse Extension Sets with Back Check Valves


Smallbore/Minibore Quadfuse Extension Set with BCV			
Part no. MZ9283	IVA		1mm
Length 18 cm (7 in) Priming vol. 1.5 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Quadfuse Extension Set Smallbore 4 MaxZero™ Bonded 4 BCV 4 SC RLL		US CE	STERILE




BD MaxZero™ In-Line Filter Extension Sets

Microbore Filter Extension Set				
Part no. MZ9289	IVA		 0.2 µm	 0.5mm
Length 15 cm (6 in)				
Priming vol. 0.6 mL				
Ctn. qty. 50				
Material PVC*	Product description BD MaxZero™ Filter Extension Set Microbore 1 MaxZero™ Removable 0.2µm Small Volume PES Filter 1 SC RLL	US	CE	STERILE
Part no. MZ9226	IVA		 0.2 µm	 0.5mm
Length 152 cm (60 in)				
Priming vol. 0.8 mL				
Ctn. qty. 50				
Material PVC*	Product description BD MaxZero™ Filter Extension Set Microbore 1 MaxZero™ Removable 0.2µm Small Volume PES Filter 1 SC RLL	US	CE	STERILE
Smallbore/Minibore Filter Extension Set				
Part no. MZ9329	IVA		 1.2 µm	 1mm
Length 198 cm (78 in)				
Priming vol. 2.2 mL				
Ctn. qty. 50				
Material PVC*	Product description BD MaxZero™ Filter Extension Set Minibore 1 MaxZero™ Bonded 1.2µm Small Volume PES Filter 1 SC RLL	US	CE	STERILE



BD MaxZero™ Bifuse Filter Extension Sets

Microbore Bifuse Filter Extension Set				
Part no. MZ9271	IVA		 0.2 µm	 0.5mm
Length 18 cm (7 in)				
Priming vol. 0.7 mL				
Ctn. qty. 50				
Material PVC*	Product description BD MaxZero™ Bifuse Filter Extension Set Microbore 2 MaxZero™ Bonded 0.2µm Small Volume PES Filter 2 SC RLL	US	CE	STERILE




Part no. MZ9294	IVA	 0.2 µm	 0.5mm
Length 18 cm (7 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Bifuse Filter Extension Set Microbore 2 MaxZero™ Removable 0.2um Small Volume PES Filter 2 SC RLL		US CE	STERILE





BD MaxZero™ Bifuse Filter Extension Sets with Back Check Valves

Microbore Bifuse Filter Extension Set with BCV			
Part no. MZ9296	IVA	 0.2 µm	 0.5mm
Length 18 cm (7 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Bifuse Filter Extension Set Microbore 2 MaxZero™ Bonded 1 BCV 0.2um Small Volume PES Filter 2 SC RLL		US CE	STERILE
Part no. MZ9293	IVA	 0.2 µm	 0.5mm
Length 19 cm (7.5 in) Priming vol. 0.9 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Bifuse Filter Extension Set Microbore 2 MaxZero™ Removable 1 BCV 0.2um Small Volume PES Filter 2 SC RLL		US CE	STERILE


BD MaxZero™ Trifuse Filter Extension Sets

Microbore Trifuse Filter Extension Set			
Part no. MZ9270	IVA	 0.2 µm	 0.5mm
Length 14 cm (5.5 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*			
Product description BD MaxZero™ Trifuse Filter Extension Set Microbore 3 MaxZero™ Bonded 0.2um Small Volume PES Filter 3 SC RLL		US CE	STERILE





Part no. MZ9274	IVA			 0.5mm
Length 15 cm (6 in) Priming vol. 1.1 mL Ctn. qty. 50 Material PVC*			US CE	STERILE
Product description BD MaxZero™ Trifuse Filter Extension Set Microbore 3 MaxZero™ Bonded 2 x 0.2um Small Volume PES Filter 3 SC RLL				
Part no. MZ9272	IVA			 0.5mm
Length 152 cm (60 in) Priming vol. 1.2 mL Ctn. qty. 50 Material PVC*			US CE	STERILE
Product description BD MaxZero™ Trifuse Filter Extension Set Microbore 3 MaxZero™ Bonded 0.2um Small Volume PES Filter 3 SC RLL				
Part no. MZ9273	IVA			 0.5mm
Length 152 cm (60 in) Priming vol. 1.5 mL Ctn. qty. 50 Material PVC*			US CE	STERILE
Product description BD MaxZero™ Trifuse Filter Extension Set Microbore 3 MaxZero™ Bonded 0.2um Small Volume PES Filter 3 SC RLL				
Part no. MZ9276	IVA			 0.5mm
Length 152 cm (60 in) Priming vol. 1.2 mL Ctn. qty. 50 Material PVC*			US CE	STERILE
Product description BD MaxZero™ Trifuse Filter Extension Set Microbore 3 MaxZero™ Bonded 1.2um Small Volume PES Filter 3 SC RLL				





Smallbore/Minibore Trifuse Filter Extension Set			
Part no. MZ9325	IVA	  	
Length 18 cm (7 in) Priming vol. 1.1 mL Ctn. qty. 50 Material PVC*			
	Product description BD MaxZero™ Trifuse Filter Extension Set Minibore 2 Open Female Luer 1 MaxZero™ Bonded 1 x 0.2µm & 1 x 1.2µm Small Volume PES Filter 3 SC RLL	US CE	STERILE

BD MaxZero™ Quadfuse Filter Extension Sets

Microbore Quadfuse Filter Extension Set			
Part no. MZ9300	IVA	 	
Length 18 cm (7 in) Priming vol. 1.1 mL Ctn. qty. 50 Material PVC*			
	Product description BD MaxZero™ Quadfuse Filter Extension Set Microbore 4 MaxZero™ Bonded 0.2µm Small Volume PES Filter 4 SC RLL	US CE	STERILE

Smallbore/Minibore Quadfuse Filter Extension Set			
Part no. MZ9324	IVA	  	
Length 18 cm (7 in) Priming vol. 1.6 mL Ctn. qty. 50 Material PVC*			
	Product description BD MaxZero™ Quadfuse Filter Extension Set Minibore 2 Open Female Luer 2 MaxZero™ Bonded 1 x 0.2µm & 1 x 1.2µm Small Volume PES Filter 4 SC RLL	US CE	STERILE

Part no. MZ9322	IVA	 	
Length 18 cm (7 in) Priming vol. 1.2 mL Ctn. qty. 50 Material PVC*			
	Product description BD MaxZero™ Quadfuse Filter Extension Set Microbore 3 Open Female Luer 1 MaxZero™ Bonded 2 BCV 0.2µm Small Volume PES Filter 4 SC RLL	US CE	STERILE



BD MaxPlus™ Extension Sets



BD MaxPlus™ Extension Sets

Smallbore/Minibore Extension Set				
Part no. MP9236-C	IVA			 1mm
Length 36 cm (14 in) Priming vol. 0.75 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Extension Set Minibore 1 MaxPlus™ Bonded 1 SC RLL		US	STERILE
Part no. MP9014-C	IVA			 1mm
Length 36 cm (14 in) Priming vol. 0.75 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Extension Set Minibore 1 MaxPlus™ Removable 1 SC RLL		US	STERILE

BD MaxPlus™ Pressure Rated Extension Sets

Smallbore/Minibore Pressure Rated (325 psi) Extension Set				
Part no. MP5305-C	IVA			 1mm
Length 19 cm (7.5 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Extension Set Smallbore 1 MaxPlus™ Removable Pressure Rated 325 PSI RLL		US	STERILE
Part no. MP5304-C	IVA			 1mm
Length 18 cm (7 in) Priming vol. 0.6 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Extension Set Minibore 1 MaxPlus™ Bonded Pressure Rated 325 PSI 1 SC RLL		US	STERILE



Part no. MP5310-C	IVA	 1mm		
Length 21.5 cm (8.5 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxPlus™ Extension Set Minibore 1 MaxPlus™ Bonded Pressure Rated 325 PSI RLL		US		

Part no. MP5302-C	IVA	 1mm		
Length 21.5 cm (8.5 in) Priming vol. 0.6 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxPlus™ Extension Set Minibore 1 MaxPlus™ Bonded Pressure Rated 325 PSI 1 SC RLL		US		

Part no. MP5313-C	IVA	 1mm		
Length 22 cm (8.5 in) Priming vol. 0.6 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxPlus™ Extension Set Minibore 1 MaxPlus™ Removable Pressure Rated 325 PSI 1 SC RLL		US		

Standardbore Pressure Rated (325 psi) Extension Set

Part no. MP5301-C	IVA	 1.5mm		
Length 21.5 cm (8.5 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxPlus™ Extension Set Standardbore 1 MaxPlus™ Bonded Pressure Rated 325 PSI 1 PC RLL		US		


Part no. MP5314-C	IVA	 1.5mm		
Length 18 cm (7 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxPlus™ Extension Set Standardbore 1 MaxPlus™ Removable Pressure Rated 325 PSI 1 SC RLL		US		




Part no. MP5303-C	IVA	 1.5mm		
Length 19 cm (7.5 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxPlus™ Extension Set Standardbore 1 MaxPlus™ Removable Pressure Rated 325 PSI 1 PC RLL		US		

Part no. MP5312-C	IVA	 1.5mm		
Length 22 cm (8.5 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxPlus™ Extension Set Standardbore 1 MaxPlus™ Bonded Pressure Rated 325 PSI RLL		US		

BD MaxPlus™ T-connector Extension Sets



Smallbore/Minibore T-connector Extension Set				
Part no. MPXT5302-C	IVA	 1mm		
Length 19 cm (7.5 in) Priming vol. 1.1 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxPlus™ Extension Set Minibore 1 MaxPlus™ Bonded 1 SC T-connector		US		

BD MaxPlus™ Pressure Rated T-connector Extension Sets

Standardbore Pressure Rated (325 psi) T-connector Extension Set				
Part no. MPXT1001-C	IVA	 1.5mm		
Length 19 cm (7.5 in) Priming vol. 0.9 mL Ctn. qty. 50 Material PVC*				
Product description BD MaxPlus™ Extension Set Standardbore 1 MaxPlus™ Removable Pressure Rated 325 PSI 1 PC T-connector		US		




BD MaxPlus™ Pressure Rated Y-site Extension Sets


Standardbore Pressure Rated (325 psi) Y-site Extension Set				
Part no. MPX5305-C	IVA			 1.5mm
Length 19 cm (7.5 in) Priming vol. 1 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Extension Set Standardbore 1 MaxPlus™ Bonded Pressure Rated 325 PSI 1 SC 1 NAC-Y RLL		US	
Part no. MPX5301-C	IVA			 1.5mm
Length 19 cm (7.5 in) Priming vol. 1 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Extension Set Standardbore 1 MaxPlus™ Removable Pressure Rated 325 PSI 1 SC 1 NAC-Y RLL		US	
Part no. MPX5300-C	IVA			 1.5mm
Length 23 cm (9 in) Priming vol. 1.2 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Extension Set Standardbore 1 MaxPlus™ Bonded Pressure Rated 325 PSI 1 PC 1 NAC-Y RLL		US	
Part no. MPX5302-C	IVA			 1.5mm
Length 23 cm (9 in) Priming vol. 1.2 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Extension Set Standardbore 1 MaxPlus™ Removable Pressure Rated 325 PSI 1 PC 1 NAC-Y RLL		US	




BD MaxPlus™ Y-site Extension Sets with BCVs


Standardbore Y-site Extension Set with BCV				
Part no. MPX9108-C	IVA		 2.7mm	
Length 32 cm (13 in) Priming vol. 2.4 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Extension Set Standardbore 1 MaxPlus™ Bonded 1 BCV 1 NAC-Y RLL		US	STERILE

BD MaxPlus™ Bifuse Extension Sets

Smallbore/Minibore Bifuse Extension Set				
Part no. MP9232-C	IVA		 1mm	
Length 16 cm (6 in) Priming vol. 1.1 mL Ctn. qty. 50/600 Material PVC*				
	Product description BD MaxPlus™ Bifuse Extension Set Minibore 2 MaxPlus™ Bonded 2 SC FLL		US	STERILE



Standardbore Bifuse Extension Set				
Part no. MP9256-C	IVA		 2.7mm	
Length 18 cm (7 in) Priming vol. 1.4 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Bifuse Extension Set Standardbore 2 MaxPlus™ Bonded 2 SC RLL		US	STERILE


BD MaxPlus™ Bifuse Extension Sets Pressure Rated

Smallbore/Minibore Pressure Rated (325 psi) Bifuse Extension Set				
Part no. MP5315-C	IVA		 1mm	
Length 21 cm (8 in) Priming vol. 1.1 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Extension Set Minibore 2 MaxPlus™ Removable Pressure Rated 325 PSI 2 SC RLL		US	STERILE



BD MaxPlus™ Trifuse Extension Sets



Smallbore/Minibore Trifuse Extension Set				
Part no. MP9233-C	IVA			 1mm
Length 18 cm (7 in) Priming vol. 1.5 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Trifuse Extension Set Minibore 3 MaxPlus™ Bonded 3 SC RLL		US	STERILE
Standardbore Trifuse Extension Set				
Part no. MP9220-C	IVA			 2.7mm
Length 18 cm (7 in) Priming vol. 1.5 mL Ctn. qty. 50 Material PVC*				
	Product description BD MaxPlus™ Trifuse Extension Set Standardbore 3 MaxPlus™ Removable 3 SC RLL		US	STERILE



BD Q-Syte™ Extension Sets




BD Q-Syte™ Extension Sets

Smallbore/Minibore Extension Set			
Part no. 385102	IVA		 0.8mm
Length 15 cm (6 in) Priming vol. 0.34 mL Ctn. qty. 200 Material PVC/DEHP			
Product description BD Q-Syte™ Extension Set Smallbore 1 Q-Syte™ Removable 1 SC FLL		US CE 	
Part no. 385151	IVA		 1mm
Length 15 cm (6 in) Priming vol. 0.25 mL Ctn. qty. 50 Material PVC*			
Product description BD Q-Syte™ Extension Set Smallbore 1 Q-Syte™ Removable 1 SC RLL		US CE	
Standardbore Extension Set			
Part no. 385101	IVA		 2.7mm
Length 15 cm (6 in) Priming vol. 1.14 mL Ctn. qty. 200 Material PVC*			
Product description BD Q-Syte™ Extension Set Standardbore 1 Q-Syte™ Removable 1 SC RLL		US CE	
Part no. 385150	IVA		 2.7mm
Length 15 cm (6 in) Priming vol. 0.6 mL Ctn. qty. 50 Material PVC*			
Product description BD Q-Syte™ Extension Set Standardbore 1 Q-Syte™ Removable 1 SC RLL		US CE	




BD Q-Syte™ Y-site Extension Sets


Standardbore Y-site Extension Set				
Part no. 385165	IVA		 2.7mm	
Length 20 cm (8 in) Priming vol. 1.4 mL Ctn. qty. 50 Material PVC*				
	Product description BD Q-Syte™ Extension Set Standardbore 1 Q-Syte™ Removable 1 SC 1 Open Port Y w Cap RLL	US	CE	STERILE

BD Q-Syte™ Bifuse Extension Sets


Smallbore/Minibore Bifuse Extension Set				
Part no. 385163	IVA		 1mm	
Length 15 cm (6 in) Priming vol. 0.45 mL Ctn. qty. 50 Material PVC*				
	Product description BD Q-Syte™ Bifuse Extension Set Smallbore 2 Q-Syte™ Removable 2 SC RLL	US	CE	STERILE


Standardbore Bifuse Extension Set				
Part no. 385161	IVA		 2.7mm	
Length 15 cm (6 in) Priming vol. 1.6 mL Ctn. qty. 50 Material PVC*				
	Product description BD Q-Syte™ Bifuse Extension Set Standardbore 2 Q-Syte™ Removable 2 SC RLL	US	CE	STERILE

BD Q-Syte™ Trifuse Extension Sets

Smallbore/Minibore Trifuse Extension Set				
Part no. 385164	IVA		 1mm	
Length 15 cm (6 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*				
	Product description BD Q-Syte™ Trifuse Extension Set Smallbore 3 Q-Syte™ Removable 3 SC RLL	US	CE	STERILE






Standardbore Trifuse Extension Set			
Part no. 385162	IVA	 2.7mm	
Length 15 cm (6 in) Priming vol. 2.25 mL Ctn. qty. 50 Material PVC*			
Product description BD Q-Syte™ Trifuse Extension Set Standardbore 3 Q-Syte™ Removable 3 SC RLL		US CE	



BD SmartSite™ Extension Sets



BD SmartSite™ Extension Sets

Smallbore/Minibore Extension Set				
Part no. 10010511	IVA			 1mm
Length 15 cm (6 in) Priming vol. 0.2 mL Ctn. qty. 100 Material PVC*				
	Product description BD SmartSite™ Extension Set Smallbore 1 SmartSite™ Bonded 1 SC FLL	US		
Part no. 10011253	IVA			 1mm
Length 20 cm (8 in) Priming vol. 0.3 mL Ctn. qty. 100 Material PVC*				
	Product description BD SmartSite™ Extension Set Smallbore 1 SmartSite™ Removable 1 SC RLL	US		
Part no. 10010983	IVA			 1mm
Length 22 cm (8.5 in) Priming vol. 0.33 mL Ctn. qty. 100 Material PVC*				
	Product description BD SmartSite™ Extension Set Smallbore 1 SmartSite™ Removable 1 SC FLL	US		
Part no. 20049E	IVA			 1mm
Length 30 cm (12.5 in) Priming vol. 0.4 mL Ctn. qty. 100 Material PVC*				
	Product description BD SmartSite™ Extension Set Smallbore 1 SmartSite™ Bonded 1 SC RLL	US		
Standardbore Extension Set				
Part no. 20059E	IVA			 2.7mm
Length 22 cm (8.5 in) Priming vol. 1 mL Ctn. qty. 100 Material PVC/DEHP				
	Product description BD SmartSite™ Extension Set Standardbore 1 SmartSite™ Bonded 1 SC RLL	US		




Part no. 10010912	IVA		2.7mm
Length 23 cm (9 in) Priming vol. 1.3 mL Ctn. qty. 100 Material PVC*			
Product description BD SmartSite™ Extension Set Standardbore 1 SmartSite™ Bonded 1 SC RLL		US	
Part no. 22059E	IVA		2.7mm
Length 23 cm (9 in) Priming vol. 1 mL Ctn. qty. 100 Material PVC/DEHP			
Product description BD SmartSite™ Extension Set Standardbore 1 SmartSite™ Removable 1 SC RLL		US	

BD SmartSite™ Pressure Rated Extension Sets


Smallbore/Minibore Pressure Rated (325 psi) Extension Set			
Part no. 20039E	IVA		1mm
Length 15 cm (6 in) Priming vol. 0.2 mL Ctn. qty. 100 Material PVC*			
Product description BD SmartSite™ Extension Set Smallbore 1 SmartSite™ Bonded Pressure Rated 325 PSI 1 SC RLL		US	
Standardbore Pressure Rated (325 psi) Extension Set			
Part no. 22003E-07	IVA		1.5mm
Length 22 cm (8.5 in) Priming vol. 0.5 mL Ctn. qty. 100 Material PVC*			
Product description BD SmartSite™ Extension Set Standardbore 1 SmartSite™ Bonded Pressure Rated 325 PSI 1 PC RLL		US	





BD SmartSite™ T-connector Extension Sets

Smallbore/Minibore T-connector Extension Set				
Part no. C25013E	IVA		 1mm	
Length 18 cm (7 in) Priming vol. 0.4 mL Ctn. qty. 100 Material PVC*				
Product description BD SmartSite™ T-connector Extension Set Smallbore 2 SmartSite™ Removable 1 SC RLL			US	

BD SmartSite™ Extension Sets with Texium™ Male Luer

Microbore Extension Set with 1 bonded SmartSite™ Connector and Texium Male Luer				
Part no. 22005E-07T	GSD		 0.5mm	
Length 157 cm (62 in) Priming vol. 1.6 mL Ctn. qty. 50 Material PVC*				
Product description BD SmartSite™ Extension Set Microbore 1 SmartSite™ Bonded 1 PC Texium™			US	

BD SmartSite™ Bifuse Extension Sets

Smallbore/Minibore Bifuse Extension Sets				
Part no. 10015488	IVA		 1mm	
Length 20 cm (8 in) Priming vol. 0.5 mL Ctn. qty. 100 Material PVC*				
Product description BD SmartSite™ Bifuse Extension Set Smallbore 2 SmartSite™ Removable 2 SC RLL			US	
Part no. 20019E	IVA		 1mm	
Length 24 cm (9.5 in) Priming vol. 0.4 mL Ctn. qty. 100 Material PVC*				
Product description BD SmartSite™ Bifuse Extension Set Smallbore 2 SmartSite™ Bonded 2 SC RLL			US	



Part no. 10012143	IVA	
Length 28 cm (11 in) Priming vol. 1.5 mL Ctn. qty. 100 Material PVC/DEHP		
Product description BD SmartSite™ Bifuse Extension Set Standardbore 1 Open Female Luer 1 SmartSite™ Removable 2 SC RLL	US	

Standardbore Bifuse Extension Sets		
Part no. 20159E	IVA	
Length 25 cm (10 in) Priming vol. 2 mL Ctn. qty. 100 Material PVC*		
Product description BD SmartSite™ Bifuse Extension Set Standardbore 2 SmartSite™ Bonded 2 SC RLL	US	

BD SmartSite™ Trifuse Extension Sets

Smallbore/Minibore Trifuse Extension Set		
Part no. 20038E	IVA	
Length 13 cm (5 in) Priming vol. 0.3 mL Ctn. qty. 100 Material PVC*		
Product description BD SmartSite™ Trifuse Extension Set Smallbore 3 SmartSite™ Bonded 3 SC RLL	US	
Part no. 20338E	IVA	
Length 20 cm (8 in) Priming vol. 0.7 mL Ctn. qty. 100 Material PVC*		
Product description BD SmartSite™ Trifuse Extension Set Smallbore 2 Open Female Luer 1 SmartSite™ Removable 4 SC RLL	US	



Standardbore Trifuse Extension Set				
Part no. 10014692	IVA			 2.7mm
Length 25 cm (10 in) Priming vol. 2 mL Ctn. qty. 100 Material PVC*				
Product description BD SmartSite™ Trifuse Extension Set Standardbore 1 Open Female Luer 2 SmartSite™ Bonded 4 SC FLL			US	

BD SmartSite™ Trifuse Extension Sets with Back Check Valves

Smallbore/Minibore Trifuse Extension Set with 3 bonded SmartSite™ Connector and 3 BCV				
Part no. 20062E	IVA			 1mm
Length 18 cm (7 in) Priming vol. 1.7 mL Ctn. qty. 100 Material PVC/DEHP				
Product description BD SmartSite™ Trifuse Extension Set Smallbore 3 SmartSite™ Bonded 1 BCV 3 SC RLL			US	

BD SmartSite™ Multiport Extension Sets

Standardbore Dual Multiport Extension Set				
Part no. 10010971	IVA			 2.7mm
Length 109 cm (43 in) Priming vol. 6 mL Ctn. qty. 100 Material PVC/DEHP				
Product description BD SmartSite™ Multiport Extension Set Standardbore 2 NF Y Multiport 2 PC 2 SmartSite™ Y RLL			US	



BD SmartSite™ Filter Extension Sets

Standardbore 0.2um Large PES Filter			
Part no. 10013902	IVA	 	
Length 41 cm (16 in) Priming vol. 4 mL Ctn. qty. 100 Material PE lined PVC/DEHP			
Product description BD SmartSite™ Filter Extension Set Standardbore 1 SmartSite™ Bonded 0.2um Large PES Filter 1 PC RLL		US 	
Macrobore Extension Sets with 1.2um Small Volume PES Filter			
Part no. 20127E	IVA	 	
Length 46 cm (18 in) Priming vol. 5 mL Ctn. qty. 100 Material PVC*			
Product description BD SmartSite™ Filter Extension Set Standardbore 1.2um Large PES Filter 1 PC 1 SmartSite™ Y RLL		US	



Extension Sets



Extension Sets

Microbore Extension Sets				
Part no. ME2020	GSD			 0.5mm
Length 152 cm (60 in) Priming vol. 0.3 mL Ctn. qty. 50/500 Material PVC*				
	Product description BD Extension Set Microbore 1 SC RLL	US		
Part no. ME2059	GSD			 0.5mm
Length 183 cm (72 in) Priming vol. 0.4 mL Ctn. qty. 50/500 Material PVC*				
	Product description BD Extension Set Microbore RLL	US		
Part no. 30910	IVA			 0.5mm
Length 94 cm (37 in) Priming vol. 0.4 mL Ctn. qty. 100 Material PVC*				
	Product description BD Extension Set Microbore 1 SC FLL	US		
Smallbore/Minibore Extension Sets				
Part no. 10373772	IVA			 1mm
Length 15 cm (6 in) Priming vol. 0.3 mL Ctn. qty. 100 Material PVC*				
	Product description BD Extension Set Smallbore 1 SC RLL	US		



Part no. 11088483		GSD		
Length	158 cm (62 in)			
Priming vol.	1.5 mL			
Ctn. qty.	100			
Material	PVC*	Product description BD Extension Set Smallbore 1 SC RLL	US	
Part no. 10896354		GSD		
Length	279 cm (110 in)			
Priming vol.	6 mL			
Ctn. qty.	100			
Material	PVC*	Product description BD Extension Set Smallbore 1 PC RLL	US	
Part no. 10031776		GSD		
Length	594 cm (234 in)			
Priming vol.	12 mL			
Ctn. qty.	50			
Material	PVC*	Product description BD Extension Set Smallbore 1 PC FLL	US	
Part no. ME2010		GSD		
Length	152 cm (60 in)			
Priming vol.	1.8 mL			
Ctn. qty.	50/500			
Material	PVC*	Product description BD Extension Set Minibore 1 SC RLL	US	
Part no. ME1048		IVA		
Length	30 cm (12 in)			
Priming vol.	0.4 mL			
Ctn. qty.	50			
Material	PVC*	Product description BD Extension Set Minibore 1 SC RLL	US	



Part no. C20016	IVA			 1.5mm
Length 76 cm (30 in) Priming vol. 2 mL Ctn. qty. 72 Material PVC*				
Product description BD Extension Set Smallbore Pressure Rated 325 PSI FLL		US		

Standardbore Extension Sets

Part no. ME2068	GSD			 2.7mm
Length 152 cm (60 in) Priming vol. 9 mL Ctn. qty. 50 Material PVC*				
Product description BD Extension Set Standardbore FLL		US		

Part no. C20022	IVA			 2.7mm
Length 51 cm (20 in) Priming vol. 3 mL Ctn. qty. 72 Material PVC*				
Product description BD Extension Set Standardbore 1 PC RLL		US		

Part no. ME1055	IVA			 2.7mm
Length 76 cm (30 in) Priming vol. 4.9 mL Ctn. qty. 50 Material PVC*				
Product description BD Extension Set Standardbore 1 SC RLL		US		

Part no. ME1217	IVA			 2.7mm
Length 76 cm (30 in) Priming vol. 5 mL Ctn. qty. 50 Material PVC*				
Product description BD Extension Set Standardbore 1 RC RLL		US		



Pressure Rated Extension Sets



Smallbore/Minibore Pressure Rated Extension Sets				
Part no. ME5302	IVA			 1mm
Length 16 cm (6.5 in) Priming vol. 0.2 mL Ctn. qty. 50/500 Material PVC*				
	Product description BD Extension Set Minibore Pressure Rated 325 PSI 1 SC RLL	US		
Part no. C30006	GSD			 1mm
Length 160 cm (63 in) Priming vol. 1.6 mL Ctn. qty. 100 Material PVC*				
	Product description BD Extension Set Smallbore Pressure Rated 325 PSI RLL	US		
Standardbore Pressure Rated (325 psi) Extension Sets				
Part no. ME5301	IVA			 1.5mm
Length 19 cm (7.5 in) Priming vol. 0.4 mL Ctn. qty. 50 Material PVC*				
	Product description BD Extension Set Standardbore Pressure Rated 325 PSI 1 PC RLL	US		

Extension Set with Back Check Valve and Smartsite™ Y-site



Extension Sets with Back Check Valve and SmartSite™ Y-sites				
Part no. 30843E	MMS			 2.7mm
Length 30 cm (12 in) Priming vol. 1 mL Ctn. qty. 100 Material PVC*				
	Product description BD Extension Set Standardbore 1 BCV 1 PC 1 SmartSite™ Y RLL	US		



T-connector Extension Sets


Smallbore/Minibore T-connector Extension Sets				
Part no. 20043E	IVA		 1mm	
Length 15 cm (6 in) Priming vol. 0.4 mL Ctn. qty. 100 Material PVC*				
Product description BD T-connector Extension Set Smallbore 1 SC SmartSite™/RLL T-connector		US		
Part no. 20051E	IVA		 1mm	
Length 15 cm (6 in) Priming vol. 0.2 mL Ctn. qty. 100 Material PVC*				
Product description BD T-connector Extension Set Smallbore 1 SC SmartSite™ T LS		US		


T-connector Extension Sets Pressure Rated

Smallbore/Minibore Pressure Rated (325 psi) T-connector Extension Sets				
Part no. MXT1003	IVA		 1mm	
Length 15 cm (6 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*				
Product description BD T-connector Extension Set Smallbore Pressure Rated 325 PSI 1 SC Nac/RLL T-connector		US		
Part no. 20041E	IVA		 1mm	
Length 15 cm (6 in) Priming vol. 0.2 mL Ctn. qty. 100 Material PVC*				
Product description BD T-connector Extension Set Smallbore Pressure Rated 325 PSI 1 SC SmartSite™/RLL T-connector		US		



J-loop Extension Sets

Smallbore/Minibore J-loop Extension Sets				
Part no. 20021E		IVA	 1mm	
Length	15 cm (6 in)			
Priming vol.	0.2 mL			
Ctn. qty.	100	Product description		US
Material	PVC*	BD J-LOOP Extension Set Smallbore 1 SmartSite™ Bonded 1 SC RLL		




Extension Sets with Anti-Siphon Valve



Extension Sets with Anti-Siphon Valves




Smallbore Extension Sets with Anti-Syphon Valve				
Part no. C25012	IVA			 1mm
Length 13 cm (5 in) Priming vol. 0.5 mL Ctn. qty. 20 Material PVC*				
	Product description BD Extension Set PVC Smallbore 1 ASV (Removable) FLL	US		
Part no. 10010525	IVA			 1mm
Length 76 cm (30 in) Priming vol. 1.2 mL Ctn. qty. 100 Material PVC*				
	Product description BD Extension Set Smallbore 1 PC 1 ASV RLL	US		




Extension Sets with Y-site Access Ports



Extension Sets with Y-sites

Smallbore/Minibore Y-site Extension Sets				
Part no. C25005E	GSD			 1mm
Length 366 cm (144 in) Priming vol. 8 mL Ctn. qty. 100 Material PVC*				
Product description BD Extension Set Smallbore 2 PC 2 SmartSite™ Y RLL			US	
Part no. 30212E	IVA			 1mm
Length 53 cm (21 in) Priming vol. 2 mL Ctn. qty. 100 Material PVC*				
Product description BD Extension Set Smallbore 1 PC 1 SmartSite™ Y RLL			US	
Part no. C25002E	GSD			 1mm
Length 599 cm (236 in) Priming vol. 12 mL Ctn. qty. 50 Material PVC*				
Product description BD Extension Set Smallbore 4 SmartSite™ Y FLL			US	
Standardbore Y-site Extension Sets				
Part no. 37262E	GSD			 2.7mm
Length 104 cm (41 in) Priming vol. 6 mL Ctn. qty. 100 Material PVC*				
Product description BD Extension Set Standardbore 2 PC 2 SmartSite™ Y FLL			US CE	



Part no. MX9060		IVA	 2.7mm
Length 20 cm (8 in) Priming vol. 1.2 mL Ctn. qty. 50 Material PVC*			
Product description BD Extension Set Standardbore 1 SC 1 NAC-Y RLL		US	
Part no. 11590099		IVA	 2.7mm
Length 25 cm (10 in) Priming vol. 1.4 mL Ctn. qty. 100 Material PVC*			
Product description BD Extension Set Standardbore 1 PC 1 SmartSite™ Y RLL		US	
Part no. C25004E		GSD	 2.7mm
Length 462 cm (182 in) Priming vol. 23 mL Ctn. qty. 100 Material PVC/DEHP			
Product description BD Extension Set Standardbore 1 PC 1 SmartSite™ Y RLL		US	 
Part no. 30202E		IVA	 2.7mm
Length 51 cm (20 in) Priming vol. 3 mL Ctn. qty. 100 Material PVC/DEHP			
Product description BD Extension Set Standardbore 1 PC 1 SmartSite™ Y RLL		US	 
Part no. MX9175		IVA	 2.7mm
Length 81 cm (32 in) Priming vol. 3.5 mL Ctn. qty. 50 Material PVC*			
Product description BD Extension Set Standardbore 1 SC 1 NAC-Y RLL		US	



Part no. MX9059	IVA		2.7mm
Length 91 cm (36 in) Priming vol. 6 mL Ctn. qty. 50 Material PVC*			
Product description BD Extension Set Standardbore 1 SC 2 NAC-Y RLL		US	STERILE

Extension Sets with Y-sites Pressure Rated


Smallbore/Minibore Pressure Rated (325 psi) Y-site Extension Sets			
Part no. MX5302	IVA		1mm
Length 15 cm (6 in) Priming vol. 0.7 mL Ctn. qty. 50 Material PVC*			
Product description BD Extension Set Minibore Pressure Rated 325 PSI 1 SC 1 NAC-Y RLL		US	STERILE

Standardbore Pressure Rated (325 psi) Y-site Extension Sets			
Part no. MX5301	IVA		1.5mm
Length 20 cm (8 in) Priming vol. 0.8 mL Ctn. qty. 50 Material PVC*			
Product description BD Extension Set Standardbore Pressure Rated 325 PSI 1 PC 1 NAC-Y RLL		US	STERILE


Extension Sets with Y-sites and Back Check Valves

Standardbore Y-site Extension Sets			
Part no. MX9178	GSD		2.7mm
Length 127 cm (50 in) Priming vol. 8.4 mL Ctn. qty. 50 Material PVC*			
Product description BD Extension Set Standardbore 1 BCV 2 NAC-Y RLL		US	STERILE




Part no. MX9182	IVA	 2.7mm		
Length 20 cm (8 in) Priming vol. 1.4 mL Ctn. qty. 50 Material PVC*				
Product description BD Extension Set Standardbore 1 BCV 1 NAC-Y RLL		US		

Extension Sets with Sideport


Extension Sets with 2 Sideport and BCV				
Part no. 30207E	IVA	 2.7mm		
Length 43 cm (17 in) Priming vol. 3 mL Ctn. qty. 100 Material PVC/DEHP				
Product description BD Bifuse Extension Set Standardbore 1 BCV 1 PC 1 SmartSite™ Y RLL		US		

Bifuse Extension Sets


Bifuse Extension Sets				
Part no. 10794983	IVA	 1mm		
Length 18 cm (7 in) Priming vol. 0.3 mL Ctn. qty. 100 Material PVC*				
Product description BD Bifuse Extension Set Smallbore 2 Open Female Luer 2 SC RLL		US		




Trifuse Extension Sets

Trifuse Extension Sets			
Part no. 10839681	IVA		 1mm
Length 18 cm (7 in) Priming vol. 0.5 mL Ctn. qty. 100 Material PVC*			
Product description BD Trifuse Extension Set Smallbore 3 Open Female Luer 4 SC RLL		US	

Quadfuse Extension Sets

Quadfuse Extension Sets			
Part no. 10015817	IVA		 1mm
Length 20 cm (8 in) Priming vol. 0.7 mL Ctn. qty. 100 Material PVC*			
Product description BD Quadfuse Extension Set Smallbore 4 Open Female Luer 5 SC RLL		US	




Specialty Tubing Extension Sets



PE Lined PVC Extension Sets

Standardbore PE Lined PVC Extension Sets				
Part no. C20014	IVA			 2.7mm
Length 76 cm (30 in) Priming vol. 4 mL Ctn. qty. 72 Material PE lined PVC*				
	Product description BD Extension Set Standardbore RLL		US	



In-Line Filter Extension Sets

0.2 μ m and 1.2 μ m



0.2 Micron Filter Extension Sets

Extension Sets with Small 0.2um Filter				
Part no. 10011865	IVA		 0.2 µm	 1mm
Length 18 cm (7 in) Priming vol. 0.5 mL Ctn. qty. 100 Material PVC*				
Product description BD Filter Extension Set Smallbore Small 0.2um Filter (+ve) 1 SC RLL			US	
Extension Sets with 0.2um Small Volume PES Filter				
Part no. 10010570	GSD		 0.2 µm	 0.5mm
Length 152 cm (60 in) Priming vol. 1 mL Ctn. qty. 100 Material PVC*				
Product description BD Filter Extension Set Microbore 0.2um Small Volume PES Filter 1 SC FLL			US	
Part no. 20029E	IVA		 0.2 µm	 1mm
Length 28 cm (11 in) Priming vol. 1 mL Ctn. qty. 100 Material PVC*				
Product description BD Filter Extension Set Smallbore 0.2um Small Volume PES Filter 1 SC 1 SmartSite™ Y RLL			US	
Part no. 20027E	IVA		 0.2 µm	 2.7mm
Length 46 cm (18 in) Priming vol. 5 mL Ctn. qty. 100 Material PVC*				
Product description BD Filter Extension Set Standardbore 0.2um Large PES Filter 1 PC 1 SmartSite™ Y RLL			US	



0.2 Micron Filter PE Lined PVC Extension Sets

Standardbore 0.2um Large PES Filter				
Part no. MX9166	IVA		0.2 μm 2.7mm	
Length 41 cm (16 in) Priming vol. 6 mL Ctn. qty. 50 Material PE lined PVC*				
	Product description BD Filter Extension Set Standardbore 0.2um Large PES Filter 1 PC 1 NAC-Y RLL	US	STERILE	
Part no. 20350E	IVA		0.2 μm 2.7mm	
Length 45 cm (18 in) Priming vol. 5 mL Ctn. qty. 100 Material PE lined PVC*				
	Product description BD Filter Extension Set Standardbore 0.2um Large PES Filter 1 PC 1 SmartSite™ Y RLL	US	STERILE	


0.2 Micron Filter PE Lined PVC Extension Sets with Texium™ Male Luer

Extension Sets with 0.2um Large Volume PES Filter and Texium™ Male Luer				
Part no. 20350ET	GSD		0.2 μm 2.7mm	
Length 51 cm (20 in) Priming vol. 6 mL Ctn. qty. 100 Material PE lined PVC/DEHP				
	Product description BD Filter Extension Set Standardbore 0.2um Large PES Filter 1 PC 1 SmartSite™ Y Texium™	US	PHT DEHP	STERILE

1.2 Micron Filter Extension Sets

Macrobores Extension Sets with 1.2um Small Volume PES Filter				
Part no. 20129E	IVA		1.2 μm 1mm	
Length 18 cm (7 in) Priming vol. 0.5 mL Ctn. qty. 100 Material PVC*				
	Product description BD Filter Extension Set Smallbore 1.2um Small Volume PES Filter 1 SC FLL	US	STERILE	




Part no. ME1225	IVA	 2.7mm
Length 38 cm (15 in) Priming vol. 4.6 mL Ctn. qty. 50 Material PVC*		
Product description BD Filter Extension Set Standardbore 1.2um Large PES Filter 1 SC RLL		US 

Part no. 10012866	IVA	 2.7mm
Length 43 cm (17 in) Priming vol. 4 mL Ctn. qty. 100 Material PVC*		
Product description BD Filter Extension Set Standardbore 1.2um Large PES Filter 1 PC RLL		US 

1.2 Micron Filter PE Lined PVC Extension Sets

Extension Sets with 1.2um Large Volume PES Filter		
Part no. MX9171	IVA	 2.7mm
Length 41 cm (16 in) Priming vol. 6 mL Ctn. qty. 50 Material PE lined PVC*		
Product description BD Filter Extension Set Standardbore 1.2um Large PES Filter 1 PC 1 NAC-Y RLL		US 



Stopcock Extension Sets





Stopcock Extension Sets

Standardbore Removable 3 Way Single Handle Stopcock Sets				
Part no. 25203	IVA		 2.7mm	
Length 51 cm (20 in) Priming vol. 3 mL Ctn. qty. 72 Material PVC*				
	Product description BD Stopcock Extension Set Standardbore Removable 3WSC HiFlow RLL	US		


Standardbore Removable 4 Way White T-handle Stopcock Sets				
Part no. 25206	IVA		 2.7mm	
Length 81 cm (32 in) Priming vol. 4 mL Ctn. qty. 72 Material PVC*				
	Product description BD Stopcock Extension Set Standardbore Removable 4WSC RLL	US		

Stopcock Extension Sets with Y-sites

Smallbore/Minibore Removable 3 Way T-handle Stopcock Sets				
Part no. MX4464	IVA		 1mm	
Length 110 cm (43.5 in) Priming vol. 7 mL Ctn. qty. 50 Material PVC*				
	Product description BD Extension Set PVC Smallbore 1 PC 1 NAC-Y 4WSC (Removable) RLL (Removable Extension 1 PC 1 NAC-Y 17cm/6.5in)	US		


Standardbore Removable 4 Way White T-handle Stopcock Sets with Y-sites				
Part no. MX4439	GSD		 2.7mm	
Length 102 cm (40 in) Priming vol. 7 mL Ctn. qty. 50 Material PVC*				
	Product description BD Stopcock Extension Set Standardbore 2 NAC-Y Removable 4WSC FLL	US		




Standardbore Bonded 4 Way White Single Handle HiFlow Stopcock Sets with Y-sites			
Part no. MX4471	IVA		 2.7mm
Length 81 cm (32 in) Priming vol. 6 mL Ctn. qty. 50 Material PVC*			
Product description BD Stopcock Extension Set Standardbore 1 NAC-Y 4WSC HiFlow RLL		US	

Standardbore Removable 4 Way White T-handle Stopcock Sets with Y-sites			
Part no. 25202E	IVA		 2.7mm
Length 94 cm (37 in) Priming vol. 5 mL Ctn. qty. 50 Material PVC/DEHP			
Product description BD Stopcock Extension Set Standardbore 1 PC 1 SmartSite™ Y 4WSC Removable RLL		US 	

Ganged 4 Way Stopcock Extension Sets


Standardbore Removable 2 Gang 4 Way White Single Handle Stopcock Sets with Y-sites			
Part no. MX4450	GSD		 1mm
Length 163 cm (64 in) Priming vol. 6.5 mL Ctn. qty. 50 Material PVC*			
Product description BD Stopcock Extension Set Standardbore 1 SC 1 NAC-Y 2 Gang 4WSC RLL		US	


Standardbore Removable 2 Gang 4 Way White T-handle Stopcock Sets with Y-sites			
Part no. 25207E	IVA		 2.7mm
Length 86 cm (34 in) Priming vol. 5 mL Ctn. qty. 50 Material PVC*			
Product description BD Stopcock Extension Set Standardbore 1 PC 1 SmartSite™ Y 2 Gang 4WSC Removable RLL		US	





Standardbore Removable 3 Gang 4 Way White T-handle Stopcock Manifold Sets with Y-sites			
Part no. 25208E	IVA		 2.7mm
Length 86 cm (34 in) Priming vol. 5 mL Ctn. qty. 50 Material PVC*			
	Product description BD Stopcock Extension Set Standardbore 1 SmartSite™ Y 3 Gang 4WSC RLL	US	

Manifold Extension Sets

Standardbore 3 Port Manifold Sets with BCV and Y-sites			
Part no. 20519E	IVA		 2.7mm
Length 109 cm (41 in) Priming vol. 8 mL Ctn. qty. 50 Material PVC/DEHP			
	Product description BD Manifold Extension Set Standardbore 1 PC 1 SmartSite™ Y 3 Port Manifold w BCV (Removable) RLL	US	 

Standardbore 3 Port Manifold Sets with BCV			
Part no. MM4021	GSD		 2.7mm
Length 112 cm (44 in) Priming vol. 7 mL Ctn. qty. 50 Material PVC*			
	Product description BD Manifold Extension Set Standardbore 3 Port Manifold w BCV (Removable) RLL	US	


Standardbore Bonded 4 Way White T-handle Manifold Sets with 2 One Way Valves			
Part no. 20511E	IVA		 2.7mm
Length 114 cm (45 in) Priming vol. 8 mL Ctn. qty. 50 Material PVC*			
	Product description BD Manifold Extension Set Standardbore 1 PC 1 SmartSite™ Y RLL (Removable Extension 4WSC 2 T Port w BCV Set 25 cm/9 in)	US	 





In-Line Flow Controller Extension Sets



Flow Controller Single Scale with Y-sites

Standardbore Flow Controller Single Scale with Y-sites				
Part no. MFS106		GSD	 2.7mm	
Length	58 cm (23 in)			
Priming vol.	3.7 mL			
Ctn. qty.	50			
Material	PVC*			
		Product description BD Flow Controller Extension Set Standardbore Single Scale Flow Controller 1 SC 1 NAC-Y RLL	US	

Gravity Sets



Sales and Customer Service


 **1-844-8BD LIFE (1-844-823-5433)**

Clinical Support

 **1-800-555-7422 Ext 5**

 **medical.information@bd.com**







Primary Gravity Sets








BD Primary Gravity Sets

Gravity Infusion Sets				
Part no. MA3131	GSD		 10 ml	 2.7mm
Length 213 cm (84 in) Priming vol. 20 mL Ctn. qty. 50 Material PVC*				
Product description BD Gravity Infusion Set Standardbore V/NV 10 DPML 1 RC 1 SC RLL			US	STERILE
Part no. C42001	GSD		 20 ml	 2.7mm
Length 229 cm (90 in) Priming vol. 19 mL Ctn. qty. 100 Material PVC*				
Product description BD Gravity Infusion Set Standardbore V/NV 20 DPML 15um 1 RC RLL			US	STERILE





BD Primary Gravity Set Y Access Basic

Gravity Infusion Sets with Y Access Basic				
Part no. 41173E	GSD		 20 ml	 2.7mm
Length 264 cm (104 in) Priming vol. 23 mL Ctn. qty. 100 Material PVC*				
Product description BD Gravity Infusion Set Standardbore V/NV 20 DPML 1 RC 1 SmartSite™ Y RLL			US	STERILE
Part no. 42133E	GSD		 60 ml	 2.7mm
Length 264 cm (104 in) Priming vol. 17 mL Ctn. qty. 100 Material PVC*				
Product description BD Gravity Infusion Set Standardbore V/NV 60 DPML 1 RC 2 SmartSite™ Y RLL			US	STERILE








Part no. 42273E		GSD	 20 ml	 2.7mm
Length	264 cm (104 in)			
Priming vol.	21 mL			
Ctn. qty.	100			
Material	PVC*			
		Product description BD Gravity Infusion Set Standardbore V/NV 20 DPML 1 RC 2 SmartSite™ Y RLL	US	
Part no. 10015171		GSD	 20 ml	 2.7mm
Length	264 cm (104 in)			
Priming vol.	31 mL			
Ctn. qty.	20			
Material	PVC*			
		Product description BD Gravity Infusion Set Standardbore Dual Drip Chamber V/NV 20 DPML 3 RC 1 SmartSite™ Y RLL	US	
Part no. MX9107		GSD	 10 ml	 2.7mm
Length	229 cm (90 in)			
Priming vol.	21 mL			
Ctn. qty.	50			
Material	PVC*			
		Product description BD Gravity Infusion Set Standardbore V/NV 10 DPML 1 RC 1 NAC-Y RLL	US	
Part no. 42103E-07		GSD	 15 ml	 2.7mm
Length	254 cm (100 in)			
Priming vol.				
Ctn. qty.	100			
Material	PVC*			
		Product description BD Gravity Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 2 SC 3 SmartSite™ Y RLL	US	
Part no. 47233E		GSD	 60 ml	 2.7mm
Length	264 cm (104 in)			
Priming vol.	22 mL			
Ctn. qty.	100			
Material	PVC*			
		Product description BD Gravity Infusion Set Standardbore V/NV 60 DPML 1 BCV 1 RC 2 SmartSite™ Y RLL	US	





Part no. 47177E		GSD	 10 ml	 2.7mm
Length	264 cm (104 in)			
Priming vol.	23 mL			
Ctn. qty.	100			
Material	PVC*			
		Product description BD Gravity Infusion Set Standardbore V/NV 10 DPML 1 BCV 1 RC 2 SmartSite™ Y RLL	US	
Part no. 42434E		GSD	 60 ml	 2.7mm
Length	264 cm (104 in)			
Priming vol.	23 mL			
Ctn. qty.	100			
Material	PVC*			
		Product description BD Gravity Infusion Set Standardbore V/NV 60 DPML 1 BCV 1 RC 3 SmartSite™ Y RLL	US	
Part no. 42493E		GSD	 20 ml	 2.7mm
Length	264 cm (104 in)			
Priming vol.	17 mL			
Ctn. qty.	100			
Material	PVC/DEHP			
		Product description BD Gravity Infusion Set Standardbore V/NV 20 DPML 1 BCV 1 RC 3 SmartSite™ Y RLL	US	 
Part no. 10802688		GSD	 10 ml	 2.7mm
Length	269 cm (106 in)			
Priming vol.	20 mL			
Ctn. qty.	100			
Material	PVC*			
		Product description BD Gravity Infusion Set Standardbore V/NV 10 DPML 1 BCV 1 RC 3 SmartSite™ Y RLL	US	
Part no. 41138E-07		GSD	 10 ml	 2.7mm
Length	345 cm (136 in)			
Priming vol.	28 mL			
Ctn. qty.	50			
Material	PVC*			
		Product description BD Gravity Infusion Set Standardbore V/NV 10 DPML 1 BCV 1 RC 3 SmartSite™ Y RLL (Removable Extension 1 SC 18 cm/7 in)	US	




Part no. 41137E	GSD	 
Length 345 cm (136 in) Priming vol. 27 mL Ctn. qty. 50 Material PVC*		
Product description BD Gravity Infusion Set Standardbore V/NV 10 DPML 1 BCV 1 RC 4 SmartSite™ Y RLL (Removable Extension 1 SC 8 cm/3 in)		US
		
Part no. MX9963	GSD	 
Length 271 cm (107 in) Priming vol. 24 mL Ctn. qty. 50 Material PVC*		
Product description BD Gravity Infusion Set Standardbore V/NV 60 DPML 1 BCV 1 RC 1 SC 2 NAC-Y RLL		US
		
Part no. MX9968	GSD	 
Length 279 cm (110 in) Priming vol. 24 mL Ctn. qty. 50 Material PVC*		
Product description BD Gravity Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 1 SC 2 NAC-Y RLL		US
		
Part no. MX9416	GSD	 
Length 279 cm (110 in) Priming vol. 25 mL Ctn. qty. 50 Material PVC*		
Product description BD Gravity Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 2 NAC-Y RLL		US
		
Part no. MX9128	GSD	 
Length 284 cm (112 in) Priming vol. 25 mL Ctn. qty. 50 Material PVC*		
Product description BD Gravity Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 1 SC 3 NAC-Y RLL		US
		




Part no. MX9433	GSD	 15 ml	 2.7mm
Length 323 cm (127 in) Priming vol. 28 mL Ctn. qty. 50 Material PVC*			
Product description BD Gravity Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 1 SC 4 NAC-Y RLL		US	STERILE

Part no. MX4457	GSD	 15 ml	 2.7mm
Length 374 cm (147 in) Priming vol. 31 mL Ctn. qty. 50 Material PVC*			
Product description BD Gravity Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 2 NAC-Y RLL (1 Removable Extensions Set 2 NAC-Y 1 SC 86 cm/34 in and 1 Extension Pressure Rated 1PC 19 cm/7.5 in)		US	STERILE




BD Primary Gravity Sets with In-Line Filter

Gravity Sets with In-Line Filters			
Part no. 42933E	GSD	 20 ml	 0.2 μm
Length 264 cm (104 in) Priming vol. 24 mL Ctn. qty. 50 Material PVC/DEHP			
Product description BD Gravity Infusion Set Standardbore V/NV 20 DPML 1 BCV 0.2um Large PES Filter 1 RC 1 PC 2 SmartSite™ Y RLL		US	PHT DEHP STERILE


BD Primary Gravity Flow Controller Sets

Gravity Flow Controller Sets			
Part no. MFS137	GSD	 20 ml	 0.2 μm
Length 188 cm (74 in) Priming vol. Ctn. qty. 50 Material PVC*			
Product description BD Gravity Flow Controller Infusion Set Standardbore V/NV 20 DPML Single Scale Flow Controller 0.2um Large PES Filter 1 SC RLL		US	STERILE







Part no. MFS140	GSD	 20 ml	 2.7mm
Length 216 cm (84 in) Priming vol. 20 mL Ctn. qty. 50 Material PVC*			
Product description BD Gravity Flow Controller Infusion Set Standardbore V/NV 20 DPML Single Scale Flow Controller 1 PC 1 NAC-Y RLL		US	STERILE
Part no. MFS102	GSD	 10 ml	 2.7mm
Length 251 cm (99in) Priming vol. 22 mL Ctn. qty. 50 Material PVC*			
Product description BD Gravity Flow Controller Infusion Set Standardbore V/NV 10 DPML Single Scale Flow Controller 1 SC 1 NAC-Y RLL		US	STERILE
Part no. MFS132	GSD	 20 ml	 2.7mm
Length 226 cm (89 in) Priming vol. 22 mL Ctn. qty. 50 Material PVC*			
Product description BD Gravity Flow Controller Infusion Set Standardbore V/NV 20 DPML 1 BCV Single Scale Flow Controller 1 PC 2 NAC-Y RLL		US	STERILE

BD Primary Gravity Burette (Add On) Sets




Gravity Burette (Add On) Sets			
Part no. 82113E	GSD	 2.7mm	
Length 56 cm (22 in) Priming vol. 33 mL Ctn. qty. 20 Material PVC*			
Product description BD Gravity Burette (Add On) Infusion Set Standardbore V/NV 1 RC 1 PC Graduated 150 mL (Top Plate V/NV 1 NF Access) Spike Adapter		US	STERILE







BD Primary Gravity Stopcock Sets

Gravity Stopcock Sets				
Part no. 42511E	GSD			
Length 307 cm (121 in) Priming vol. 23 mL Ctn. qty. 50 Material PVC*		Product description BD Gravity Infusion Set Standardbore V/NV 20 DPML 1 BCV 1 RC 1 SmartSite™ Y 4WSC 2 T Access w BCV RLL (Removable Extension 1 PC 1 SmartSite™ Y 36 in)		US 
Part no. 42502E	GSD			
Length 312 cm (123 in) Priming vol. 26 mL Ctn. qty. 50 Material PVC/DEHP		Product description BD Gravity Infusion Set Standardbore V/NV 10 DPML 1 BCV 1 RC 1 SmartSite™ Y 4WSC w BCV RLL (2 Removable Extensions 2 T Port Extension cm/8 in and 1 SmartSite Y 1 PC 96 cm/38 in)		US  
Part no. 1001265	GSD			
Length 320 cm (126 in) Priming vol. 24 mL Ctn. qty. 50 Material PVC*		Product description BD Gravity Manifold Infusion Set Standardbore V/NV 10 DPML 1 BCV 1 RC 2 GANG 3WSC(Removable) RLL (Removable Extension 1 SmartSite™ Y 1 PC 84 cm/33 in)		US 
Part no. 42504E	GSD			
Length 295 cm (116 in) Priming vol. 24 mL Ctn. qty. 50 Material PVC*		Product description BD Gravity Manifold Infusion Set Standardbore V/NV 10 DPML 1 BCV 1 RC 1 SmartSite™ Y 3 GANG 4WSC(Removable) RLL (Removable Extension 1 SmartSite™ Y 1 PC 34 in/13 in)		US 





Part no. CS42522E-07	GSD	 	
Length 378 cm (149 in) Priming vol. 30 mL Ctn. qty. 20 Material PVC*			
Product description BD Gravity Manifold Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 1 SmartSite™ Y 3 Gang SmartSite™ 4WSC RLL (1 Removable Extension Set 1 PC 127 cm/50 in and 1 Removble Extension Set 1 SmartSite™ Y 1 SC 20 cm/8 in)		US	
Part no. 10050187	GSD	 	
Length 188 cm (74 in) Priming vol. 17 mL Ctn. qty. 20 Material PVC/DEHP			
Product description BD Gravity Manifold Infusion Set Standardbore V/NV 60 DPML 1 RC 2 SmartSite™ Y 4 GANG 4WSC RLL (1 Removable Extension Set 1 PC 46 cm/18 in and 1 Removable Extension 15 cm/6 in)		US	 
Part no. 42521E	GSD	 	
Length 346 cm (136 in) Priming vol. 28 mL Ctn. qty. 20 Material PVC*			
Product description BD Gravity Manifold Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 2 SmartSite™ Y 2 GANG 4WSC (Removable) RLL (Removable Extension 1 SmartSite™ Y 1 SC 58 cm/23 in)		US	
Part no. C42006E	GSD	 	
Length 351 cm (138 in) Priming vol. 26 mL Ctn. qty. 20 Material PVC/DEHP			
Product description BD Gravity Manifold Infusion Set Standardbore V/NV 10 DPML 1 BCV 1 RC 1 PC 2 SmartSite™ Y 2 GANG 4WSC RLL (Removable Extension 1 SmartSite™ Y 1 PC 58 cm/23 in)		US	 




Part no. 11059727	GSD	 
Length 259 cm (102 in) Priming vol. 21 mL Ctn. qty. 20 Material PVC/DEHP		
Product description BD Gravity Stopcock Infusion Set Standardbore V/NV 60 DPML 1 BCV 1 RC 2 SmartSite™ Y 4WSC RLL (1 Removable Extension Sets 1 SmartSite™ Y 89 cm/35 in and 1 Removable 1 SmartSite™ Y 1 SC 20 cm/8 in)		US  
Part no. 11059728	GSD	 
Length 259 cm (102 in) Priming vol. 21 mL Ctn. qty. 20 Material PVC/DEHP		
Product description BD Gravity Infusion Set Standardbore V/NV 20 DPML 1 BCV 1 RC 2 SmartSite™ Y 4WSC RLL (1 Removable Extension Sets 1 SmartSite™ Y 65 cm/25 in and 1 Removable 1 SmartSite™ Y 1 SC 20 cm/8 in)		US  
Part no. C42014E-07	GSD	 
Length 297 cm (117 in) Priming vol. 24 mL Ctn. qty. 20 Material PVC*		
Product description BD Gravity Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 1 SC 3 SmartSite™ Y 4WSC(Removable) RLL (Removable Extension 1 SmartSite™ Y 1 SC 76 cm/30 in)		US 
Part no. C42010E	GSD	 
Length 360 cm (142 in) Priming vol. 25 mL Ctn. qty. 20 Material PVC/DEHP		
Product description BD Gravity Infusion Set Standardbore V/NV 10 DPML 1 BCV 1 RC 2 SmartSite™ Y 4WSC RLL (Removable Extension 1 SmartSite™ Y 1 PC 84 cm/33 in)		US  




Part no. 42100E-07	GSD	 
Length 332 cm (131 in) Priming vol. 26 mL Ctn. qty. 20 Material PVC*		
Product description BD Gravity Manifold Infusion Set Standardbore V/NV 10 DPML 1 BCV 1 RC 1 SC 2 SmartSite™ Y 2 GANG 3WSC(Removable) RLL (1 Removable Extension Set 2 SmartSite™ Y 1 SC 97 cm/38 in and 1 Removable Set 1 SC 15 cm/6 in)		US 


Part no. MX4452	GSD	 
Length 355 cm (140 in) Priming vol. Ctn. qty. 25 Material PVC*		
Product description BD Gravity Manifold Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 2 NAC-Y 2 GANG 4WSC(Removable) RLL (1 Removable Extension 1 SC 1 NAC-Y 135 cm/53 in)		US 


Part no. MX4465	GSD	 
Length 265 cm (106 in) Priming vol. Ctn. qty. 25 Material PVC*		
Product description BD Gravity Stopcock Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 1 NAC-Y 2 GANG 4WSC (Removable) RLL (1 Removable Extension Set 1 BCV 56 cm/22 in and 1 Removable Extension Set 1 BCV 1 SC 1 NAC-Y 18 cm/7 in)		US 

BD Primary Gravity Manifold Sets



Gravity Manifold Sets		
Part no. 42519E	GSD	 
Length 300 cm (118 in) Priming vol. 23 mL Ctn. qty. 20 Material PVC/DEHP		
Product description BD Gravity Manifold Infusion Set Standardbore V/NV 20 DPML 1 BCV 1 RC 1 PC 2 SmartSite™ Y 3 Port Manifold w BCV (Removable) RLL		US  




Part no. CM42500E-07	GSD	 15 ml	 2.7mm
Length 335 cm (132 in) Priming vol. 28 mL Ctn. qty. 20 Material PVC*			
Product description		US	
BD Gravity Manifold Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 2 SmartSite™ Y 3 Port SmartSite™ Manifold (Removable) RLL (Removable Extension 1 PC 1 SmartSite™ Y 40 cm/16 in)			


Part no. MX4436	GSD	 15 ml	 2.7mm
Length 421 cm (166 in) Priming vol. 34 mL Ctn. qty. 25 Material PVC*			
Product description		US	
BD Gravity Manifold Infusion Set Standardbore V/NV 15 DPML 1 BCV 1 RC 2 NAC-Y 3 Port Manifold (Removable) 4WSC (Removable) RLL (Removable Extension 1 SC 1 NAC-Y 43 in/110 cm)			

BD Primary Gravity Blood Filter Sets

Gravity Blood Filter Sets			
Part no. 10387509	GSD	 20 ml	 2.7mm
Length 267 cm (105 in) Priming vol. 27 mL Ctn. qty. 50 Material PVC*			
Product description		US	
BD Gravity Blood Infusion Set Standardbore NV 20 DPML 200um 1 RC RLL			
Part no. 42081E	GSD	 15 ml	 2.7mm
Length 264 cm (104 in) Priming vol. 30 mL Ctn. qty. 20 Material PVC/DEHP			
Product description		US	
BD Gravity Blood Infusion Set Dual Spike Standardbore NV 15 DPML 180um 3 RC 1 SmartSite™ Y RLL			



Part no. 10010903		GSD	 15 ml	 2.7mm
Length	264 cm (104 in)			
Priming vol.	67 mL			
Ctn. qty.	20			
Material	PVC*			
		Product description BD Gravity Blood Infusion Set Dual Spike Standardbore NV 15 DPML 180um 3 RC 1 SmartSite™ Y Blood Handpump RLL	US	CE
				



Secondary Administration Sets





BD Secondary Gravity Sets

Secondary Infusion Set				
Part no. 10015501	GSD		 0.5mm	
Length 70 cm (28 in) Priming vol. 1.3 mL Ctn. qty. 50 Material PVC*				
	Product description BD Secondary Syringe Adapter Infusion Set Smallbore V/NV 1 PC RLL	US		
Part no. 72215N	GSD		 15 ml	 2.7mm
Length 102 cm (40 in) Priming vol. 13 mL Ctn. qty. 100 Material PVC*				
	Product description BD Secondary Infusion Set Standardbore V/NV 15 DPML 1 RC RLL	US CE		
Part no. MT6000	GSD		 2.7mm	
Length 56 cm (22 in) Priming vol. Ctn. qty. 50 Material PVC*				
	Product description BD Secondary Infusion Set (Transfer Set) Standardbore V/NV 1 RC Luer Slip	US		
Part no. 10016073	GSD		 20 ml	 2.7mm
Length 79 cm (31 in) Priming vol. 11 mL Ctn. qty. 100 Material PVC*				
	Product description BD Secondary Infusion Set Standardbore V/NV 20 DPML 1 RC RLL	US CE		
Part no. 11610341	GSD		 20 ml	 2.7mm
Length 79 cm (31 in) Priming vol. 12 mL Ctn. qty. 100 Material PVC*				
	Product description BD Secondary Infusion Set Standardbore NF Container Access 20 DPML 1 RC FLL	US		



Part no. 10014881	GSD			
Length 80 cm (31 in) Priming vol. 12 mL Ctn. qty. 100 Material PE Lined PVC*				
	Product description BD Secondary Infusion Set Standardbore NF Container Access 20 DPML 1 RC FLL	US		
Part no. C70101	GSD			
Length 92 cm (36 in) Priming vol. 12 mL Ctn. qty. 30 Material PVC*				
	Product description BD Secondary Infusion Set Standardbore V/NV 20 DPML 15um 1 RC RLL	US		

BD Secondary Texium™ Sets

Secondary Infusion Set (Texium™)				
Part no. 70001B-07T	GSD			
Length 80 cm (31 in) Priming vol. 12 mL Ctn. qty. 100 Material PE Lined PVC/DEHP				
	Product description BD Secondary Infusion Set Standardbore NF Container Access 20 DPML 1 RC Texium™	US		
Part no. 40000-07T	GSD			
Length 81 cm (32 in) Priming vol. 12 mL Ctn. qty. 100 Material PVC*				
	Product description BD Secondary Infusion Set Standardbore V/NV 20 DPML 1 SC Texium™	US		
Part no. 4030B-07T	GSD			
Length 81 cm (32 in) Priming vol. 14 mL Ctn. qty. 100 Material PE lined PVC/DEHP				
	Product description BD Secondary Infusion Set Standardbore NF Container Access 20 DPML 0.2um Large PES Filter 1 RC 1 PC Texium™	US		



BD Secondary Transfer Set

Secondary Transfer Set				
Part no. 2309E	GSD		 2.7mm	
Length 13 cm (5 in)				
Priming vol. 0.93 mL				
Ctn. qty. 100	Product description BD Container Access Device Standardore NV 1BCV 1 Smartsite™ Bonded		US	
Material PVC*				

BD Secondary Blood Infusion Set

Secondary Transfusion Set				
Part no. 10016242	GSD		 20 ml	 2.7mm
Length 81 cm (32 in)				
Priming vol. 16 mL				
Ctn. qty. 50	Product description BD Secondary Blood Infusion Set Standardbore NV 20 DPML 200um 1 RC RLL		US CE	
Material PVC*				

BD Alaris™ Infusion System - Infusion Sets



Sales and Customer Service


 **1-844-8BD LIFE (1-844-823-5433)**

Clinical Support

 **1-800-555-7422 Ext 5**

 **medical.information@bd.com**












BD Alaris™ Pump Infusion Sets




BD Alaris™ Pump Infusion Sets


BD Alaris™ Pump Infusion Sets			
Part no. 2204-0007	MMS	 20 ml	 2.7mm
Length 119 cm (47 in) Priming vol. 14 mL Ctn. qty. 20 Material PVC/DEHP			
Product description		US	 
BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 RC RLL			
Part no. 10942011	MMS	 20 ml	 2.7mm
Length 277 cm (109 in) Priming vol. 23 mL Ctn. qty. 20 Material PVC*			
Product description		US	
BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 RC RLL			
Part no. 10863358	MMS	 15 µm	 20 ml
Length 296 cm (117 in) Priming vol. 25 mL Ctn. qty. 20 Material PVC/DEHP			
Product description		US	 
BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML 15um Alaris™ Pump Segment 1 RC RLL			
BD Alaris™ Pump Infusion Sets with Smartsite™ Y-sites and Back Check Valves			
Part no. 24001-0007	MMS	 20 ml	 2.7mm
Length 120 cm (47 in) Priming vol. 14 mL Ctn. qty. 20 Material PVC/DEHP			
Product description		US	 
BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 1 SmartSite™ Y RLL			




Part no. 2420-0007		MMS	 20 ml	 2.7mm	
Length	297 cm (117 in)				
Priming vol.	25 mL				
Ctn. qty.	20				
Material	PVC*				
Product description		BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 2 SmartSite™ Y RLL	US	CE	
Part no. 10013186		MMS	 60 ml	 2.7mm	
Length	297 cm (117 in)				
Priming vol.	25 mL				
Ctn. qty.	20				
Material	PVC/DEHP				
Product description		BD Alaris™ Pump Infusion Set Standardbore V/NV 60 DPML Alaris™ Pump Segment 1 BCV 2 RC 3 SmartSite™ Y RLL	US		
Part no. 2426-0007		MMS	 20 ml	 2.7mm	
Length	320 cm (126 in)				
Priming vol.	26 mL				
Ctn. qty.	20				
Material	PVC*				
Product description		BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 1 PC 3 SmartSite™ Y RLL	US		
Part no. 11426965		MMS	 20 ml	 2.7mm	
Length	349 cm (137 in)				
Priming vol.	28 mL				
Ctn. qty.	20				
Material	PVC*				
Product description		BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 1 SC 5 SmartSite™ Y RLL	US		
Part no. 2452-0007		MMS	 10 ml	 2.7mm	
Length	282 cm (111 in)				
Priming vol.	25 mL				
Ctn. qty.	20				
Material	PVC*				
Product description		BD Alaris™ Pump Infusion Set Standardbore V/NV 10 DPML Alaris™ Pump Segment 2 BCV 1 RC 3 SmartSite™ Y RLL	US		




Part no. 11171447	MMS	 20 ml	 2.7mm
Length 320 cm (126 in) Priming vol. 26 mL Ctn. qty. 20 Material PVC/DEHP			
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 2 BCV 1 RC 1 PC 1 SC 3 SmartSite™ Y RLL		US	 

Part no. 10561554	MMS	 15 μm	 20 ml	 2.7mm
Length 321 cm (126 in) Priming vol. 26 mL Ctn. qty. 20 Material PVC/DEHP				
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML 15um Alaris™ Pump Segment 1 BCV 1 RC 1 PC 3 SmartSite™ Y RLL		US	 	


BD Alaris™ Pump PE Lined PVC Infusion Sets


BD Alaris™ Pump PE Lined PVC Infusion Sets with no Y-sites				
Part no. 11426964	MMS	 20 ml	 2.7mm	
Length 118 cm (47 in) Priming vol. 15 mL Ctn. qty. 20 Material PE lined PVC/DEHP				
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 RC RLL		US	 	
Part no. 2260-0500	MMS	 20 ml	 2.7mm	
Length 286 cm (113 in) Priming vol. 23 mL Ctn. qty. 20 Material PE lined PVC/DEHP				
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 RC RLL		US	 	




BD Alaris™ Pump PE Lined PVC Infusion Sets with Smartsite™ Y-sites and Back Check Valves			
Part no. 24600-0007	MMS	 20 ml	 2.7mm
Length 304 cm (120 in) Priming vol. 25 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 3 SmartSite™ Y RLL		US	 

BD Alaris™ Pump 0.2µm Filter Infusion Sets


BD Alaris™ Pump Filter Infusion Set with 0.2µm Inline Filter, Smartsite™ Y-sites and Back Check Valves			
Part no. 2434-0007	MMS	 20 ml	 0.2 µm
		 2.7mm	
Length 293 cm (115 in) Priming vol. 27 mL Ctn. qty. 20 Material PVC*			
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 0.2µm Large PES Filter 1 RC 1 SC 1 SmartSite™ Y RLL		US	

Part no. 2432-0007	MMS	 20 ml	 0.2 µm
		 2.7mm	
Length 299 cm (118 in) Priming vol. 28 mL Ctn. qty. 20 Material PVC*			
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 0.2µm Large PES Filter 1 RC 1 SC 3 SmartSite™ Y RLL		US	


BD Alaris™ Pump PE Lined PVC 0.2µm Filter Infusion Sets

BD Alaris™ Pump PE Lined PVC Filter Infusion Set with 0.2µm Inline Filter			
Part no. 10010454	MMS	 20 ml	 0.2 µm
		 2.7mm	
Length 291 cm (115 in) Priming vol. 26 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 0.2µm Large PES Filter 1 RC 1 PC 1 SmartSite™ Y RLL		US	 




BD Alaris™ Pump PE Lined PVC Filter Infusion Set with 0.2um Inline Filter, Smartsite™ Y-sites and Back Check Valves			
Part no. 11532269	MMS	  	
Length 318 cm (125 in) Priming vol. 29 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
	Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 0.2um Large PES Filter 1 RC 1 PC 2 SmartSite™ Y RLL	US 	

BD Alaris™ Pump 1.2µm Filter Infusion Sets


BD Alaris™ Pump Filter Infusion Set with 1.2um Inline Filter			
Part no. 2202-0007	MMS	  	
Length 309 cm (122 in) Priming vol. 27 mL Ctn. qty. 20 Material PVC/DEHP			
	Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1.2um Large PES Filter 1 RC RLL	US 	

BD Alaris™ Pump PE Lined PVC 1.2µm Filter Infusion Sets


BD Alaris™ Pump PE Lined PVC Filter Infusion Set with 1.2um Inline Filter and SmartSite™ Y-sites			
Part no. 10010453	MMS	  	
Length 291 cm (115 in) Priming vol. 26 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
	Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1.2um Large PES Filter 1 RC 1 SC 1 SmartSite™ Y RLL	US 	




BD Alaris™ Pump Burette Infusion Sets

BD Alaris™ Pump Burette Infusion Set with SmartSite™ Y-sites			
Part no. 2441-0007	MMS	 	
Length 330 cm (130 in) Priming vol. 20 mL Ctn. qty. 10 Material PVC/DEHP			
Product description BD Alaris™ Pump Burette Infusion Set Standardbore V/NV 60 DPML Alaris™ Pump Segment 2 RC 1 SC 3 SmartSite™ Y Graduated 150 mL (Top Plate V/NV 1 NF Access) RLL		   	





BD Alaris™ Pump PE Lined PVC Burette 0.2µm Filter Infusion Sets

BD Alaris™ Pump PE Lined PVC Burette Infusion Set with 0.2µm Inline Filter and SmartSite™ Y-sites			
Part no. 10015012	MMS	  	
Length 317 cm (125 in) Priming vol. 26 mL Ctn. qty. 10 Material PE lined PVC/DEHP			
Product description BD Alaris™ Pump Burette Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 0.2µm Large PES Filter 2 RC 2 PC 1 SmartSite™ Y Graduated 150 mL (Top Plate V/NV 1 NF Access) RLL		  	

BD Alaris™ Pump Stopcock/Manifold Infusion Sets

BD Alaris™ Pump Stopcock/Manifold Infusion Set with SmartSite™ Y-sites and Back Check Valves			
Part no. 10015896	MMS	 	
Length 318 cm (125 in) Priming vol. 22 mL Ctn. qty. 10 Material PVC*			
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 10 DPML Alaris™ Pump Segment 1 BCV 1 RC 1 SC 2 SmartSite™ Y 2 GANG 4WSC (Removable) RLL (Removable Extension 1 SmartSite™ Y 1 SC Set 76 cm/30 in)		 	



Part no. 2423-0007	MMS	 
Length 345 cm (136 in) Priming vol. 28 mL Ctn. qty. 10 Material PVC*		
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 10 DPML Alaris™ Pump Segment 1 BCV 1 RC 1 SC 2 SmartSite™ Y 2 GANG 4WSC (Removable) RLL (Removable Extension Set 2 SmartSite™ Y 1 SC 102 cm/40 in)		US
		
Part no. 2422-0007	MMS	 
Length 401 cm (158 in) Priming vol. 31 mL Ctn. qty. 20 Material PVC*		
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 2 SmartSite™ Y 4WSC (Removable) RLL (Removable Extension Set 2 SmartSite™ Y 2 PC 97 cm/38 in)		US
		
Part no. 11419365	MMS	 
Length 403 cm (159 in) Priming vol. 31 mL Ctn. qty. 20 Material PVC*		
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 1 PC 3 SmartSite™ Y Manifold and 3WSC (Removable) RLL (Removable Extension Set 5 SmartSite™ Y 1 PC 1 SC 69 cm/27 in)		US
		
Part no. 10813621	MMS	 
Length 404 cm (159 in) Priming vol. 31 mL Ctn. qty. 20 Material PVC*		
Product description BD Alaris™ Pump Infusion Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 1 SC 3 SmartSite™ Y 3 GANG 4WSC (Removable) RLL (Removable Extension Set 1 SmartSite™ Y 71 cm/28 in)		US
		



BD Alaris™ Pump Blood Infusion Sets

BD Alaris™ Pump Blood Single Spike Infusion Sets			
Part no. 2278-0500	MMS		
Length 307 cm (121 in) Priming vol. 33 mL Ctn. qty. 20 Material PVC/DEHP			
Product description BD Alaris™ Pump Blood Infusion Set Standardbore NV 20 DPML 200um Alaris™ Pump Segment 1 RC RLL			

BD Alaris™ Pump Blood Dual Spike Infusion Sets

BD Alaris™ Pump Blood Dual Spike Infusion Sets			
Part no. 10062818	MMS		
Length 308 cm (121 in) Priming vol. 33 mL Ctn. qty. 10 Material PVC/DEHP			
Product description BD Alaris™ Pump Blood Infusion Set Standardbore NV 15 DPML 180um Alaris™ Pump Segment 1 RC 2 SC RLL			

BD Alaris™ Pump Blood Dual Spike Infusion Sets with SmartSite™ Y-sites			
Part no. 2477-0007	MMS		
Length 308 cm (122 in) Priming vol. 33 mL Ctn. qty. 10 Material PVC*			
Product description BD Alaris™ Pump Blood Infusion Set Standardbore NV 15 DPML 180um Alaris™ Pump Segment 1 RC 2 SC 1 SmartSite™ Y RLL			





BD Alaris™ Pump Infusion Sets with SmartSite™ Bag Access Port and Texium™ Closed Male Luer

BD Alaris™ Pump Infusion Sets with SmartSite™ Bag Access Port and Texium™ Closed Male Luer			
Part no. 22000-B007T	MMS	 20 ml	 2.7mm
Length 115 cm (45 in) Priming vol. 14 mL Ctn. qty. 20 Material PVC/DEHP			
	Product description BD Alaris™ Pump Set Standardbore NF Container Access 20 DPML Alaris™ Pump Segment 1 RC Texium™	US	 

BD Alaris™ Pump Infusion Sets PE Lined PVC with SmartSite™ Bag Access Port and Texium™ Closed Male Luer


BD Alaris™ Pump PE Lined PVC Infusion Sets with SmartSite™ Bag Access Port and Texium™ Closed Male Luer			
Part no. 22603-B007T	MMS	 20 ml	 2.7mm
Length 120 cm (47 in) Priming vol. 14 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
	Product description BD Alaris™ Pump Set Standardbore NF Container Access 20 DPML Alaris™ Pump Segment 1 RC Texium™	US	 

Part no. 22602-B007T	MMS	 20 ml	 2.7mm
Length 300 cm (118 in) Priming vol. 24 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
	Product description BD Alaris™ Pump Set Standardbore NF Container Access 20 DPML Alaris™ Pump Segment 1 RC Texium™	US	 

Part no. 24601-B007T	MMS	 20 ml	 2.7mm
Length 300 cm (118 in) Priming vol. 24 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
	Product description BD Alaris™ Pump Set Standardbore NF Container Access 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 2 SmartSite™ Y Texium™	US	 




BD Alaris™ Pump Infusion Sets with Texium™ Closed Male Luer

BD Alaris™ Pump Infusion Sets with Texium™ Closed Male Luer			
Part no. 24010-0007T	MMS	 20 ml	 2.7mm
Length 321 cm (126 in) Priming vol. 26 mL Ctn. qty. 20 Material PVC*			
Product description BD Alaris™ Pump Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 1 PC 3 SmartSite™ Y Texium™		US	

BD Alaris™ Pump Infusion Sets PE Lined PVC with Texium™ Closed Male Luer



BD Alaris™ Pump PE Lined PVC Infusion Sets with Texium™ Closed Male Luer			
Part no. 22600-0007T	MMS	 20 ml	 2.7mm
Length 120 cm (47 in) Priming vol. 14 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
Product description BD Alaris™ Pump Set Standardbore V/NV 20 DPML Alaris™ Pump Segment 1 RC Texium™		US	 

BD Alaris™ Pump Infusion Sets with SmartSite™ Bag Access Port


BD Alaris™ Pump Infusion Sets with SmartSite™ Bag Access Port			
Part no. 10013361T	MMS	 20 ml	 2.7mm
Length 297 cm (117 in) Priming vol. 25 mL Ctn. qty. 20 Material PVC/DEHP			
Product description BD Alaris™ Pump Set Standardbore NF Container Access 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 2 SmartSite™ Y Texium™		US	 



BD Alaris™ Pump Infusion Sets PE Lined PVC with SmartSite™ Bag Access Port


BD Alaris™ Pump PE Lined PVC Infusion Sets with SmartSite™ Bag Access Port			
Part no. 10015862	MMS	 20 ml	 2.7mm
Length 115 cm (45 in) Priming vol. 14 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
Product description BD Alaris™ Pump Set Standardbore NF Container Access 20 DPML Alaris™ Pump Segment 1 RC RLL		US	 
Part no. 2466-0007	MMS	 20 ml	 2.7mm
Length 220 cm (87 in) Priming vol. 21 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
Product description BD Alaris™ Pump Set Standardbore NF Container Access 20 DPML Alaris™ Pump Segment 1 RC RLL		US	 
Part no. 10015861A	MMS	 20 ml	 2.7mm
Length 323 cm (127 in) Priming vol. 25 mL Ctn. qty. 20 Material PE lined PVC/DEHP			
Product description BD Alaris™ Pump Set Standardbore NF Container Access 20 DPML Alaris™ Pump Segment 1 BCV 1 RC 1 PC 3 SmartSite™ Y RLL		US	 

BD Alaris™ Pump Infusion Sets PE Lined PVC Filter Infusion Sets with SmartSite™ Bag Access Port and Texium™ Closed Male Luer

BD Alaris™ Pump PE Lined PVC Filter Infusion Sets with 0.2um Inline Filter, SmartSite™ Bag Access Port, and Texium™ Closed Male Luer			
Part no. 24301-0007T	MMS	 20 ml	 0.2 μm
Length 292 cm (115 in) Priming vol. 27 mL Ctn. qty. 20 Material PE lined PVC/DEHP			 2.7mm
Product description BD Alaris™ Pump Set Standardbore NF Container Access 20 DPML Alaris™ Pump Segment 0.2um Large PES Filter 1 RC 1 PC 1 SmartSite™ Y Texium™		US	 




BD Alaris™ Pump Infusion Sets PE Lined PVC Filter Infusion Sets with SmartSite™ Bag Access Port

BD Alaris™ Pump PE Lined PVC Filter Infusion Sets with 0.2um Inline Filter and SmartSite™ Bag Access Port			
Part no. 2465-0007	MMS	 20 ml	 0.2 μm
Length 294 cm (116 in)			 2.7mm
Priming vol. 27 mL	Product description BD Alaris™ Pump Set Standardbore NF Container Access 20 DPML Alaris™ Pump Segment 0.2um Large PES Filter 1 RC 1 PC 1 SmartSite™ Y RLL		  
Ctn. qty. 20			
Material PE lined PVC/DEHP			

BD Alaris™ Pump Syringe Adapter Infusion Sets

BD Alaris™ Pump Syringe Adapter Infusion Sets			
Part no. 10010483	MMS		 1mm
Length 208 cm (82 in)			
Priming vol. 4 mL	Product description BD Alaris™ Pump Syringe Adapter Set Smallbore V/NV Alaris™ Pump Segment 1 SC RLL		  
Ctn. qty. 20			
Material PVC/DEHP			

BD Alaris™ Pump Infusion Half Sets

BD Alaris™ Pump Infusion Half Sets with SmartSite™ Y-sites			
Part no. 2403-0007	MMS		 2.7mm
Length 220 cm (87 in)			
Priming vol. 13 mL	Product description BD Alaris™ Pump Half Set Standardbore Alaris™ Pump Segment 1 RC 1 SmartSite™ Y RLL		  
Ctn. qty. 20			
Material PVC/DEHP			



BD Alaris™ Syringe Infusion Sets



BD Alaris™ Syringe Infusion Sets with Pressure Sensing Disc

BD Alaris™ Syringe Infusion Sets with Pressure Sensing Disc				
Part no. 10014914	MMS			 0.5mm
Length 152 cm (60 in) Priming vol. 0.7 mL Ctn. qty. 100 Material PVC*				
Product description BD Alaris™ Syringe Infusion Set Microbore PSD 1 SC FLL			US	
Part no. 10014917	MMS			 1mm
Length 152 cm (60 in) Priming vol. 1.7 mL Ctn. qty. 100 Material PVC*				
Product description BD Alaris™ Syringe Infusion Set Smallbore PSD 1 SC FLL			US	

BD Alaris™ Syringe PE Lined PVC Infusion Sets with Pressure Sensing Disc

BD Alaris™ Syringe PE Lined PVC Infusion Sets with Pressure Sensing Disc				
Part no. 10014918	MMS			 0.5mm
Length 152 cm (60 in) Priming vol. 1.2 mL Ctn. qty. 100 Material PE lined PVC*				
Product description BD Alaris™ Syringe Infusion Set Microbore PSD 1 PC FLL			US	

BD Alaris™ Syringe Filter Infusion Sets with Pressure Sensing Disc

BD Alaris™ Syringe Filter Infusion Sets with 0.2um Inline Filter and Pressure Sensing Disc				
Part no. 10014916	MMS		 0.2 μm	 0.5mm
Length 152 cm (60 in) Priming vol. 1 mL Ctn. qty. 100 Material PVC*				
Product description BD Alaris™ Syringe Infusion Set Microbore PSD 0.2um Small Volume PES Filter 1 SC FLL			US	




BD Alaris™ Syringe Infusion Sets

BD Alaris™ Syringe Infusion Sets with Anti-Siphon Valve, Back Check Valve and Pressure Sensing Disc				
Part no. 10798703	MMS		20 ml	1mm
Length 244 cm (96 in) Priming vol. 11 mL Ctn. qty. 50 Material PVC*				
	Product description BD Alaris™ Syringe Extension Set Smallbore V/NV 20 DPML PSD 1 BCV 1 ASV Side Arm w NFC RLL		US	

BD Alaris™ Syringe Infusion Set

BD Alaris™ Syringe Infusion Sets with no PSD				
Part no. 10689602	MMS		20 ml	0.5mm
Length 241 cm (95 in) Priming vol. 13 mL Ctn. qty. 50 Material PVC*				
	Product description BD Alaris™ Syringe Extension Set Smallbore 1 SmartSite™ Bonded V/ NV 20 DPML 1 BCV 1 ASV Side Port w SmartSite™ FLL		US	




BD Alaris™ PCA Infusion Sets



Alaris™ PCA Infusion Set Kit

Alaris™ PCA Infusion Sets with Anti-Siphon Valves and Monoject™ Plunger Rod				
Part no. 10800176	MMS		 0.5mm	
Length 234 cm (92 in) Priming vol. 2.1 mL Ctn. qty. 50 Material PVC*				
Product description			US	
BD Alaris™ PCA Extension Set PVC Microbore 1 ASV 1 PC 1 FL w ASV FLL (w Plunger Rod)				

Alaris™ PCA Infusion Set with Open Y-site

Alaris™ PCA Infusion Sets with Open Y-site, Back Check Valve, and Anti-Siphon Valve				
Part no. 30852	MMS		 1mm	
Length 177 cm (69.5 in) Priming vol. 1.8 mL Ctn. qty. 100 Material PVC*				
Product description			US	
BD Alaris™ PCA Extension Set Smallbore 1BCV 1ASV 2 PC Sidearm Access FLL				

Part no. 30873	MMS		 0.5mm	
Length 229 cm (90 in) Priming vol. 3 mL Ctn. qty. 50 Material PVC*				
Product description			US	
BD Alaris™ PCA Extension Set Microbore 1BCV 1ASV 2 PC Sidearm Access FLL				

Alaris™ PCA Infusion Sets with Open Y-site, Back Check Valve, and Anti-Siphon Valve and Monoject™ Plunger Rod				
Part no. 10800175	MMS		 0.5mm	
Length 229 cm (90 in) Priming vol. 3 mL Ctn. qty. 50 Material PVC*				
Product description			US	
BD Alaris™ PCA Extension Set Microbore 1BCV 1ASV 2 PC Sidearm Access PLUNGER ROD FLL				

Glossary

Add-in/Add-on device

Refers to an infusion extension set or specialty component that can be manually added to any infusion set. Add-in/add-on devices are used to add functionality required for a specific clinical situation.¹

Burette/Buretrol

A burette is an inline metered chamber designed to support fluid volume management and allow for compounding of intermittent drug delivery.¹

Change interval

Infusion sets should be replaced according to CDC, National Guidelines, Hospital protocols and drug manufacturers recommendations.¹

Di (2-ethylhexyl) phthalate (DEHP)

A plasticizer used in some polyvinyl chloride medical devices is used to yield clarity and soft pliable tubing. DEHP has been shown to produce a range of adverse effects in experimental animals, notably liver toxicity, and testicular atrophy. Although the toxic and carcinogenic effects of DEHP have been well established in experimental animals, the ability of this compound to produce adverse effects in humans is controversial. There is no evidence that neonates, infants, pregnant and breast-feeding women exposed to DEHP experience any related adverse effects. However, a lack of evidence of causation between DEHP-PVC and any disease or adverse effect does not mean that there are no risks. Many manufacturers utilize alternate phthalates that have not been shown to cause harm.^{1,2}

Drip chamber (DC)

Drip chambers have two primary functions.

- Allows an estimate of the rate at which fluid is administered. Flow rate can be calculated with the help of the observations from the drip chamber drop former and its drop factor.
- Acts as an air trap. Allows gas (such as air) to rise out from a fluid so that it is not passed downstream.¹

DP/mL: Drops per mL or drop factor

The number of drops delivered equals 1 mL. Required to manually calculate flow rates. Common drop factors are¹



10 DP/mL



15 DP/mL



20 DP/mL



60 DP/mL

Auto-vent (AV)

Characteristic of the spike that enters the solution container. Designed for collapsible and non-collapsible containers. The auto vent automates venting without clinical intervention, allowing air to enter the container to allow fluid to flow, in addition it filters air, aiding in blocking micro-organism ingress.¹

Non-vented (NV)

Characteristic of the spike that enters the solution container. Designed for collapsible plastic containers as no air is required to enter the container for fluid to flow.¹

Vented (V/NV)

Characteristic of the spike that enters the solution container. Designed for containers that do not collapse when emptying (glass or semi-rigid plastic). Vent allows air to enter the container to allow fluid to flow, in addition filters air, aiding in blocking micro-organism ingress.¹

Extension sets

A collection of tubing, connectors, and components intended to establish an extension of tubing where the standard length and or components of the tubing in an intravenous (infusion) administration set is insufficient.¹

Filters

0.2-micron inline filter—an inline filter that removes particulate matter larger than 0.2 microns.¹

1.2-micron inline filter—an inline filter that removes particulate matter larger than 1.2 microns.¹

Drip chamber 15-micron screen—a screen in the base of designated drip chambers.

Supor / Polyether sulfone (PES) filter membrane—filter membranes which have low protein and drug binding characteristics.

Blood filter—filters containing a large screen to remove coarse fibrin and byproducts of stored blood. They range in size from 170–260um.¹

Glossary (continued)

Flow controller (FC)

A calibrated and scaled device used to regulate the flow of fluid.¹

Hand pump (gravity transfusion sets only)

A cylindrical chamber that when filled allows for infusate to be manually pumped (hand squeezed). Used for situations where rapid infusion is required.¹

Internal diameter (ID)

Internal lumen size of the infusion set tubing. The tubing ID dictates the calculated infusion set priming volume.¹



0.5 mm

Microbore tubing 0.5 mm internal diameter



1 mm

Smallbore tubing 1 mm internal diameter



2.1 mm

Standardbore tubing 2.1 mm internal diameter



3 mm

Macrobore tubing 3 mm internal diameter

In-line

Refers to a specialty component or feature that is an integral part of the infusion set.¹

Light protective administration set

Specialized tubing and components designed to provide protection from ultraviolet light when administering medications that degrade if exposed to light.^{1,3}

Luer connections (male/female)

Universal male and female connections (fittings) that make a correct fit for standard intravenous sets to be connected.¹

Manifolds

Allows for multiple infusions to flow into one tubing lumen.¹

Multiway adapters

Allows for multiple infusions to flow into one tubing lumen.¹

Needle-free Connectors (needle free-Y and needle free-T)

Designed to provide swabable, needle free, closed access to an infusion device.¹

Non-PVC

Tubing either manufactured completely without polyvinyl chloride or which has an inert lining that helps prevent exposure of infused fluids to polyvinyl chloride.¹

Non-latex

Natural rubber latex is not part of the material formulation.¹

One way valves

Anti siphon valve (ASV)

Aids in preventing fluid syphoning during container exchanges or line disconnections. Anti-siphon valves also prevent against back flow and free flow in the infusion line.¹

Back check valve (BCV)

A one-way valve that allows the solution to flow in one direction. BCVs are commonly found in infusion sets with needle-free Y access ports (located above a needle-free connector [NFC] to prevent back flow of medication during administration of an infusion push medication or secondary infusion).¹

Pinch and slide clamps (PC/SC)

Clamps designed to provide on/off fluid control.¹

Primary set

Typically, the main administration set used to deliver the infusate from the container to the patient.¹

Protective cap

Single-use component protective caps support contamination free handling during setup.

Prime cap/auto-prime cap single-use, tapered seal, male luer, cap with hydrophobic membrane designed to allow air to escape, but not fluid during initial set priming.¹

Priming volume (PV)

The amount of fluid required to fill the entire length of the infusion set eliminating the air in line.¹

Secondary set

A secondary administration set is defined as an infusion set attached to a primary administration set. Secondary sets are widely used to administer intermittent medications, such as antibiotics and are commonly referred to as a “piggyback.” (Primary intermittent sets are primary sets that are used for infusions that are disconnected and reconnected during use.)¹

Spike

Device used to access infusion containers.¹

Double spike infusion sets

Primary sets with two spikes. For use when rapid infusion or dual administration is required.¹

Stopcocks and stopcock gangs

A manually manipulated device used to direct the flow of solution. Additionally allows for multiple infusions to flow into one tubing lumen.¹

Roller clamp (RC)

Clamp allowing for manual regulation of flow. Tubing is incrementally occluded as the roller clamp is tightened.¹

Needle-free access ports/Needle-free connectors

Serves as an access to the tubing for medication administration, located at various locations along an infusion set— swabable needle-free access ports align with infection prevention and needlestick injury reduction initiatives.¹

Transfusion/Blood administration set

Transfusion/Blood infusion sets are designed to accommodate the viscous properties of blood and allow for rapid transfusion (if clinically required). These sets contain a large screen filter (170 micron to 240 micron) to filter coarse fibrin and clots.¹

Y-Access port/s

Allows for additional infusions to flow into one tubing lumen. Y-access may be open or contain a needle-free access device depending on the tubing manufacturers offerings.¹

References

1. Infusion Nurses Society, Corrigan, A., Gorski, L., Hankins, J., Perucca, R., & Alexander, M. (2014). *Infusion Nursing: an Evidence-Based Approach*. Elsevier Health Sciences
2. Guidance for Industry Limiting the Use of Certain Phthalates as Excipients in CDER-Regulated Products. (2012). https://www.fda.gov/files/drugs/published/Limiting-the-Use-of-Certain-Phthalates-as-Excipients-in-CDER-Regulated-Products_pdf.pdf
3. Tønnesen, H. (2004). *Photostability of Drugs and Drug Formulations* (2nd ed.). CRC Press. <https://doi.org/10.1201/9780367801021>

Symbols

This symbols glossary can be found online and not all symbols are applicable to products in this catalog.

Graphic	Symbol title	Explanatory text and standard
Medical devices — Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: general requirements		
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. ISO 15223-1:2016
	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized. ISO 15223-1:2016
	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide. ISO 15223-1:2016
	Sterile	Indicates a medical device that has been subjected to a sterilization process. ISO 15223-1:2016
	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation. ISO 15223-1:2016
	Sterile fluid path	Indicates the presence of a sterile fluid path within the medical device in cases when other parts of the device, including the exterior, might not be supplied sterile. ISO 15223-1:2016
	Example, sterile fluid path sterilized by ethylene oxide	Indicates the medical device contains a sterile fluid path that has been sterilized using ethylene oxide. ISO 15223-1:2016
	Example, sterile fluid path sterilized using irradiation	Indicates the medical device contains a sterile fluid path that has been sterilized using irradiation. ISO 15223-1:2016
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223-1:2016
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified. ISO 15223-1:2016
	Serial number	Indicates the manufacturer's serial number so that a specific medical device can be identified. ISO 15223-1:2016
	Date of manufacture	Indicates the date when the medical device was manufactured. ISO 15223-1:2016
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. ISO 15223-1:2016
	Use-by-date	Indicates the date after which the medical device is not to be used. ISO 15223-1:2016
	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened. ISO 15223-1:2016

Symbols (continued)

This symbols glossary can be found online and not all symbols are applicable to products in this catalog.

Graphic	Symbol title	Explanatory text and standard
Medical devices — Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: general requirements		
	Contains or presence of natural rubber latex	Indicates the presence of dry natural rubber or dry natural rubber latex as a material of construction within the medical device or the packaging of a medical device. ISO 15223-1:2016
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed. ISO 15223-1:2016
	Keep dry	Indicates a medical device that needs to be protected from moisture. ISO 15223-1:2016
	Consult instructions for use	Indicates the need for the user to consult the instructions for use. ISO 15223-1:2016
	Nonpyrogenic	Indicates a medical device that is nonpyrogenic. ISO 15223-1:2016
	Drops per milliliter	Indicates the number of drops per milliliter. ISO 15223-1:2016
	Liquid filter with pore size	Indicates an infusion or transfusion system of the medical device that contains a filter of a particular nominal pore size. ISO 15223-1:2016
	Open here	Indicates tab for package opening ISO 7000
	Inner diameter	Indicates inner diameter
	MR Safe	An item which is marked MR safe poses no known hazards in all MR environments
	MR Conditional	An item with demonstrated safety in the MR environment within defined conditions ASTM. F2503
Symbol for use in the labeling of medical devices — Requirements for labeling of medical devices containing phthalates.		
	Contains or presence of phthalate: bis (2-ethylhexyl) phthalate (DEHP)	Indicates the medical device is manufactured with products containing the phthalate: bis (2-ethylhexyl) phthalate (DEHP). BS EN 15986:2011
CE marking of conformity or “CE marking” means a marking by which a manufacturer indicates that a device is in conformity with the applicable requirements set out in this regulation and other applicable union harmonisation legislation providing for its affixing		
	CE marking	Signifies European technical conformity. EU 2017-745 EU 2017-746
Graphical symbols for use on equipment		
	Country of manufacture	Indicates the country of manufacture. IEC 60417



Components details

Male luer locks			
	Fixed male luer lock (FLL)	Rotating male luer lock (RLL)	Rotating male luer lock bifuse
	Rotating male luer lock trifuse	Rotating male luer lock quadfuse	Male luer cap
	Fixed male luer lock w back check valve	Fixed male luer lock w anti syphon valve	Rotating male luer lock w back check valve
	Rotating male luer lock w anti syphon valve	Male luer w flash blub	

Female Luer Locks			
	Female luer lock winged (removable components)	Female luer lock no wings (bonded components)	Female luer lock w back check valve
	Female luer lock w anti syphon valve		





Needle-free connectors			
	SmartSite™ Needle-Free connector	MaxPlus™ Needle-Free connector	MaxZero™ Needle-Free connector
	Q-Syte™ Needle-Free connector	Texium™ with Priming Cap	




Check valves	
	Anti siphon valve (ASV)

Stopcocks		
	Stopcock (blue tap)	Stopcock (white 'T' tap)
	Stopcock injection port (white 'T' tap)	Stopcock smartsite (white 'T' tap)










Components details (continued)

Stopcocks gangs/Manifolds			
	3 port manifold		3 Port Manifold with Smartsite™
	Stopcock 2 gang		Stopcock 3 Gang
	3 Stopcock manifold Smartsite™		3 Stopcock manifold

T-connectors			
	T-connector w injection port / rotating luer		T-connector w SmartSite™ / rotating luer
			T-connector w SmartSite™ / luer slip




Y - access			
	SmartSite™-Y		Injection port-Y
			Nac-Y
	Multiport 2 SmartSite™-Y		Multiport 4 SmartSite™-Y
			T access port
	Open Y access port		



Clamps			
	Slide clamp		Pinch clamp
			Roller clamp


Filter sets			
	0.2um Small volume PES filter		1.2um Small volume PES filter
			0.2um Small volume endotoxin retentive filter (+ve)
	1.2um Small volume endotoxin retentive filter (+ve)		0.2um Medium volume PES filter
			1.2um Medium volume PES filter
	Large inline 0.2 Micron filter generic		Large inline 1.2 Micron filter generic









Components details (continued)


Tube fuse junction		
	2 - 1 Tube fuse junction	
		
		2 - 1 Y tube fuse junction

Flow controllers	
	Flow controller - single scale
	Flow controller - double scale

Alaris™ Pump Segment	GP pump segment
	
Alaris™ System Pump Segment	Gp pump segment







GW pump segment	Pressure sensing disc
	
GW pump segment	Pressure sensing disc (PSD) Alaris™ Syringe Pump


Spikes	
	Spike non-vented
	Spike autovent
	Spike blunt blood non-vented
	Spike vented / non-vented
	Spike needle-free container access


In-line drip chambers	
	In-line drip chamber 20DP
	In-line drip chamber 20DP w 15um filter
	In-line drip chamber 20DP w 200um filter





Components details (continued)

Drip chambers			
	Drip chamber vented/ non-vented 60DP		Drip chamber vented/ non-vented 60DP 15um filter
	Drip chamber vented/ non-vented 60DP ARD		Drip chamber vented/non-vented 10DP/15DP/20DP
	Drip chamber vented/non-vented 10DP/15DP/20DP 15um filter		Drip chamber non-vented 10DP/15DP/20DP 15um filter
	Drip chamber vented/non-vented 10DP/15DP/20DP ARD		Drip chamber autovent 20DP
	Drip chamber with nf container access 20DP		

Transfusion blood drip chambers			
	Blood drip chamber 180um/200um non-vented		Blood drip chamber 180um/200um non-vented ball valve shut-off

Burette sets	
	Burette/buretrol set

Blood hand pump	
	Blood hand pump

Hanger	
	Iv bag/ container hanger

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

[bd.com](https://www.bd.com)



BD, the BD Logo, Alaris, MaxPlus, MaxZero, Q-Syte, SmartSite, Texium, Connecta and NRFit are trademarks of Becton Dickinson and Company or its affiliates. © 2024 BD. All rights reserved. BD-116883 (08 / 24)