

# Prevue™ II Peripheral Vascular Access System Voice Recognition



## Voice Commands (if enabled)

- Prevue brightness up
  - Prevue brightness down
  - Prevue capture image
  - Prevue scan image
  - Prevue save image
  - Prevue change depth
  - Prevue volume up
  - Prevue volume down
- Prevue evaluate vein
  - Prevue continue (*To dismiss Cue™ magnetizer reminder screen*)
  - Prevue yes (*To respond to "Accept Evaluation?" popup during automatic evaluation*)
  - Prevue no (*To respond to "Accept Evaluation?" popup during automatic evaluation*)



**For example, you could say: "Prevue save image."**