

BD PhaSeal™ vs BD PhaSeal™ Optima

Connection technique for administration



The BD PhaSeal™ System is a closed system drug transfer device (CSTD) that has been demonstrated to prevent exposure to hazardous drugs, from drug preparation to IV administration.

Airtight and leakproof, and built on the legacy of the BD PhaSeal™ CSTD, the BD PhaSeal™ Optima System mechanically prohibits the transfer of environmental contaminants into the system and the escape of drug vapor concentrations outside the system, thereby minimizing individual and environmental exposure to drug vapor, aerosols and spills.

BD PhaSeal™ System



BD PhaSeal™ Optima System



Make sure you hear the audible click when engaging the Connector/Injector
Remember to check your connections