



BD PyxisTM MedBankTM MedPass

A **new medication management software** module displaying medication orders with due times over a 24 hour period- **helping nurses dispense medications to the right patient at the right time.**



BD Pyxis™ MedBank™ MedPass

Software features include:

Patient Grid – List of active patients

Medication Grid – List of active medication orders by patient

Timing Grid – 24-hour timing window for orders to be dispensed

Icons – Color-coded, interactive icons help aid the nurse during the dispensing process

Alerts – Help prevent removal of medications outside the scheduled due time

MedPass
5/18/2023 8:14:12 AM
Logged in User: Williams, Theresa

Resolve Discrepancy Acknowledge Discrepancy

Search - Patient Search - Medication ☒ Stat/Once

Unit Patients All Patients

South Wing - Cab 01

Thursday, May 18, 2023

Patient ID / Name / DOB	Medications	Available	0800	1000	1200	1400
4527853 Johnson, Alexander 11/02/2021 Allergies	Alprazolam 1MG Tab Dose: N/A Qty 1 Freq: Override Start: 5/18/2023 - 12:00 AM Practitioner: Last: 11/14/2022 - 3:50 PM Instructions					
169604 Johnson, Carl 05/05/2021 Allergies	Clarithromycin 500MG Tab Dose: 500 MG Freq: Two times a day @0900, 2000 Start: 9/9/2022 - 12:00 AM Practitioner: Wilson, Greg Last: 9/10/2022 - 8:16 PM Instructions					
955564 Johnson, Clarissa 09/05/1955 Allergies	Metformin 850MG Tab Dose: 850 MG Freq: Every 12 hours Start: 7/1/2022 - 8:00 AM Practitioner: Gregory, Mark Last: 5/18/2023 - 8:13 AM Instructions					
TEMP742385 Johnson, Dominique 06/14/1976 Allergies	Metoprolol Succ ER 100MG Tab Dose: 100 MG Freq: Every Day @1200 Start: 7/1/2022 - 12:00 AM Practitioner: Wilson, Greg Last: 9/10/2022 - 12:16 PM Instructions					
3445464 Johnson, John 06/14/1976 Allergies	Acetaminophen 500MG Tab Dose: 500 MG Freq: PRN Start: 7/1/2022 - 12:00 AM Practitioner: Smith, Samantha Last: 5/18/2023 - 8:13 AM Instructions					
169605 Johnson, Enrique 08/10/1993 Allergies	Hydrocortisone Cream Dose: 2 MG Freq: PRN Start: 5/16/2023 - 12:00 AM Practitioner: Smith, Samantha Last: NONE Instructions					

Admin Inventory Stocking Add Patient Return/Waste Load/Remove MedPass Tray Override Dispense Medications

Data displayed in this material is fictitious and not intended to represent actual users, patients, or sites and is meant for illustrative purposes only.



Helps dispense
the correct medication



Supports safe and
secure medication storage



Automates medication
inventory management



Helps support on-time
medication administration

MedPass is the latest addition to the cloud-enabled BD Pyxis™ MedBank™ portfolio

MedBank™ Platform Features :

- Enterprise reporting and program management
- Web-based interfaces with supported 3rd party systems
- Flexible BD Pyxis™ CUBIE™ pocket options for customizable solutions

MedPass is available on the BD Pyxis™ MedBank™ **Mini**, **Tower** and **MedFlex** cabinets. More options are available for additional inventory tracking and visibility.

Minimum Win10 OS is required



Contact us to get started

BD—Canada, Mississauga, ON, L5N 0B3

bd.com

BD, the BD Logo, MedBank and Pyxis are trademarks of Becton, Dickinson and Company or its affiliates.
© 2025 BD. All rights reserved. BD-146438 (0325)

