

Monitor Accessory Compatibility Card

Alaris™ EtCO₂ Module Model 8300



WARNING:

Use only approved Microstream™ breath sampling accessories with the Alaris™ EtCO₂ Module. All accessories are single-patient use only and latex-free. Using accessories other than those specified may cause improper module performance and result in inaccurate readings.

NOTE:

Not all sets are available in all countries.

Note: “H” denotes enhanced humidity control for longer term use.

Administration Set Description
Microstream Non-Intubated Patients, Long-Term Use
Adult Smart CapnoLine® H Plus: Pain Management, Surgical GCF, Critical Care
Pediatric Smart CapnoLine H: Pain Management, Surgical GCF, Critical Care
CapnoLine H: When nasal sampling is preferred
Microstream Non-Intubated Patients, Short-Term Use
Adult Smart CapnoLine Plus: Procedural Sedation, Pain Management, OR, EMS, ED, Rapid Response
Pediatric Smart CapnoLine: Procedural Sedation, EMS, ED, Rapid Response
Smart CapnoLine Guardian System: Procedural Sedation, including Upper Endoscopy, Bronchoscopy
O ₂ /CO ₂ Nasal FilterLine®: When nasal sampling is preferred
Nasal/NIV Line™: EMS, ED
Microstream Intubated Patients, Long-Term Use
Adult/Pediatric FilterLine H Set: Critical Care, humidified environments
Adult/Pediatric VitaLine™ H Set: Critical Care, high ambient humidity
FilterLine H Set for Infant/Neonates: Critical Care, humidified environments
VitaLine H Set for Infant/Neonates: Critical Care, humidity controlled incubators
Microstream Intubated Patients, Short-Term Use
FilterLine Set for Infant/Neonates: OR, EMS, ED, Rapid Response, Transport Team

DME 10000300399-00 PN P00000198

Alaris is a registered trademark of CareFusion or one of its affiliates. All other trademarks are the property of their respective owners. Smart CapnoLine, FilterLine and NIV-Line are registered trademarks of Oridion or one of its subsidiaries. © 2008, 2016 CareFusion or one of its affiliates. All rights reserved.