



# BD Pyxis™ MedBank™ MedPass

A new medication management software module displaying medication orders with due times over a 24 hour period- **helping nurses dispense medications to the right patient at the right time.**



# BD Pyxis™ MedBank™ MedPass

Software features include:

**Patient Grid** – List of active patients

**Medication Grid** – List of active medication orders by patient

**Timing Grid** – 24-hour timing window for orders to be dispensed

**Icons** – Color-coded, interactive icons help aid the nurse during the dispensing process

**Alerts** – Help prevent removal of medications outside the scheduled due time

The screenshot displays the MedPass interface. At the top, it shows the user is logged in as Theresa Williams on 5/18/2023 at 8:14:12 AM. The interface is divided into two main sections: a Patient Grid on the left and a Medication Grid on the right.

**Patient Grid:** Lists patients with their ID, Name, and DOB. The selected patient is Johnson, Enrique (ID: 169605, DOB: 08/10/1993).

**Medication Grid:** Shows medication orders for South Wing - Cab 01. The grid has columns for 'Available' and time slots (0800, 1000, 1200, 1400). Medications listed include Alprazolam 1MG Tab, Clarithromycin 500MG Tab, Metformin 850MG Tab, Metoprolol Succ ER 100MG Tab, Acetaminophen 500MG Tab, and Hydrocortisone Cream. Each medication entry includes dose, frequency, start/end dates, and practitioner information. Color-coded icons (plus sign, green checkmark, orange circle) are used to indicate medication status.

At the bottom of the interface, there are several action buttons: Admin, Inventory, Stocking, Add Patient, Return/Waste, Load/Remove MedPass Tray, Override, and Dispense Medications.

Data displayed in this material is fictitious and not intended to represent actual users, patients, or sites and is meant for illustrative purposes only.



Helps dispense the correct medication



Supports safe and secure medication storage



Automates medication inventory management



Helps support on-time medication administration

# MedPass is the latest addition to the cloud-enabled BD Pyxis™ MedBank™ portfolio

## MedBank™ Platform Features :

- Enterprise reporting and program management
- Web-based interfaces with supported 3<sup>rd</sup> party systems
- Flexible BD Pyxis™ CUBIE™ pocket options for customizable solutions

MedPass is available on the BD Pyxis™ MedBank™ **Mini** and BD Pyxis™ MedBank™ **Tower**. More options are available for additional inventory tracking and visibility.

Minimum Win10 OS is required



### Contact us to get started

Email: [bd.mms@bd.com](mailto:bd.mms@bd.com)

Manufactured by CareFusion 303, Inc., 10020 Pacific Mesa Blvd. San Diego, CA 92121 United States.

[bd.com](http://bd.com)

BD, the BD Logo, MedBank and Pyxis are trademarks of Becton, Dickinson and Company or its affiliates. © 2023 BD. All rights reserved. BD-105036 November 2023

