

# Security Patches: BD FocalPoint™ Review Station

April 2023

BD has identified patches from Microsoft® that have been identified as critical or security related up to January 2023. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

## Microsoft® patches

Patch name	Description	Patch ID	Notes
Servicing stack update for Windows 8.1	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	<a href="#">KB5018922</a>	This update must be installed prior to all other updates
Cumulative update for .NET Framework 3.5	Includes cumulative reliability improvements in .NET framework 3.5.	<a href="#">KB5020862</a>	
Cumulative update for .NET 4.6.2, 4.7, 4.7.1, 4.7.2	Includes cumulative reliability improvements in .NET framework 4.6.2, 4.7, 4.7.1, 4.7.2	<a href="#">KB5020868</a>	
Cumulative update for .NET 4.8	Includes cumulative reliability improvements in .NET framework 4.8.	<a href="#">KB5020878</a>	

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

**bd.com**

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.





## Microsoft® patches

Patch name	Description	Patch ID	Notes
Cumulative Security update for .NET Framework 3.5	Includes cumulative security-only improvements in .NET framework 3.5.	<a href="#">KB5020897</a>	
Cumulative Security update for .NET 4.6.2, 4.7, 4.7.1, 4.7.2	Includes cumulative security-only improvements in .NET framework 4.6.2, 4.7, 4.7.1, 4.7.2	<a href="#">KB5020899</a>	
Cumulative Security update for .NET 4.8	Includes cumulative security-only improvements in .NET framework 4.8.	<a href="#">KB5020902</a>	
Cumulative Security update for Windows 8.1	Fixes a memory leak in the Local Security Authority Subsystem Service (LSASS.exe) in Windows Domain Controllers.	<a href="#">KB5021296</a>	
Cumulative Security update for Windows 8.1	Resolves multiple security issues in Windows 8.1	<a href="#">KB5022346</a>	
Cumulative Security update for Windows 8.1	Resolves multiple security issues in Windows 8.1	<a href="#">KB5022352</a>	

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

**bd.com**

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.