

Security Patches: BD FocalPoint™ Large Lab Server

April 2023

BD has identified patches from Microsoft® that have been identified as critical or security related up to January 2023. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

Microsoft® patches

Patch name	Description	Patch ID	Notes
Servicing stack update for Windows Server 2012	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	KB5018922	This update must be installed prior to all other updates
Security update for .NET Framework 4.8	Includes cumulative reliability improvements in .NET framework 4.8.	KB5020608	
Security update for .NET Framework 3.5	Includes cumulative reliability improvements in .NET framework 3.5.	KB5020862	
Security update for .NET Framework 4.6.2, 4.7, 4.7.1 and 4.7.2	Includes cumulative reliability improvements in .NET framework 4.6.2, 4.7, 4.7.1 and 4.7.2.	KB5020868	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.





Microsoft® patches

Patch name	Description	Patch ID	Notes
Security update for .NET Framework 4.8	Includes cumulative security improvements in .NET framework 4.8.	KB5020878	
Security update for .NET Framework 3.5	Includes cumulative security improvements in .NET framework 3.5.	KB5020897	
Security update for .NET Framework 4.6.2, 4.7, 4.7.1 and 4.7.2	Includes cumulative reliability improvements in .NET framework 4.6.2, 4.7, 4.7.1 and 4.7.2.	KB5020899	
Security update for .NET Framework 4.8	Includes cumulative security improvements in .NET framework 4.8.	KB5020902	
Security update for Windows Server 2012 R2	Fixes a memory leak in the Local Security Authority Subsystem Service (LSASS.exe) for domain controllers.	KB5021296	
Security update for Windows Server 2012 R2	Fixes a memory leak in the Local Security Authority Subsystem Service (LSASS.exe) for domain controllers.	KB5022346	
Security update for Windows Server 2012 R2	Fixes various security issues for Windows Server 2012 R2.	KB5022352	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.