

# Security Patches: BD FocalPoint™ Large Lab Server

September 2023

BD has identified patches from Microsoft® that have been identified as critical or security related up to August 2023. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

## Microsoft® patches

Patch name	Description	Patch ID	Notes
Servicing stack update for Windows Server 2012	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	<a href="#">KB5029368</a>	This update must be installed prior to all other updates
Monthly rollup for Windows Server 2012 R2	Fixes various security issues for Windows Server 2012 R2.	<a href="#">KB5029312</a>	
Security update for .NET Framework 3.5	Includes cumulative reliability improvements in .NET framework 3.5.	<a href="#">KB5027116</a>	
Security update for .NET Framework 4.6.2, 4.7, 4.7.1 and 4.7.2	Includes cumulative reliability improvements in .NET framework 4.6.2, 4.7, 4.7.1 and 4.7.2.	<a href="#">KB5027112</a>	

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

[bd.com](#)

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.





## Microsoft® patches

Patch name	Description	Patch ID	Notes
Security update for Windows Server 2012 R2	Fixes various security issues for Windows Server 2012 R2.	<a href="#">KB5027282</a>	
Security update for .NET Framework 3.5	Includes cumulative security improvements in .NET framework 3.5.	<a href="#">KB5028970</a>	
Security update for .NET Framework 4.6.2, 4.7, 4.7.1 and 4.7.2	Includes cumulative reliability improvements in .NET framework 4.6.2, 4.7, 4.7.1 and 4.7.2.	<a href="#">KB5028962</a>	
Security update for Windows Server 2012 R2	Fixes various security issues for Windows Server 2012 R2.	<a href="#">KB5028223</a>	
Security update for .NET Framework 4.6.2, 4.7, 4.7.1 and 4.7.2	Includes cumulative reliability improvements in .NET framework 4.6.2, 4.7, 4.7.1 and 4.7.2.	<a href="#">KB5028977</a>	
Security update for .NET Framework 3.5	Includes cumulative security improvements in .NET framework 3.5.	<a href="#">KB5028982</a>	
Security update for Windows Server 2012 R2	Fixes various security issues for Windows Server 2012 R2.	<a href="#">KB5029304</a>	

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

**bd.com**

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.