

# Security Patches: BD FocalPoint™ Small Lab Server

December 2022

BD has identified patches from Microsoft® that have been identified as critical or security related up to November 2022. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

## Microsoft® patches

Patch name	Description	Patch ID	Notes
Servicing stack update for Windows 8.1	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	<a href="#">KB5018922</a>	This update must be installed prior to all other updates
Cumulative Security update for .NET Framework 3.5	Includes cumulative reliability improvements in .NET framework 3.5.	<a href="#">KB5016268</a>	
Cumulative Security update for Windows 8.1	Addresses miscellaneous security improvements to internal OS functionality.	<a href="#">KB5017365</a>	
Cumulative Security update for Windows 8.1	Addresses an issue that leads to User Datagram Protocol (UDP) packet drops from Linux virtual machines.	<a href="#">KB5018476</a>	

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

[bd.com](#)

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.





## Microsoft® patches

Patch name	Description	Patch ID	Notes
Cumulative Security update for Internet Explorer	Includes cumulative security updates for Microsoft Internet Explorer.	<a href="#">KB5019958</a>	
Cumulative Security update for Windows 8.1	Addresses security vulnerabilities in the Kerberos and Netlogon protocols.	<a href="#">KB5020010</a>	
Cumulative Security update for Windows 8.1	Addresses miscellaneous security improvements to internal OS functionality.	<a href="#">KB5020023</a>	
Cumulative Security update for .NET Framework 4.6.2, 4.7, 4.7.1 and 4.7.2	Includes cumulative reliability and security improvements in .NET frameworks.	<a href="#">KB5020611</a>	
Cumulative Security update for .NET Framework 4.6.2, 4.7, 4.7.1 and 4.7.2	Includes cumulative reliability and security improvements in .NET frameworks.	<a href="#">KB5020629</a>	

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

**bd.com**

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.