

Security Patches: BD FocalPoint™ Small Lab Server

April 2023

BD has identified patches from Microsoft® that have been identified as critical or security related up to January 2023. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

Microsoft® patches

Patch name	Description	Patch ID	Notes
Servicing stack update for Windows 8.1	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	KB5018922	This update must be installed prior to all other updates
Cumulative Security update for Windows 8.1	Fixes a memory leak in the Local Security Authority Subsystem Service (LSASS.exe) in Windows Domain Controllers.	KB5021296	
Cumulative Security update for Windows 8.1	Resolves multiple security issues in Windows 8.1	KB5022346	
Cumulative Security update for Windows 8.1	Resolves multiple security issues in Windows 8.1	KB5022352	

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

