

# Security Patches: BD FocalPoint™ Small Lab Server

June 2022

BD has identified patches from Microsoft® that have been identified as critical or security related for April 2022. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s)

## Microsoft® patches

Patch name	Description	Patch ID	Notes
Windows 8.1 Cumulative security update.	Addresses an issue in which Event ID 37 might be logged during certain password change scenarios.  Addresses an issue in Windows Media Center where some users might have to reconfigure the application on each start.	<a href="#">KB5012670</a>	
Service Stack Update	This update makes quality improvement to the servicing stack, which is the component that installs Windows updates.	<a href="#">KB5001403</a>	

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

**bd.com**

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

