

BD FocalPoint™ Guided Screen Review Station

Date of Critical or Security Patches: December 2021

Critical Security Patches: October 2021-November 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for October 2021 through November 2021. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Notes
Security and Quality Update for .NET Framework.	This security update resolves a vulnerability in Microsoft .NET Framework that may cause an information disclosure that allows bypassing Cross-origin Resource Sharing (CORS) configurations.	KB4480054	
Security Update for Internet Explorer	Cumulative security update for Internet Explorer	KB5006671	
Security Update for Windows 8.1	Monthly security rollup for Windows 8.1	KB5007247	

Notes

