

# BD Kiestra™ Inoqua +

## July 2020

Critical and/or Security Patches for April – June 2020

### Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for June 2020. These patches were not found to adversely affect this BD product and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Month
Security Update for Microsoft Windows 10 (1809) for X64-based Systems	Security improvements and fixes published by Microsoft between April and June 2020 for Windows 10 (1809).	<a href="#">4549949</a>	April
		<a href="#">4551853</a>	May
		<a href="#">4561608</a>	June

### Notes

1. Customers using the System Control Unit (SCU) server system: Your regional Customer Support Center will contact you to make an appointment to install the patches

