

Security Patches: BD Kiestra™ ReadA

June 2023

BD Kiestra™ ReadA

Critical and/or Security Patches: April – June 2023.

BD has identified patches from Microsoft® that have been identified as critical or security related for April – June 2023. These patches were not found to adversely affect BD products and will be applied according to customers' service agreements.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Microsoft® patches

Patch Name	Patch Description	Patch ID	Month
Security Update for Microsoft Windows 10 (1809) for X64-based Systems	Security improvements and fixes published by Microsoft for Windows 10 (1809).	KB5025229	APR-2023
		KB5026362	MAY-2023
		KB5027222	JUNE-2023
Service Stack Update (SSU) for Microsoft Windows 10 (1809) for X64-based Systems	Servicing stack updates fixes and improves the reliability of the update process to mitigate potential issues while installing.	KB5005112	APR-2023
		KB5005112	MAY-2023
		KB5005112	JUNE-2023
Security Update for Microsoft Windows Server 2016 (1607) for x64-based Systems	Security improvements and fixes published by Microsoft for Windows Server 2016.	KB5025228	APR-2023
		KB5026363	MAY-2023
		KB5027219	JUNE-2023
Service Stack Update (SSU) Update for Microsoft Windows Server 2016 (1607) for x64-based Systems	Security improvements and fixes published by Microsoft for Windows Server 2016.	KB5023788	APR-2023
		KB5023788	MAY-2023
		KB5023788	JUNE-2023

Notes

1. Customers using the System Control Unit (SCU) server system:
Your regional Customer Support Center will contact you to make an appointment to install the patches.

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

