

BD Product Name: BD Pyxis ES System

Date of Critical or Security Patches: July 2019

Abstract: Microsoft Patches report as failed but systems report fully patched

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for July 2019.

The patches include updates that may adversely affect BD Pyxis ES System on Server 2008 R2. These patches will not be and should not be installed on BD Pyxis ES System at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

Patch Name	Patch Description	Patch ID	Notes
Security Update for Microsoft Windows 2008 Server R2	<p>This security update includes quality improvements.</p> <ul style="list-style-type: none"><li>Key changes include: Security updates to Windows Server, Microsoft Graphics Component, Windows Storage and Filesystems, Windows Shell, Windows Input and Composition, and Windows Kernel.</li></ul>	<u>4507411</u>	1
Security Only Update	<p>This security update includes quality improvements. Key changes include:</p> <ul style="list-style-type: none"><li>Security updates to Windows Server, Microsoft Graphics Component, Windows Storage and Filesystems, Windows Shell, Windows Input and Composition, and Windows Kernel.</li></ul>	<u>4507456</u>	1

Notes

1. Patch fails install after deployment and installation, manual installation required.

